


NIST Standard Reference Materials® Catalog

NAT'L INST OF STAND & TECH



A11106 926652

NIST
PUBLICATIONS

NIST SP 260
JANUARY 2004

ENGINEERING MATERIALS

FOOD & AGRICULTURE

HEALTH & CLINICAL

FORENSICS

ENVIRONMENTAL

HIGH PURITY MATERIALS

INDUSTRIAL MATERIALS

PHYSICAL PROPERTIES

RADIOACTIVITY


INDUSTRIAL HYGIENE

To order:
www.nist.gov/srm
Phone: 301-975-6776
Fax: 301-948-3730
Email: srminfo@nist.gov

SRM

NIST

National Institute of
Standards and Technology
Technology Administration
U.S. Department of Commerce



QC
100
457
#260
2004
c.2

MAIL ORDERS

Mail Orders (in English) for all NIST SRMs/RMs should be directed to:

Standard Reference Materials Program
National Institute of Standards and Technology
100 Bureau Drive, Stop 2322
Gaithersburg, MD 20899-2322
USA

Telephone: (301) 975-6776
Fax: (301) 948-3730
E-Mail: srminfo@nist.gov
www.nist.gov/srm

Each purchase order should give the number of units, SRM number, and name of each reference material requested.

Example:

1 unit, SRM 930e Glass Filters for Spectrophotometry

The following information must be included with each order:

- name of customer
- shipping address
- billing address
- telephone number
- fax number
- purchase order number
- a customer identification number, i.e., a social security number (SSN) for consumer customers, tax identification number (TIN) for commercial customers, or agency code (ALC) for U.S. Government customers

Note: NIST SRMs/RMs are only distributed in the units of issue listed in this catalog and its supplement (price list). All purchase orders must be in English.

Receipt of an order does not imply acceptance of provisions set forth in the order that are contrary to the policies, practices, or regulations of the National Institute of Standards and Technology or the United States Government.

NIST SP 260 - 2004

Standard Reference Materials® Catalog

January 2004

Joylene W.L. Thomas, Editor

Yvonne A. Branden, Co-Editor

Robbin D. Howard, Co-Editor

Revised by: Regina R. Montgomery and

Joan C. Sauerwein

Standard Reference Materials Program

Technology Services

National Institute of Standards and Technology

Gaithersburg, MD 20899-2320



U.S. Department of Commerce

Donald L. Evans, Secretary

Technology Administration

Phillip J. Bond

Under Secretary of Commerce for Technology

National Institute of Standards and Technology

Arden L. Bement, Jr., Director

Please visit our website

www.nist.gov/srm

Certain commercial entities, equipment, or materials may be identified in this document in order to describe an experimental procedure or concept adequately. Such identification is not intended to imply recommendation or endorsement by the National Institute of Standards and Technology, nor is it intended to imply that the entities, materials, or equipment are necessarily the best available for the purpose.

National Institute of Standards and Technology
Special Publication 260
Supersedes NIST SP 260, 2003
138 pages (January 2004)
CODEN: NSPUE2

U.S. GOVERNMENT PRINTING OFFICE
WASHINGTON: 2004

For sale by the Superintendent of Documents,
U.S. Government Printing Office
Internet: bookstore.gpo.gov
Phone: (202) 512-1800
Fax: (202) 512-2250
Mail: Stop SSOP, Washington, DC 20402-0001

NIST Standard Reference Materials® (SRMs®) are used by industry, government, and academia to ensure the highest quality measurements. This catalog lists over 1300 individual reference materials produced and sold by NIST, each with carefully assigned values for chemical composition and physical properties.

SRMs find use in calibrating instruments and in assuring the long-term integrity of quality assurance programs. They are also key mechanisms for verifying important measurement results and in developing new measurement methods. SRMs provide users with tools to assist in establishing traceability of measurement results to NIST.

Each SRM comes carefully packaged with documentation containing assigned values with stated uncertainties and a material safety data sheet, if applicable. Details on use, stability, and NIST analytical methods are also included.

For further information and prices, contact us at:

Telephone: (301) 975-6776

Fax: (301) 948-3730

E-mail: srminfo@nist.gov

www.nist.gov/srm

■ Engineering Materials

- 1 SIZING
 - 1 Particle Size
 - 2 Cement Turbidity and Fineness
 - 2 Specific Surface Area of Powders
 - 2 Mercury Porosimetry Standards
 - 2 Particle Count Materials
- 3 SURFACE FINISH
 - 3 Abrasive Wear
 - 3 Surface Roughness
- 3 FIRE RESEARCH
 - 3 Surface Flammability
 - 4 Smoke Density Chamber
 - 4 Smoke Toxicity
 - 4 Flooring Radiant Panel
- 5 NONDESTRUCTIVE EVALUATION
 - 5 Artificial Flaw for Eddy Current NDE
- 5 PERFORMANCE ENGINEERING MATERIALS
 - 5 Fracture Toughness of Steels (Charpy V-Notch Test Blocks)
 - 5 Rockwell Hardness
 - 6 Microindentation Hardness (Knoop and Vickers Test Blocks)
 - 6 Coordinate Measuring Machine (CMM) Probe Performance
 - 6 Tape Adhesion Testing
 - 7 Bleached Kraft Pulp
 - 7 Secondary Ferrite Number (FN) Materials
 - 7 Fracture Toughness of Ceramics
 - 7 Magnetic Moment Standards

■ Food & Agriculture

- 9 Trace Elements in Food and Dairy Products
- 9 Wheat Hardness
- 10 Nutrition Composition
- 11 Trace Elements in Botanicals
- 11 Fertilizers
- 11 Whole Biomass Feedstock

■ Health & Clinical

- 13 Pure, Crystalline Standards
- 13 Biological Buffer Systems
- 14 Human Serum
- 14 Bovine Serum
- 15 DNA Profiling
- 15 Biomaterials
- 15 Toxic Substances in Urine
- 15 Miscellaneous Health-Related Materials

■ Forensics

- 17 Ethanol Solutions
- 17 DNA Profiling
- 17 Drugs of Abuse in Human Hair
- 17 Drugs of Abuse in Urine

■ Environmental

- 19 ORGANICS
 - 19 Gas Chromatography/Mass Spectrometry (GC/MS) and Characterizing Liquid Chromatography (LC) System Performance
 - 20 Organic Contaminant Calibration Solutions
 - 21 Organic Contaminants in Natural Matrix Materials
 - 22 EPA: Organic Compounds Related to Water Analysis
- 23 INORGANICS
 - 23 Metal Constituents in Natural Matrices: Air Particulate, Indoor Dust, Sediment, Mine Waste, Sludge, Soil, and Water
 - 25 Carbon Modified Silica
 - 25 Used Auto Catalysts
 - 25 Primary Gas Mixtures
- 29 FOSSIL FUELS
 - 29 Metal Constituents in Fossil Fuels
 - 29 High Purity Liquids for Fuel Rating
 - 29 Trace Elements in Coals and Coke
 - 30 Alcohols and Ethers [Oxygenates] in Reference Fuels
 - 31 Sulfur in Fossil Fuels
 - 32 Moisture in Oils and Alcohols

33 GEOLOGICAL MATERIALS AND ORES

- 33 Ores
- 34 Ore Bioleaching Substrate
- 34 Chinese Ores
- 34 Clays
- 34 Rocks and Minerals
- 35 Refractories

35 MICROANALYSIS

- 35 Metals
- 35 Synthetic Glasses
- 36 Thin Film for Transmission Electron Microscope

36 ENGINE WEAR MATERIALS

- 36 Metallo-Organic Compounds
- 37 Lubricating Base Oils
- 38 Catalyst Characterization Material
- 38 Wear-Metals in Oil

38 INDUSTRIAL HYGIENE

- 38 Materials on Filter Media
- 39 Trace Constituent Elements in Blank Filters
- 39 Respirable Silica
- 40 Lead in Paint, Dust, and Soil
- 41 Asbestos

■ High Purity Materials

- 43 Elemental Composition in High Purity Metals
- 44 Fine Gold Standards
- 44 Stoichiometric Standards
- 45 Microchemistry
- 46 Spectrometric Single Element Solutions
- 48 Anion Chromatography Solutions
- 48 Stable Isotopic Materials
- 49 Light Stable Isotopic Materials

■ Industrial Materials

51 FERROUS METALS

- 51 Steels
- 51 Plain Carbon Steels
- 52 Low Alloy Steels
- 54 Special Low Alloy Steels
- 54 High Alloy Steels
- 55 Stainless Steels
- 56 Specialty Steels
- 56 Tool Steels
- 57 Cast Steels, White Cast Irons, and Ductile Irons
- 57 Steelmaking Alloys
- 58 Cast Irons
- 59 High Temperature Alloys
- 59 Gases in Metals: Iron and Steel

TABLE OF CONTENTS

60 NONFERROUS METALS

- 60 Aluminum Base Alloys
- 60 Cobalt Base Alloys
- 61 Copper "Benchmark"
- 62 Copper Base Alloys
- 62 Lead Base Alloys
- 63 Lead Base Materials
- 63 Nickel Oxides
- 63 Nickel Base Alloys
- 64 Trace Elements in Nickel Base Superalloys
- 64 Tin Base Alloys
- 64 Titanium Base Alloys
- 65 Hydrogen in Titanium
- 65 Zirconium Base Alloys
- 65 Zinc Base Alloys

66 CERAMICS AND GLASSES

- 66 Carbides
- 66 Cemented Tungsten Carbides
- 67 Glasses
- 67 Trace Elements

68 CEMENTS

- 68 Portland Cements
- 68 Portland Cement Clinkers

69 LUBRICANTS

- 69 Metallo-Organic Compounds

Physical Properties

71 ION ACTIVITY

- 71 pH Calibration
- 72 Biological Buffer Systems
- 72 pD Calibration
- 72 Ion-Selective Electrode Calibration
- 73 Electrolytic Conductivity
- 73 Positive Electrophoretic Mobility

74 POLYMERIC PROPERTIES

- 74 Molar Mass/Molecular Weight
- 75 Melt Flow Rate
- 75 Viscosity
- 75 Biomaterials

76 THERMODYNAMIC PROPERTIES

- 76 Calorimetry - Combustion
- 76 Calorimetry - Solution
- 76 Enthalpy and Heat Capacity
- 77 Differential Scanning Calorimetry
- 77 Differential Thermal Analysis
- 77 Defining Fixed Points, International Temperature Scale of 1990, ITS-90
- 78 Reference Points
- 78 Freezing Point, Melting Point, and Triple Point Cells

- 78 Thermal Expansion of Metal and Glass

- 78 Thermal Resistance of Glass, Silica, and Polystyrene

- 79 Vapor Pressure of Metals

- 79 Thermal Conductivity of Graphite and Iron

- 79 Laboratory Thermometer

- 79 Thermocouple Material, Platinum

80 OPTICAL PROPERTIES

- 80 Molecular Transmittance and Absorbance

- 81 Transmittance Wavelength Standards

- 81 Fluorescence

- 81 Specular Spectral Reflectance

- 81 Near Infrared Reflectance Wavelength Standard

- 82 Optical Rotation

- 82 Liquid Refractive Index

- 82 X-ray and Photographic Imaging

83 ELECTRICAL PROPERTIES

- 83 Electrical Resistivity and Conductivity of Electrolytic Iron and Graphite

- 83 Electrical Resistivity and Conductivity of Silicon

84 OPTOELECTRONICS

84 METROLOGY

- 84 Optical Microscope Linewidth Measurement

- 85 Scanning Electron Microscope (SEM)

- 85 Depth Profiling

- 85 Solder Thickness for X-ray Fluorescence

- 86 Coating Thickness

- 86 Ellipsometry

- 87 Oxygen Concentration in Silicon

- 87 Superconducting Critical Current

87 CERAMICS AND GLASSES

- 87 Chemical Resistance [Durability]

- 87 Electrical Properties

- 88 Viscosity

- 88 Viscosity Fixpoints

- 88 Relative Stress Optical Coefficient

- 89 Density (glass & liquid)

- 89 Glass Liquidus Temperature

89 X-RAY SPECTROMETRY

- 89 X-ray Diffraction

- 89 X-ray Stage Calibration

Radioactivity

- 91 Radioactive Solutions
- 93 Radioactive Point Sources
- 93 Radiopharmaceuticals
- 94 Beryllium Isotopic Ratio Standard
- 94 Carbon-14 Dating
- 95 Natural Matrix Materials
- 95 Neutron Density Monitor Wire
- 95 Fission Track Glass

Industrial Hygiene

- 97 Materials on Filter Media
- 97 Trace Constituent Elements in Blank Filters
- 97 Respirable Silica
- 98 Lead in Paint, Dust, and Soil
- 99 Asbestos

Subject Index 100

Numeric Index 116



ENGINEERING MATERIALS

- 1 Sizing
- 3 Surface Finish
- 3 Fire Research
- 5 Nondestructive Evaluation
- 5 Performance Engineering Materials





SIZING

Particle Size

These SRMs are used for particle size measuring instruments, including light scattering, electrical zone flow-through counters, optical and scanning electron microscopes, sedimentation systems, and wire cloth sieving devices.

SRM	Particle Diameter (Mesh Size)	Unit Size (g)
Glass Beads, Soda Lime		
1021	2 μm to 12 μm	4
1003c	20 μm to 50 μm (No. 635 to No. 325)	28
1004b	53 μm to 125 μm (No. 270 to No. 120)	43
1017b	106 μm to 355 μm (No. 140 to No. 45)	70
1018b	250 μm to 710 μm (No. 60 to No. 25)	87
1019b	850 μm to 2000 μm (No. 20 to No. 10)	200
Sand		
RM 8010	(No. 30 to No. 325)	3 \times 150 g
Silicon Nitride (equiaxed)		
659	0.2 μm to 10 μm	5 \times 2.5 g
Zirconium Oxide (irregular)		
1978	0.2 μm to 10 μm	5
1982	10 μm to 150 μm	10
Tungsten Carbide/Cobalt (spheroidal)		
1984	9 μm to 30 μm	14
1985	18 μm to 55 μm	14
Polystyrene Spheres		
<i>Unit Size: 5 mL vial (unless otherwise noted)</i>		
1690 (0.5 % in H ₂ O)	0.895 μm	
1691 (0.5 % in H ₂ O)	0.269 μm	
1692 (0.25 % in H ₂ O)	2.982 μm	
1960* (0.4 % in H ₂ O)	9.89 μm	
1961* (0.5 % in H ₂ O)	29.64 μm	
1963** (0.5 % in H ₂ O)	0.1007 μm	
1965 (Slide Mounted: 1 slide)	9.94 μm (hexagonal array) 9.89 μm (unordered clusters)	

*Developed in cooperation with NASA

**This SRM is limited to the calibration of electron microscope and surface scanning inspection systems (not suitable for applications where monosize, unagglomerated spheres are necessary).

Cement Turbidity and Fineness

This SRM is suitable for use with ASTM C 430-92, C 115-93, and C 204-92.

SRM	Description	Properties Certified	Value	Unit Size
114p	Portland Cement	Sieve Residue (45 μm (No. 325) Sieve)	8.24 %	20 pouches \times 10 g
		Specific Surface Area (Wagner Turbidimeter)	2086 $\text{cm}^2 \cdot \text{g}^{-1}$	
		Specific Surface Area (Blaine Air Permeability)	3774 $\text{cm}^2 \cdot \text{g}^{-1}$	

Specific Surface Area (SSA) of Powders (Brunauer, Emmett, and Teller Method)

SRM	Description	Surface Area (m^2/g)		Unit Size (g)
		Multi-point	Single Point	
1897	SSA (Silica Alumina)	258.32	253.08	7
1899	SSA (Silicon Nitride)	10.67	10.52	4
1900	SSA (Silicon Nitride)	2.85	2.79	4

Mercury Porosimetry Standards


SRM	Description	Unit Size (g)
1917	Mercury Porosimetry Standard (Alumina Beads)	10
1918	Mercury Porosimetry Standard (Extruded Silica-Alumina)	12



Particle Count Materials

These SRMs are suitable for use with particle sizing instrumentation, including optical counters, in accordance with National Fluid Power Association (NFPA) T2.9.6 R2-1998 and ISO/DIS 11171.

SRM	Description	Particle Concentration	Unit Size
2806	Medium Test Dust in Hydraulic Fluid	2.8 mg/L	400 mL
RM 8631	Medium Test Dust	1 μm to 50 μm	20 g
RM 8632	Ultrafine Test Dust	1 μm to 20 μm	20 g



SURFACE FINISH

Abrasive Wear

This SRM is suitable for use with ASTM G 65, Procedure A.

SRM	Description	Unit Size
1857	D-2 Tool Steel	2 blocks: 0.78 cm × 2.5 cm × 7.6 cm

Surface Roughness

Unit Size: 25 mm × 34 mm × 12 mm

These SRMs are used for calibrating stylus instruments that measure surface roughness. These electroless-nickel coated steel blocks have a sinusoidal roughness profile machined on the top surface.

SRM	Roughness, R_a (μm)	Wavelength, D (μm)
<i>Sinusoidal Roughness (Knoop Hardness 500)</i>		
2071b	0.3137	100
2073a	0.034	100
2074	0.025	40
2075	0.012	800

FIRE RESEARCH

Surface Flammability

This SRM is suitable for checking the operation of radiant panel test equipment in accordance with ASTM E 162-78.

SRM	Description	Certification	Unit Size (cm)
1002d	Hardboard Sheet	Flame Spread Index, I = 203 Heat Evolution Factor, Q = 42.0	4 sheets: 15.2 × 45.7 × 0.6



Smoke Density Chamber

These SRMs are suitable for use with National Fire Protection Agency (NFPA) 258-1998. SRM 1006d is also suitable for use with ASTM E 662-95.

SRM	Description	Maximum Specific Optical Density (D_m (corr.))	Unit Size (cm)
1006d	Non-Flaming Exposure Condition (paper)	193	9 sheets: 17.2 × 25.4 × 0.165
1007b	Flaming Exposure Condition (plastic)	388 to 512	1 sheet: 25.4 × 25.4 × 0.076



Smoke Toxicity

SRM	Description	Combustion on Mode	Observation Time	Values		Unit Size
				LC ₅₀	N-Gas	
1048	Cup Furnace Smoke Toxicity Method Standard (ABS copolymer)	Flaming	WE*	27	1.4	8 sheets: (16 × 16 × 0.76) mm
			WE & PE**	25	1.5	
		NonFlaming	WE*	58	1.2	
			WE & PE**	53	1.4	
1049	University of Pittsburgh I Smoke Toxicity Method Standard (Nylon 6/6)		30 min exposure, plus 10 min post-exposure	4.4		150 g

*WE = within 30 minutes

**WE & PE = 30 minutes + 14 days

Flooring Radiant Panel

This SRM is suitable for use with ASTM E 648-78 and NFPA 253-1978.

SRM	Description	Critical Radiant Flux	Unit Size (cm)
1012	Flooring Radiant Panel (Kraft Paperboard)	0.36 W/cm ²	3 sheets: 104.1 × 25.4 × 0.305



NONDESTRUCTIVE EVALUATION

Artificial Flaw for Eddy Current NDE

RM	Description	Flaw Size	Unit Size
8458	Artificial Flaw (Aluminum Alloy)	3.0 mm × 0.1 mm	7 cm × 7 cm × 2 cm

PERFORMANCE ENGINEERING MATERIALS

Fracture Toughness of Steels (Charpy V-Notch Test Blocks)

Unit Size: set of 10 mm × 10 mm × 54 mm specimens

These SRMs are suitable for use with ASTM E 23 and ISO/DIS 12736.

SRM	Description	Energy Range (J)
2092	Low Energy (4340 Alloy Steel)	13 to 20
2096	High Energy (4340 Alloy Steel)	88 to 136
2098	Super High Energy (Maraging Steel)	176 to 244



Rockwell Hardness

Unit size: 60 mm diameter × 15 mm

SRM	Description	Nominal Hardness (HRC)
2810	Rockwell C Scale Hardness - Low Range	25
2811	Rockwell C Scale Hardness - Mid Range	45
2812	Rockwell C Scale Hardness - High Range	62

Microindentation Hardness (Knoop and Vickers Test Blocks)

Unit Size: 1.15 cm × 1.15 cm (unless otherwise noted)

These SRMs are suitable for use with ASTM E 384.

SRM	Description	Load (N)	Hardness (kg/mm ²)
Copper, Bright			
1893	Knoop	0.245, 0.49, 0.98	125
Nickel, Bright			
1894a	Vickers	0.245, 0.49, 0.98	125
1895	Knoop	0.245, 0.49, 0.98	600
1896a	Vickers	0.245, 0.49, 0.98	600
1905	Knoop	2.943	600
1906	Knoop	4.905	600
1907	Knoop	9.81	600
1908	Vickers	2.943	500
1909	Vickers	9.81	500
2798a	Vickers	4.905	600
Silicon Nitride, Ceramic			
2830 (22 mm diameter × 9.54 mm)	Knoop	19.6	1500
2831 (25 mm diameter × 9.5 mm)	Vickers	9.8	1530

Coordinate Measuring Machine (CMM) Probe Performance

These SRMs are suitable for use with ANSI/ASME B89.4.1.

SRM	Description	Unit Size
2084	Tungsten Carbide Sphere	10 mm diameter (stem mounted with stand)
2084R	Tungsten Carbide Sphere	10 mm diameter (stem mounted)
2085	Stainless Steel Sphere	25 mm diameter (stem mounted)

Tape Adhesion Testing


This SRM is suitable for use with ASTM D 2860 and ASTM D 3654.

SRM	Description	Unit Size
1810a	Linerboard for Tape Adhesion Testing	50 sheets: 21.6 cm × 28 cm



Bleached Kraft Pulps

These RMs are intended primarily for use in fundamental studies on the physical properties of fibers and paper sheets. No extensive property measurements have been made on these materials beyond ensuring that they were within the control limits of the normal production run.

RM	Description	Unit Size	
8495*	Northern Softwood	10 standard lap sheets: 0.5 kg each	
8496*	Eucalyptus Hardwood	10 standard lap sheets: 0.5 kg each	

*Developed in cooperation with the Pulp Material Research Committee

Secondary Ferrite Number (FN) Materials

The RMs are suitable for use with ANSI/AWS A4.2 and ISO 8249.

RM	Ferrite Number	Unit Size (mm)
8480	0 to 30	10 × 12 × 20
8481	30 to 120	10 × 12 × 20

Fracture Toughness of Ceramics

Unit Size: 3 mm × 4 mm × (45 to 47) mm

SRM	Description	Fracture Toughness (MPa · m ^{1/2})	No. of Specimens
2100	Silicon Nitride Flexure Specimens	4.57	5

Magnetic Moment Standards

SRM	Description	Certified Property	Unit Size
762	Nickel Disk	Specific Magnetization	disk: 6 mm diameter × 0.13 mm
772a	Nickel Sphere	Magnetic Moment	sphere: 2.383 mm diameter sphere
2853	Yttrium Garnet Sphere	Magnetic Moment	sphere: 1 mm diameter (2.8 mg)

FOOD & AGRICULTURE

- 9 Trace Elements in Food and Dairy Products
- 9 Wheat Hardness
- 10 Nutrition Composition
- 11 Trace Elements in Botanicals
- 11 Fertilizers
- 11 Whole Biomass Feedstock





Trace Elements in Food and Dairy Products



FOOD & AGRICULTURE

SRM	Description	Unit Size (g)
1577b	Bovine Liver	50
RM 8414*	Bovine Muscle Powder	50
RM 8413*	Corn Kernel	47
RM 8412*	Corn Stalk	34
RM 8436*	Durum Wheat Flour	50
RM 8437*	Hard Red Spring Wheat Flour	50
1549	Non-fat Milk Powder	100
1566b	Oyster Tissue	25
1568a	Rice Flour	80
RM 8438*	Soft Winter Wheat Flour	50
1570a	Spinach Leaves	60
1548a	Typical Diet	2 × 6.5 g
1567a	Wheat Flour	80

* Developed by Agriculture Canada in cooperation with NIST

Wheat Hardness


Unit Size: 50 × 20 g

RM	Description	Wheat Numbers
8441*	Wheat Hardness	Hard-1 through Hard-5 Soft-1 through Soft-5

Nutrition Composition

Please visit our website to view the relevant certificate or report of investigation for available certified and non-certified values.

NIST Food-Matrix SRMs and RMs



1. SRM 1563
2. SRM 2384
3. SRM 2387
4. SRM 1546
RM 8415
5. SRM 2383
RM 8432
RM 8433
RM 8436
6. SRM 1846
RM 8435
SRM 1548a
SRM 1544
7. SRM 1566b
SRM 1570a
SRM 2385
9. SRM 1946
SRM 1947
SRM 1974a
RM 8418

SRM	Description	Certified Constituents **	Unit Size (g)
1544	Fatty Acids and Cholesterol in Frozen Diet Composite	Cholesterol, Fatty Acids, Proximates	4 × 15 g
1546	Meat Homogenate	Cholesterol, Fatty Acids, Proximates, Vitamins, Minerals	4 × 85 g
1548a	Typical Diet	Proximates, Trace Elements, Total Dietary Fiber	2 × 6.5 g
1563	Cholesterol and Fat-Soluble Vitamins in Coconut Oil	Cholesterol, Ergocalciferol, dl-α-Tocopheryl Acetate	10 ampoules: 5 fortified, 5 natural
1589a	PCBs, Pesticides, and Dioxins/Furans in Human Serum	Cholesterol, Triglycerides	5 × 10 mL
1845	Whole Egg Powder	Cholesterol	35
1846	Infant Formula (milk-based)	Minerals, Proximates, Vitamins, Fatty Acids	10 × 30 g
2383	Baby Food Composite	Carotenoids, Cholesterol, Minerals, Proximates, Vitamins	4 × 70 g
RM 8415*	Whole Egg Powder	Fatty Acids, Minerals, Proximates, Vitamins	35
RM 8418*	Wheat Gluten	Fatty Acids, Minerals, Proximates, Vitamins	50
RM 8432*	Corn Starch	Fatty Acids, Minerals, Proximates, Vitamins	50
RM 8433*	Corn Bran	Fatty Acids, Minerals, Proximates, Vitamins	50
RM 8435*	Whole Milk Powder	Fatty Acids, Minerals, Proximates, Vitamins	40
RM 8436*	Durum Wheat Flour	Fatty Acids, Minerals, Proximates, Vitamins	50
1570a	Spinach Leaves	Fatty Acids, Trace Elements, Proximates, Total Dietary Fibers	60
2384	Baking Chocolate	Fat, Fatty Acids, Calcium, Iron, Caffeine, Theobromine, Catechins	5 × 91 g
1566b	Oyster Tissue	Fatty Acids, Nitrogen, Proximates, Total Dietary Fiber, Trace Elements, Mercury, Methylmercury	25
1974a	Organics in Mussel Tissue (<i>Mytilus Edulis</i>)	Selected Proximate	3 × 15 g
2385	Slurried Spinach	Calcium, Iron, Magnesium, Manganese, Phosphorous, Potassium, Zinc, Lutein, Beta Carotene	3 × 70 g
2387	Peanut Butter	Fat, Fatty Acids, Elements, Tocopherols	170 g
3240	Ephedra-Related Dietary Supplement	Alkaloids	5 g × 10
3241	Ephedra-Related Dietary Supplement	Alkaloids	1.2 g × 10
3242	Ephedra-Related Dietary Supplement	Alkaloids	1.2 g × 10
3243	Ephedra-Related Dietary Supplement	Alkaloids	2.5 g × 10
3244	Ephedra-Related Dietary Supplement	Alkaloids	12 g × 10



* Developed by Agriculture Canada in cooperation with NIST

** Proximates are provided as reference values.



Trace Elements in Botanicals

SRM	Description	Unit Size (g)
1515	Apple Leaves	50
1547	Peach Leaves	50
1570a	Spinach Leaves	60
1575a	Pine Needles	50
1573a	Tomato Leaves	50
2695*	Fluoride in Vegetation	2 × 25 g
RM 8412	Corn Stalk (Zea Mays)	34
RM 8413	Corn Kernel (Zea Mays)	47



* Developed in cooperation with Aluminum Association, Inc.

Fertilizers

Unit Size: 90 g

SRM	Description	Certified Constituents
120c	Phosphate Rock (Florida)	Minerals
193	Potassium Nitrate	N, K
194	Ammonium Dihydrogen Phosphate	N, P
200a	Potassium Dihydrogen Phosphate	K, P
694	Phosphate Rock (Western)	Minerals

Whole Biomass Feedstock*


Unit Size: 5 × 10 g


RM	Description	Reference Constituents
8491	Sugarcane Bagasse	Ash, Ethanol Extractives, Acid Soluble Lignin,
8492	Populus Deltoides	Acid Insoluble Lignin, Total Lignin, Glucuronic Acid,
8493	Monterey Pine	Arabinan, Xylan, Mannan, Galactan, Glucan
8494	Wheat Straw	

* Developed by the International Energy Agency (IEA) Biomass Annex, NREL, and NIST

HEALTH & CLINICAL

- 13 Pure, Crystalline Standards
- 13 Biological Buffer Systems
- 14 Human Serum
- 14 Bovine Serum
- 15 DNA Profiling
- 15 Biomaterials
- 15 Toxic Substances in Urine
- 15 Miscellaneous Health-Related Materials





Pure, Crystalline Standards

SRM	Description	Purity (%)	Unit Size (g)
998	Angiotensin I (Human)	94.1	0.5
916a	Bilirubin	98.3	0.1
915a	Calcium Carbonate	99.9	20
911b	Cholesterol	99.8	2
921	Cortisol (Hydrocortisone)	98.9	1
914a	Creatinine	99.7	10
917b	D-Glucose (Dextrose)	99.7	50
920	D-Mannitol	99.8	50
937	Iron Metal (Clinical)	99.90	50
928	Lead Nitrate	100.00	30
924a	Lithium Carbonate	99.867	30
929	Magnesium Gluconate Dihydrate	5.403 Mg	5
918a	Potassium Chloride	99.9817	30
919a	Sodium Chloride	99.89	30
910	Sodium Pyruvate	98.7	25
1595	Tripalmitin	99.5	2
912a	Urea	99.9	25
913a	Uric Acid	99.6	10
925	VMA (4-hydroxy-3-methoxy-DL-mandelic acid)	99.4	1

**Values in parentheses are not certified and are given for information only.*

Biological Buffer Systems


Unit Size: 60 grams

SRM	Description	pH(S) Values (at 37 °C)	
		0.05 molal	0.08 molal
2181	HEPES Free Acid	7.364*	7.373*
2182	NaHEPESate		
2183	MOPS0 Free Acid	6.699*	6.694*
2184	NaMOPS0ate		



**This pH results only when the two SRMs listed are used as an admixture in solution.*

Human Serum



SRM	Description	Certified Constituents	Reference Values Constituents	Form	No. of Levels
1599	Anticonvulsant Drug Level Assay (Valproic Acid and Carbamazepine)	PCB Congeners (16), Chlorinated Pesticides (5), Total Cholesterol	PCB Congeners (9), Chlorinated Pesticides (5), Total Cholesterol, Triglycerides, "Free" Cholesterol, Phospholipids	Lyophilized	1
900	Antiepilepsy Drug Level Assay	Antiepileptics (4)	---	Lyophilized	3
970	Ascorbic Acid in Frozen Human Serum	Total Ascorbic Acid	---	Frozen	2
1952a	Cholesterol in Human Serum	Cholesterol	---	Lyophilized	3
956a	Electrolytes in Frozen Human Serum	Total Ca, Li, Mg, K, Na	Ionized Ca, Cl	Frozen	3
968c	Fat-Soluble Vitamins, Carotenoids, and Cholesterol in Human Serum	Vitamins (4), Cholesterol, Carotenoids (4)	Carotenoids (8), Vitamin D	Lyophilized	2
965a	Glucose in Human Serum	Glucose	---	Frozen	3
909b	Human Serum	Organics (6), Inorganics (6)	Bilirubin	Lyophilized	2
1951b	Lipids in Frozen	Total Cholesterol, Total Glycerides Triglycerides	HDL-, LDL-, and Total Cholesterol, Triglycerides, Free Glycerol		
1589a	PCBs, Pesticides, Dioxins/Furans in Serum	PCB Congeners (16), Chlorinated Pesticides (5), Total Cholesterol	PCB Congeners (9), Chlorinated Pesticides (5), Total Cholesterol (5), Triglycerides, "Free" Cholesterol, Phospholipids	Lyophilized	1

Bovine Serum

SRM	Description	Certified Constituents	Reference Constituents	Form	No. of Levels
927c	Bovine Serum Albumin (7 % Solution)	Protein Concentration	---	Solution	1
1598	Inorganic Constituents in Bovine Serum	Elements (13)	---	Frozen	1
955b	Lead in Bovine Blood	Pb	---	Frozen	4
966	Toxic Elements in Bovine Blood	Pb, Cd	Pb, Cd, Total Hg, Inorganic Hg	Frozen	2



DNA Profiling

SRM 2392 is intended to provide quality control when performing the polymerase chain reaction (PCR) and sequencing of human mitochondrial DNA (mtDNA) for medical diagnosis, or mutation detection. It may also be used as a control when amplifying (PCR) and sequencing any DNA.


SRM	Description	Unit Size
2392	Mitochondrial DNA Sequencing (Human)	3 components



Biomaterials

SRM	Description	Certified Properties	Reference Properties	Unit Size
2910	Calcium Hydroxyapatite	Calcium Phosphorus Ca/P Molar Ratio Specific Surface Area Solubility Product		5 g (powder)
RM 8456	Ultra High Molecular Weight Polyethylene		Young's Modulus Yield Strength Ultimate Strength Elongation	3 in diameter × 60 in (bar) (7.62 cm diameter × 152.4 cm)

Toxic Substances in Urine



SRM	Description	No. of Levels	Unit Size
2670a	Toxic Elements in Urine	2	4 × 20 mL
2671a	Fluoride	2	4 × 20 mL
2672a	Mercury	2	4 × 20 mL

Miscellaneous Health-Related Materials

SRM	Description	Certified Constituents	Form	Unit Size
2389	Amino Acids in 0.1 mol/L HCl	Amino Acids (17)	Solution	5 ampoules
2921	Cardiac Troponin	Troponin I	Solution	—
1400	Bone Ash	Elements (8)	Powdered	50 g
1486	Bone Meal	Elements (8)	Powdered	50 g

FORENSICS

17 Ethanol Solutions

17 DNA Profiling

17 Drugs of Abuse
in Human Hair

17 Drugs of Abuse
in Urine



Ethanol Solutions

This SRM is for use in the calibration of instruments and techniques for the determination of ethanol (ethyl alcohol) in breath and blood.




SRM	Description	Ethanol Mass Fraction (%)	Unit Size
1828b	Ethanol - Water Solutions (Set: 5 ampoules)	95.60	1 × 20 mL
		1.9957	2 × 20 mL
		0.09480	1 × 5 mL
		0.02186	1 × 5 mL

DNA Profiling

SRMs 2390, and 2391a are intended for use in the standardization of forensic and paternity quality assurance procedures and instructional law enforcement or non-clinical research purposes.

SRM	Description	Unit Size
2390	DNA Profiling Standard - RFLP	20 components
2391b	PCR-Based DNA Profiling Standard	12 components



Drugs of Abuse in Human Hair

SRM	Description	Certified Constituents
2379	Drugs of Abuse in Human Hair I	6
2380	Drugs of Abuse in Human Hair II	4

Drugs of Abuse in Urine

SRM	Description	Certified Constituents	Reference Constituent	Form	Unit Size
1508	Cocaine Metabolites in Urine	Benzoylcegonine		Lyophilized	3 levels, plus 1 blank
RM 8444	Cotinine in Urine		Cotinine (nicotine metabolite)	Lyophilized	2 levels, plus 1 blank
1507b	Marijuana Metabolites in Urine	TH-9-COOH		Lyophilized	3 levels, plus 1 blank
2381	Morphine and Codeine in Urine	Morphine and Codeine		Lyophilized	3 levels, plus 1 blank
2382	Morphine Glucuronide in Urine	Free Morphine		Lyophilized	3 levels, plus 1 blank
1511	Multi Drugs of Abuse in Urine	Drugs of Abuse (5)		Lyophilized	1 level



ENVIRONMENTAL

19 Organics

23 Inorganics

29 Fossil Fuels


33 Geological Materials and Ores

35 Microanalysis

36 Engine Wear Materials

38 Industrial Hygiene





ORGANICS


Gas Chromatography/Mass Spectrometry (GC/MS) and Characterizing Liquid Chromatography (LC) System Performance


SRM	Description	Certified Constituents	Solvent	No. of Levels	Unit Size
1543	GC/MS System	Methyl Stearate, Benzophenone	Hexane	2	4 × 1 mL
RM 8443	Consists of 5 units of SRM 1543				
877	LC Chiral Selectivity	various Chiral components	Ethanol	—	5 × 1 mL
870	LC Performance	Silanol Activity, Trace Metal Activity, Hydrophobic Retention, Methylene Selectivity	Methanol	1	5 × 1 mL
869a	LC Selectivity	Shape Selectivity: PAHs (3)	Acetonitrile	1	5 × 1 mL



Organic Contaminant Calibration Solutions

SRM	Description	Certified Constituents	Non-Certified Constituents	Unit Size
RM 8467	4,4'-DDE (neat)	—	—	Vial: 100 mg
RM 8469	4,4'-DDT (neat)	—	—	Vial: 100 mg
RM 8466	g-HCH (Lindane) (neat)	—	—	Vial: 100 mg
1491	Aromatic Hydrocarbons in Hexane/Toluene	PAHs (23)	PAHs (1)	5 ampoules
2260a	Aromatic Hydrocarbons in Toluene	PAHs (23)	PAHs (1)	5 ampoules
1493	Chlorinated Biphenyl Congeners in 2,2,4-Trimethylpentane	PCBs (18)	PCBs (2)	5 ampoules
2262	Chlorinated Biphenyl Congeners in 2,2,4-Trimethylpentane	PCBs (25)	PCBs (4)	5 ampoules
2275	Chlorinated Pesticide Solution-II	Pesticides (9)	—	5 ampoules
1492	Chlorinated Pesticides in Hexane	Pesticides (15)	—	5 ampoules
2261	Chlorinated Pesticides in Hexane	Pesticides (15)	—	5 ampoules
2273	DDTs and Metabolites in Solution	DDTs, Metabolites (7)	5 ampoules	
1596	Dinitropyrene Isomers and 1-Nitropyrene in Methylene Chloride	Nitro-PAHs (4)	—	5 ampoules
1614	Dioxin (2,3,7,8-TCDD) in Iso-octane	Dioxins (2)	Dioxins (2)	6 ampoules
1639	Halocarbons (in Methanol) for Water Analysis	Halocarbons (7)	—	5 ampoules
1586	Isotopically Labeled and Unlabeled Priority Pollutants in Methanol	Priority pollutants (10)	—	6 ampoules
1587	Nitrated PAHs in Methanol	Nitro-PAHs (6)	Nitro-PAHs (1)	4 ampoules
2274	PCB Congener Solution-II	PCBs (11)		5 ampoules
2269	Perdeuterated PAH-I	Perdeuterated PAHs (5)	—	5 ampoules
2270	Perdeuterated PAH-II	Perdeuterated PAHs (6)	—	5 ampoules
1647d	Priority Pollutant PAHs (in Acetonitrile)	PAHs (16)	—	5 ampoules
1584	Priority Pollutant Phenols in Methanol	Phenols (10)	Phenols (1)	5 ampoules
2276	Three Planar PCBs in Solution	PCBs (3)	—	5 ampoules





Organic Contaminants in Natural Matrix Materials

SRM	Description	Certified Constituents	Non-Certified Constituents	Unit Size
1597	Complex Mixture of PAHs from Coal Tar	PAHs (12)	PAHs/PASH/PANH (18)	4 ampoules
1975	Diesel Particulate Extract	PAHs (8)	PAHs (29), Nitro-PAHs	4 ampoules
1650b	Diesel Particulate Material	PAHs (5), Nitro-PAHs (1)	PAHs (6), Nitro-PAHs (3), PAQ (1)	0.01 g
2975	Diesel Particulate Matter	PAHs (11) (Industrial Forklift)	PAHs (28), Total Extractable Mass, Particle Size Distribution	1 g
2978	Mussel Tissue (Organic Contaminants - Raritan Bay, NJ)	PAHs (7), PCB Congeners (22), Pesticides (12)	PAHs (20), PCBs (2)	10 g
2977	Mussel Tissue	PAHs (14), PCB Congeners (25), Pesticides (7), Trace Elements (6), Methylmercury	PAHs (16), Trace Elements (9)	10 g
2976	Mussel Tissue	Methylmercury, Total Mercury, Trace Elements (7)	Trace elements (20)	25 g
1941b	Organics in Marine Sediment	PAHs (24), PCBs (29), Pesticides (7)	PAHs (43), PCBs (13), Pesticides (2), Tin Species (3) and Total Tin	50 g
1944	New York/New Jersey Waterway Sediment	PAHs (24), PCBs (35), Pesticides (4), Trace Elements (9)	PAHs (32), Pesticides (7), Trace Elements (20), PCDDs/PCDFs (17), Particle Size, Total Organic Carbon	50 g
1946	Lake Superior Fish Tissue	PCBs (30), Pesticides (15), Fat and Fatty Acids (14), Total Mercury, Methylmercury, Arsenic, Iron	PCBs (12), Pesticides (2), Fatty Acids (12), Proximates, Caloric Content, Trace Elements (9)	5 × 7 – 9 g
1588a	Organics in Cod Liver Oil	PCBs (24), Pesticides (4)	PCDDs/PCDFs (7), PCBs (34), Pesticides (3)	5 ampoules
1974a	Organics in Mussel Tissue (<i>Mytilus Edulis</i>) (Frozen)	PAHs (15), PCBs (20), Pesticides (7), Total Mercury, Methylmercury	Aliphatics (16), Trace Elements (32), PAHs (18), PCBs (4), Pesticides (4), Proximates, Calories	3 × 15 g
1580	Organics in Shale Oil	PAHs (5), Phenols (3), PANH (1)	Phenols (6), PANH (1)	5 ampoules
1945	Organics in Whale Blubber (Frozen)	PCBs (27), Pesticides (15)	PCBs (2), Pesticides (2)	2 bottles
1589a	PCBs, Pesticides, and Dioxins/Furans in Human Serum	PCBs Congeners (16), Pesticides (5), Total Cholesterol	CDC Lipid Laboratory: PCB Congeners (9), Pesticides (5), Total Cholesterol, Triglycerides, "Free" Cholesterol, Phospholipids, PCDDs, PCDFs, non-ortho, PCBs	5 × 10 mL
1582	Petroleum Crude Oil	PAHs (5), PASH (1)	PAHs (5), Phenols (2), PANH (1)	5 ampoules

(continued)


Organic Contaminants in Natural Matrix Materials (continued)

SRM	Description	Certified Constituents	Non-Certified Constituents	Unit Size
1939a	Polychlorinated Biphenyls (Congeners) in River Sediment	PCBs (20), Pesticides (3)	PCBs (4)	50 g
1649a	Urban Dust	PAHs (22), PCBs (35), Pesticides (8)	PAHs (22), Pesticide (1, Mutagenic Activity), PCDD/PCDFs (17), Trace Elements (32), Particle size, Total Organic Carbon	2.5 g
1648	Urban Particulate Matter	Trace Elements (9)	Trace Elements (25), PAH (13)	2 g

EPA: Organic Compounds Related to Water Analysis

These SRMs are intended primarily for the calibration of instrumentation and validation of methods for volatile or semi-volatile organic compound determinations. Because of its miscibility with water, each SRM can also be used to fortify aqueous samples with known amounts of the organic compound. These SRMs were developed by the NIST Analytical Chemistry Division (ACD) primarily to support the Chemical Calibration Providers of the Proficiency Testing Program with support by the U.S. Environmental Protection Agency (EPA).

SRM	Description	Unit Size
3000	Benzene in Methanol	2 × 2.5 mL
3001	Toluene in Methanol	2 × 2.5 mL
3002	Ethylbenzene in Methanol	2 × 2.5 mL
3003	o-Xylene in Methanol	2 × 2.5 mL
3004	m-Xylene in Methanol	2 × 2.5 mL
3005	p-Xylene in Methanol	2 × 2.5 mL
3006	Carbon Tetrachloride in Methanol	2 × 2.5 mL
3008	Methylene Chloride in Methanol	2 × 2.5 mL
3009	1,2-Dichloropropane in Methanol	2 × 2.5 mL
3010	Tetrachloroethene (Tetrachloroethylene) in Methanol	2 × 2.5 mL
3011	1,1,1-Trichloroethane in Methanol	2 × 2.5 mL
3012	1,2-Dichloroethane in Methanol	2 × 2.5 mL
3014	1,2,3-Trichloropropane in Methanol	2 × 2.5 mL
3015	Isopropylbenzene in Methanol	2 × 2.5 mL
3016	sec-Butylbenzene in Methanol	2 × 2.5 mL
3063	2, 3, 7, 8-Tetrachlorodibenzo-p-dioxin (2, 3, 7, 8-TCDD) in Methanol	5 × 1.2 mL
3064	Endothall in Water	5 × 1.2 mL



EPA: Organic Compounds Related to Water Analysis (continued)

SRM	Description	Unit Size
3067	Toxaphene in Methanol	5 × 1.2 mL
3068	Chlordane in Methanol	5 × 1.2 mL
3071	Glyphosate in Water	5 × 1.2 mL
3072	Diquat Dibromide Monohydrate in Water	5 × 1.2 mL
3075	Aroclor 1016 in Transformer Oil	5 × 1.2 mL
3076	Aroclor 1232 in Transformer Oil	5 × 1.2 mL
3077	Aroclor 1242 in Transformer Oil	5 × 1.2 mL
3078	Aroclor 1248 in Transformer Oil	5 × 1.2 mL
3079	Aroclor 1254 in Transformer Oil	5 × 1.2 mL
3080	Aroclor 1260 in Transformer Oil	5 × 1.2 mL
3081	Aroclor 1016 in Methanol	5 × 1.2 mL
3082	Aroclor 1232 in Methanol	5 × 1.2 mL
3083	Aroclor 1242 in Methanol	5 × 1.2 mL
3084	Aroclor 1248 in Methanol	5 × 1.2 mL
3085	Aroclor 1254 in Methanol	5 × 1.2 mL
3086	Aroclor 1260 in Methanol	5 × 1.2 mL
3090	Aroclors in Transformer Oil (set SRMs 3075-3080)	6 × 1.2 mL
3091	Aroclors in Methanol (set SRMs 3081 - 3086)	5 × 1.2 mL

INORGANICS

Metal Constituents in Natural Matrices: Air Particulate, Indoor Dust, Sediment, Mine Waste, Sludge, Soil, and Water

SRM	Description	Elements	Unit Size
<i>Air Particulate</i>			
2783	Air Particulate on Filter Media	18 certified 9 reference	2 filters, plus 2 blanks
1648	Urban Particulate Matter	15 certified	2 g
<i>Indoor Dust, Trace Elements in</i>			
2583	Nominal 90 mg/kg Lead	5 certified	8 g
2584	Nominal 1 % Lead	5 certified 10 reference	8 g



(continued)

Metal Constituents in Natural Matrices: Air Particulate, Indoor Dust, Sediment, Mine Waste, Sludge, Soil, and Water (continued)

Sediment

RM 8704	Buffalo River Sediment	25 reference	50 g
1646a	Estuarine Sediment	20 certified	70 g
1944	New York/New Jersey Waterway Sediment	72 certified 78 reference	50 g
1946	Lake Superior Fish Tissue	3 certified 9 reference	5 × 7 - 9 g
2702	Marine Sediment	25 certified 8 reference	50 g
2703	Sediment for Solid Sampling	—	5 g

Mine Waste and Sludge

2780	Hard Rock Mine Waste	12 certified 7 reference	50 g
2781	Domestic Sludge	10 certified	40 g
2782	Industrial Sludge	10 certified 16 reference	70 g
RM8785	Particulate matter on Filters	3 reference	3 filters/unit


Soil, Trace Elements in

2710	Montana Soil Highly Elevated Trace Element Concentrations		21 certified 50 g
2711	Montana Soil Moderately Elevated Trace Element Concentrations		24 certified 50 g
2709	San Joaquin Soil	26 certified	50 g
2586	Nominal 500 mg/kg Lead	4 certified 18 reference	8 g
2587	Nominal 3000 mg/kg Lead	4 certified 14 reference	8 g
2780	Hard Rock Mine Waste	12 certified 7 reference	50 g

Water

1641d	Mercury in Water	1 certified	10 × 10 mL
1640	Natural Water	17 certified 10 reference	250 mL
1643e	Trace Elements in Water	—	250 mL





Carbon Modified Silica

Unit Size: 3 × 1 g


This SRM is chemically modified microparticulate silica intended for the calibration of instruments used to measure total carbon.

SRM	Description	Bottle	Mass Fraction (%)
1216	Carbon Modified Silica	I	0.70
		II	9.06
		III	17.04

Used Auto Catalysts

Unit Size: 70 g

SRM	Description	Elemental Composition
2557	Recycled Monolith	Pt, Pd, Rh, Pb
2556	Recycled Pellet	



Primary Gas Mixtures

These SRMs are supplied in a DOT 3AL specification aluminum (6061 alloy) cylinder with a nominal pressure exceeding 12.4 MPa that provides the user with approximately 0.73 m³ of usable mixture.

SRM	Nominal Amount-of-Substance ($\mu\text{mol/mol}$)
<i>Ambient Non-Methane Organics in Nitrogen (15 components in large cylinder)</i>	
1800	5 nmol/mol
1804c	5 nmol/mol
<i>Carbon Dioxide in Air (Certified for CO₂)</i>	
1671a	340
1672a	350
1676	365


(continued)


Primary Gas Mixtures (continued)

SRM	Nominal Amount of Substance Fraction ($\mu\text{mol/mol}$)
Carbon Monoxide in Air (Certified for CO)	
2612a	10
2613a	20
2614a	45
Carbon Dioxide in Nitrogen (Certified for CO₂)	
1674b*	7 mol %
1675b*	14 mol %
2619a	0.5 mol %
2620a	1.0 mol %
2621a	1.5 mol %
2622a	2.0 mol %
2623a	2.5 mol %
2624a	3.0 mol %
2625a*	3.5 mol %
2626a	4.0 mol %
2745*	16 mol %
Carbon Monoxide in Nitrogen (Certified for CO)	
1677c*	10
1678c*	50
1679c*	100
1680b*	500
1681b*	1000
2635a*	25
2636a*	250
2637a*	2500
2638a*	5000
2639a	1 mol %
2640a	2 mol %
2641a	4 mol %
2642a*	8 mol %

*Available as a NIST Traceable Reference Material (NTRM); from commercial suppliers.
A suppliers list is available on our website.

(continued)





Primary Gas Mixtures (continued)

SRM	Nominal Amount of Substance Fraction ($\mu\text{mol/mol}$)
Carbon Monoxide in Nitrogen (Certified for CO) continued	
2740a	10 mol %
2741a	13 mol %
Hydrogen Sulfide in Nitrogen (Certified for H₂S)	
2730	5
2731	20
Methane in Air (Certified for CH₄)	
1658a	1
1659b	10
1660a (also certified for C ₃ H ₈)	4 (methane) 1 (propane)
2750	50
2751	100
Nitric Oxide in Nitrogen (Certified for NO)	
1683b*	50
1684b*	100
1685b*	250
1686b*	500
1687b*	1000
2629a*	20
2630*	1500
2631a*	3000
2735	800
2736a	2000
2737	0.5
2738	1.0
Oxides of Nitrogen in Air (Certified for NO_x)	
2660a*	100

*Available as a NIST Traceable Reference Material (NTRM); from commercial suppliers.
A suppliers list is available on our website.

(continued)

Primary Gas Mixtures (continued)



SRM	Nominal Amount of Substance Fraction ($\mu\text{mol/mol}$)
<i>Oxygen in Nitrogen (Certified for O_2)</i>	
2657a*	2 mol %
2658a*	10 mol %
2659a*	21 mol %
<i>Propane in Air (Certified for CH_4)</i>	
1660a (also certified for C_3H_8)	4 (methane) 1 (propane)
1665b	3
1666b	10
1667b	50
1668b*	100
1669b	500
2764	0.25
<i>Propane in Nitrogen (Certified for C_3H_8)</i>	
2643a	100
2644a	250
2645a	500
2646a	1000
2647a	2500
2648a	5000
<i>Sulfur Dioxide in Nitrogen (Certified for SO_2)</i>	
1661a*	500
1662a*	1000
1663a*	1500
1664a*	2500
1693a*	50
1694a*	100
1696a*	3500

*Available as a NIST Traceable Reference Material (NTRM); from commercial suppliers. A suppliers list is available on our website.

The gas NTRM program was established in 1992 in partnership with the U.S. EPA and specialty gas companies as a means for providing end users with the wide variety of certified gas standards needed to implement the Emissions Trading Provision of the 1990 Clean Air Act.



FOSSIL FUELS

Metal Constituents in Fossil Fuels

SRM	Pb Concentration	Unit Size
2713	Lead in Reference Fuels (19.4 µg/g Pb)	6 × 20 mL
2714	Lead in Reference Fuels (28.1 µg/g Pb)	6 × 20 mL
1634c	Trace Elements in Fuel Oil "No. 6" (As, Co, Ni, Pb, S, Se, V)	100 mL
RM 8505	Vanadium in Crude Oil	250 mL

High Purity Liquids for Fuel Rating

Unit Size: 100 mL

SRM	Description	Purity (%)
1816a	Isooctane (2,2,4-Trimethylpentane)	99.987
1815a	n-Heptane	99.987



Trace Elements in Coals and Coke

SRM	Description	Constituents	Unit Size (g)
2719	Calcined Petroleum Coke	6 certified, 2 reference	50
1632c	Coal (Bituminous)	15 certified, 25 reference	50
1635	Coal (Subbituminous)	16 certified	75
1633b	Coal Fly Ash	23 certified	75
2689	Coal Fly Ash	11 certified	3 × 10 g
2690	Coal Fly Ash	11 certified	3 × 10 g
2691	Coal Fly Ash	11 certified	3 × 10 g
2718	Green Petroleum Coke	6 certified, 2 reference	50



Alcohols and Ethers [Oxygenates] in Reference Fuels



SRM	Description	Constituents	Unit Size
Alcohols in Gasoline			
1829	Alcohols (t-Butanol, Ethanol, Methanol)	4 certified	6 × 20 mL
1838	Ethanol	1 certified	5 × 20 mL
2285	Arson Text Mixture	15 certified	5 × 1.2 mL
2286	Ethanol	2 certified	3 × 20 mL
2287	Ethanol	2 certified	3 × 20 mL
1839	Methanol	1 certified	5 × 20 mL
1837	Methanol	2 certified	5 × 20 mL
Ethers in Gasoline			
<i>Unit Size: 3 × 20 mL</i>			
2288	t-Amyl Methyl Ether	2 certified	
2289	t-Amyl Methyl Ether	2 certified	
2290	Ethyl t-Butyl Ether ETBE	2 certified	
2291	Ethyl t-Butyl Ether ETBE	2 certified	
2292	Methyl t-Butyl Ether MTBE	2 certified	
2293	Methyl t-Butyl Ether MTBE	2 certified	
Ethers and Ethanol in Reformulated Gasoline			
<i>Unit Size: 2 × 20 mL</i>			
2294	11 % MTBE	4 certified 26 reference	
2295	15 % MTBE	4 certified 26 reference	
2296	13 % ETBE	4 certified 26 reference	
2297	10 % Ethanol	4 certified 26 reference	



Sulfur/Mercury in Fossil Fuels

SRM	Description	% S	Hg (µg/kg)
Coke Foundry			
<i>Unit Size: 50 g</i>			
2775	Foundry Coke	0.5816	—
2776	Foundry Coke	0.825	—
Diesel Fuel Oil			
<i>Unit Size: 10 × 10 mL</i>			
2723a	Sulfur in Diesel Fuel Oil	0.00110	—
2724b	Sulfur in Diesel Fuel Oil	0.04304	—
Gasolines			
2294	Reformulated Gasoline (nominal 11 % MTBE) (2 × 20 mL)	0.00409	—
2295	Reformulated Gasoline (nominal 15 % MTBE) (2 × 20 mL)	0.0308	—
2296	Reformulated Gasoline (nominal 13 % ETBE) (2 × 20 mL)	0.00400	—
2297	Reformulated Gasoline (nominal 10 % Ethanol) (2 × 20 mL)	0.03037	—
2298	Reformulated Gasoline (5 × 20 mL)	0.00047	—
2299	Gasoline (High Octane) (5 × 20 mL)	0.00136	—
Kerosine			
<i>Unit Size: 100 mL</i>			
1616a	Sulfur in Kerosine	0.01462	—
1617a	Sulfur in Kerosine	0.17307	—
Petroleum Coke			
<i>Unit Size: 50 g</i>			
2719	Trace Elements in Calcined Petroleum Coke	0.8877	—
2718	Trace Elements in Green Petroleum Coke	4.7032	—



Sulfur/Mercury in Fossil Fuels (continued)

SRM	Description	% S	Hg (µg/kg)
Residual Fuel Oil			
<i>Unit Size: 100 mL</i>			
1619b	Sulfur in Residual Fuel Oil	0.6960	—
1620c	Sulfur in Residual Fuel Oil	4.561	—
1621e	Sulfur in Residual Fuel Oil	0.9480	—
1622e	Sulfur in Residual Fuel Oil	2.1468	—
1623c	Sulfur in Residual Fuel Oil	0.3806	—
2717a	Sulfur in Residual Fuel Oil	2.9957	—
Crude Oil			
<i>Unit Size: 5 × 10 mL</i>			
2721	Crude Oil	1.5832	0.0525
2722	Crude Oil	0.21037	0.1441
Coals			
<i>Unit Size: 50 g (unless otherwise noted)</i>			
2683b	Sulfur and Mercury in Coal	1.955	90.0
2684b	Sulfur and Mercury in Coal	3.076	97.4
2685b	Sulfur and Mercury in Coal	4.730	146.2
2692b	Sulfur and Mercury in Coal	1.170	133.3
2693	Low Sulfur/Mercury Coal	0.4567	37.3
2682b	Sulfur and Mercury in Coal (Subbituminous)	0.4917	108.8
1632c	Trace Elements in Coal Bituminous	1.462	93.8
1635	Trace Elements in Coal (Subbituminous) (75 g)	0.3616	10.9

Moisture in Oils and Alcohols

SRM	Description	Unit Size (mL)
RM 8509	Methanol	5 mL
RM 8507	Mineral Oil	10 mL
RM 8510	Moisture in Methanol	5 mL
RM 8506a	Transformer Oil	5 × 9.5 mL
2890	Water Saturated 1-Octanol	5 × 2 mL

GEOLOGICAL MATERIALS AND ORES

Elements in Ores

MINING



SRM	Description	Certified Constituents	Unit Size (g)
699	Alumina (Reduction Grade)	13	60
69b	Bauxite, Arkansas	15	60
697	Bauxite, Dominican	15	60
698	Bauxite, Jamaican	15	60
696	Bauxite, Surinam	15	60
1835	Borate Ore	15	60
330	Copper Ore Mill Heads	3	100
331	Copper Ore Mill Tails	3	100
79a	Fluorspar, Customs Grade	1	120
180	Fluorspar, High Grade	1	120
886	Gold Ore, Refractory	2	200
670	Iron Ore, Canada	6	90
690	Iron Ore, Canada	11	100
692	Iron Ore, Labrador	11	100
693	Iron Ore, Nimba	11	100
691	Iron Oxide, Reduced	13	100
182	Lithium Ore (Petalite)	1	45
181	Lithium Ore (Spodumene)	1	45
183	Lithium Ore (Lepidolite)	1	45
25d	Manganese Ore	8	60
120c	Phosphate Rock, Florida	14	90
694	Phosphate Rock, Western	13	90
600	Rutile Ore	16	90
2430	Scheelite Ore	6	100
277	Tungsten Concentrate	1	45
113b	Zinc Concentrate	9	100

ENVIRONMENTAL

Ore Bioleaching Substrate

This RM is for use as a bioleaching substrate and for testing bioleaching rates.

RM	Description	Unit Size (g)
8455	Pyrite Ore	100

Elements in Chinese Ores

Unit Size: 100 g

These RMs are a well characterized series (more than 50 elements and minerals) of skarn deposit ores developed and certified by the Hubei Geological Research Laboratory, Hubei Province, China.

RM	Description
8600	Copper
8601	Copper
8602	Lead
8603	Lead
8605	Molybdenum
8606	Molybdenum
8607	Tungsten
8608	Tungsten
8604	Zinc



COPPER WIRE

Elements in Clay

SRM	Description	Certified Constituents	Unit Size (g)
679	Brick Clay	12	75
97b	Flint Clay	12	60
98b	Plastic Clay	12	60

Elements in Rocks and Minerals

SRM	Description	Certified Constituents	Unit Size (g)
688	Basalt Rock	12	60
70a	Feldspar, Potash	10	40
99a	Feldspar, Soda	11	40

(continued)



Elements in Rocks and Minerals (continued)



SRM	Description	Certified Constituents	Unit Size (g)
81a	Glass Sand	5	75
165a	Glass Sand (Low Iron)	4	75
1413	Glass Sand (High Alumina)	9	75
1c	Limestone, Argillaceous	12	50
88b	Limestone, Dolomite	11	75
278	Obsidian Rock	12	35

Elements in Refractories

SRM	Description	Certified Constituents	Unit Size (g)
76a	Burnt Refractory (Al ₂ O ₃ -40 %)	11	75
77a	Burnt Refractory (Al ₂ O ₃ -60 %)	11	75
78a	Burnt Refractory (Al ₂ O ₃ -70 %)	11	75
198	Silica Brick	11	45
199	Silica Brick	11	45
154c	Titanium Dioxide	1	90

MICROANALYSIS

Elements in Metals

SRM	Description	Certified Constituents	Unit Size
482	Gold-Copper Wires for Microprobe Analysis	2	wires: 6
481	Gold-Silver Wires for Microprobe Analysis	2	wires: 6
480	Tungsten-20 % Molybdenum Alloy Electron Microprobe Standard	2	rod: 1

Elements in Synthetic Glasses

SRM	Description	Certified Constituents	Unit Size
1873	Barium-Zinc-Silicate Glasses for Microanalysis (K-458, K-489, K-963)	2	rod: 2 mm × 2 mm × 20 mm
2066	Glass Microspheres (K-411)	4 certified 1 reference	glass microspheres: 50 mg
1872	Lead-Germanate Glasses for Microanalysis (K-453, K-491, K968)	2	rod: 2 mm × 2 mm × 20 mm

Thin Film for Transmission Electron Microscope

SRM	Description	Certified Element	Unit Size
2063a	Microanalysis Thin Film Mineral Glass	Ar, Ca, Fe, Mg, O, Si	1 glass film



ELECTRON MICROSCOPE

ENGINE WEAR MATERIALS

Metallo-Organic Compounds

Unit Size: 5 g

These SRMs are for preparing solutions in oils of known and reproducible concentrations of metals.

SRM	Description	Elemental Composition
1075a	Aluminum 2-Ethylhexanoate	8.07 Al
1051b	Barium Cyclohexanebutyrate	28.7 Ba
1080a	Bis (1-phenyl-1,3-butanediono)copper (II)	16.37 Cu
1052b	Bis(1-phenyl-1,3-butanediono)oxovanadium (IV)	13.01 V
1053a	Cadmium Cyclohexanebutyrate	24.8 Cd
1057b	Dibutyltin bis (2-ethylhexanoate) (tin)	22.95 Sn
1059c	Lead Cyclohexanebutyrate	37.5 Pb
1060a	Lithium Cyclohexanebutyrate	4.1 Li
1065b	Nickel Cyclohexanebutyrate	13.89 Ni
1066a	Octaphenylcyclotetrasiloxane	14.14 Si
1077a	Silver 2-Ethylhexanoate	42.60 Ag

(continued)

**Metallo-Organic Compounds** (continued)

SRM	Description	Elemental Composition
1069b	Sodium Cyclohexanebutyrate	12.0 Na
1070a	Strontium Cyclohexanebutyrate	20.7 Sr
1071b	Triphenyl Phosphate	9.48 P
1078b	Tris (1-phenyl-1,3-butanediono)chromium (III)	9.6 Cr
1079b	Tris (1-phenyl-1,3-butanediono)iron (III)	10.45 Fe
1073b	Zinc Cyclohexanebutyrate	16.66 Zn

**Lubricating Base Oils**

These SRMs are for determining the concentrations of a single element in lubricating base oil. SRMs 1818a and 1819a consist of five bottles, approximately 20 g of liquid each; SRM 1836 consists of four sets of four ampoules, each ampoule containing approximately 4 g of liquid.

SRM	Description	Elemental Composition (mg/kg)				
		I	II	III	IV	V
1818a	Total Chlorine	31.6	60.0	78.2	154.4	234.0
1836	Total Nitrogen	9.0	50.9	113.3	166.2	
1819a	Total Sulfur	423.5	741.1	4022	4689	6135

Catalyst Characterization Material

This RM is for determining the activity of FCC Catalysts by Microactivity Test and is distributed by NIST in cooperation with ASTM.

RM	Description	Unit Size
8590	High Sulfur Gas Oil Feed	946 mL

Wear-Metals in Oil

SRM	Description	Unit Size
1848	Lubricating Oil Additive Package	100 mL
1084a	Wear-Metals	5 × 1.6 g
1085b	Wear-Metals	5 × 1.2 g
1083	Wear-Metals (Base Oil)	150 mL

INDUSTRIAL HYGIENE

Materials on Filter Media

These SRMs consist of potentially hazardous materials deposited on filters to be used to determine the levels of these materials in industrial atmospheres.

SRM/RM	Description	Set Size	Elemental Composition	Diameter (mm)	Pore Size (μm)
2679a	Quartz on Filter Media	2 × 3 levels, plus 2 blanks	Quartz, Clay	47	0.45
2783	Air Particulate on Filter	2 filters, plus 2 blanks	18 certified values 9 reference values	47	0.4
8785	Particulate Matter on Filters	3 filters	1 reference value 2 information values	37	—



Trace Constituent Elements in Blank Filters

SRMs 2678 and 2681 are for use in evaluating the performance of air sampling filter methods with either certified values (in μg) or limits of detection (X_D) for each of 30 constituent elements, as well as six leachable anions and cations.

SRM	Description	Diameter (mm)	Pore Size (μm)	Filter Weight (g)
2678	Cellulose Acetate Membrane	47	0.45	0.09
2681	Ashless Blank Filter	42.5	—	0.14

Respirable Silica

These SRMs are intended for use in determining, by X-ray diffraction, the levels of respirable silica in an industrial atmosphere according to the National Institute for Occupational Safety and Health (NIOSH) Analytical Method 7500 or equivalent methods.

SRM	Description	Mass Fraction/Mass Loading	Unit Size
1878a	Respirable Alpha Quartz	100.00% \pm 0.21%	5 g
1879a	Respirable Cristobalite	95.6% \pm 0.4%	5 g
2950	Respirable Alpha Quartz on Filter Media	(10, 20, 50, 100, 250, 500) $\mu\text{g}/\text{filter}$	set SRMs 2952-57
2951	Respirable Alpha Quartz on Filter Media	5 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2952	Respirable Alpha Quartz on Filter Media	10 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2953	Respirable Alpha Quartz on Filter Media	20 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2954	Respirable Alpha Quartz on Filter Media	50 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2955	Respirable Alpha Quartz on Filter Media	100 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2956	Respirable Alpha Quartz on Filter Media	250 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2957	Respirable Alpha Quartz on Filter Media	500 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2958	Respirable Alpha Quartz on Filter Media	1000 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2960	Respirable Alpha Cristobalite on Filter Media	(5, 10, 20, 50, 100, 250) $\mu\text{g}/\text{filter}$	set SRMs 2961-66
2961	Respirable Alpha Cristobalite on Filter Media	5 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2962	Respirable Alpha Cristobalite on Filter Media	10 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2963	Respirable Alpha Cristobalite on Filter Media	20 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2964	Respirable Alpha Cristobalite on Filter Media	50 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2965	Respirable Alpha Cristobalite on Filter Media	100 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2966	Respirable Alpha Cristobalite on Filter Media	250 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)
2967	Respirable Alpha Cristobalite on Filter Media	500 $\mu\text{g}/\text{filter}$	5 filters (5 blanks)

Lead in Paint, Dust, and Soil

These SRMs and RM have been developed in conjunction with the U.S. EPA to monitor paint, dust, and soil sources of lead.

SRM	Lead Concentration	Unit Size
Paint Film		
2570	<0.001 mg/cm ²	1 blank film
2571	3.58 mg/cm ²	1 film, plus 1 blank
2572	1.527 mg/cm ²	1 film, plus 1 blank
2573	1.040 mg/cm ²	1 film, plus 1 blank
2574	0.714 mg/cm ²	1 film, plus 1 blank
2575	0.307 mg/cm ²	1 film, plus 1 blank
2579a (Set of 6: SRMs 2570 to 2575)	0.307 to 3.58 mg/cm ²	5 films, plus 1 blank
2576 (High Level)	5.59 mg/cm ²	1 film, plus 1 blank
Powdered Paint		
2580	4.34 %	30 g
2581	0.449 %	35 g
2582	209.8 mg/kg	20 g
2589	9.99 %	35 g
Indoor Dust, Trace Elements in (As, Cd, Cr, Hg, Pb)		
2583	85.9 mg/kg	8 g
2584	9761 mg/kg	8 g
Soil, Trace Elements in		
2586	432 mg/kg	50 g
2587	3242 mg/kg	50 g
Paint on Fiberboard		
RM 8680	1 to 2 mg/cm ²	1 sheet: (10.2 × 15.2 × 1.3) cm



Asbestos

SRM	Description	Asbestos Type	Unit Size
1866b	Common Commercial Asbestos	chrysotile grunerite (Amosite) riebeckite (Crocidolite)	3 × 4 g
1868	Quantitative Asbestos in Building Material		
1876b	Chrysotile Asbestos for TEM	—	10 sections: 3 mm × 3 mm
RM 8411	Mixed Asbestos Research Filter	chrysotile asbestos grunerite (Amosite)	1 cm ²



HIGH PURITY MATERIALS


- 43 Elemental Composition
in High Purity Metals
- 44 Fine Gold Standards
- 44 Stoichiometric Standards
- 45 Microchemistry
- 46 Spectrometric Single
Element Solutions
- 48 Anion Chromatography
Solutions
- 48 Stable Isotopic Materials
- 49 Light Stable Isotopic Materials





HIGH PURITY MATERIALS

Elemental Composition in High Purity Metals



SRM	Description	Unit Size	
685R	High Purity Gold	rod:	5.9 mm diameter × 25 mm
685W	High Purity Gold	wire:	1.4 mm diameter × 102 mm
680a(L1)	High Purity Platinum	wire:	0.51 mm diameter × 10 cm
680a(L2)	High Purity Platinum	wire:	0.51 mm diameter × 1 m
682	High Purity Zinc	semicirc:	57 mm
885	Refined Copper	pin:	200 g
726	Selenium, Intermediate Purity	shot:	450 g
683	Zinc Metal	semicirc:	57 mm
728	Zinc, Intermediate Purity	shot:	450 g

Fine Gold Standards

These RMs are a series of fine gold and gold bullion products developed and certified by the Royal Canadian Mint (RCM), Ottawa, Canada and distributed by NIST. The fine gold RMs are primarily intended for use as calibration standards for the determination of trace elements by solid sample spectrometric methods; the gold bullion RMs are primarily intended for use as quality control check standards for fire assay. There are five sets of RMs in the gold bullion series (RMs 8068-8082) available in three forms: disc (25 mm diameter \times 20 mm); wire (2 mm diameter); and foil (35 mm \times 40 mm \times 1 mm). There are six sets of RMs in the fine gold series (RMs 8050-8067) available in three forms: block (25 mm \times 25 mm \times 2.5 mm); wire (2 mm diameter); and turnings (25 g).

Stoichiometric Standards

These SRMs are defined as primary, working, and secondary standards in accordance with recommendations of the Analytical Chemistry Section of the International Union of Pure and Applied Chemistry [Ref. Analyst 90, 251 (1965)]. These definitions are as follows:

- Primary Standard: a commercially available substance of purity $100\% \pm 0.02\%$ (Purity 99.98+ %)
- Working Standard: a commercially available substance of purity $100\% \pm 0.05\%$ (Purity 99.95+ %)
- Secondary Standard: a substance of lower purity which can be standardized against a primary grade standard

SRM/RM	Description	Certified Use	Stoichiometric Purity (%)	Unit Size (g)
951	Boric Acid	Acidimetric and Boron Isotopic Value	100.00	100
84k	Potassium Hydrogen Phthalate	Acidimetric Standard	99.9911	60
350a	Benzoic Acid	Acidimetric Standard	99.9958	30
351	Sodium Carbonate	Acidimetric Standard	99.9796	50
723d	Tris(hydroxymethyl)aminomethane	Acidimetric Standard	99.924	50
987	Strontium Carbonate	Assay and Isotopic Values	99.98	1
999a	Potassium Chloride	Assay Values for: 1. Potassium Chloride 2. Potassium 3. Chloride	99.9817 52.4354 47.5463	60
136e	Potassium Dichromate	Oxidimetric Standard	99.984	60
17e	Sucrose	Polarimetric Standard	99.950	60
917b	D-Glucose (Dextrose)	Polarimetric Standard	99.7	50
8040	Sodium Oxalate	Reductometric Standard	99.972	60
83d	Arsenic Trioxide	Reductometric Standard	99.9926	60

Microchemistry

Unit Size: 2 g



SRM	Description	Certified Component
141d	Acetanilide	C, H, N, O
142	Anisic Acid	CH ₃ O-
143d	Cystine	C, H, N, S, O
2144	m-Chlorobenzoic Acid	Cl
148	Nicotinic Acid	C, H, N
2143	p-Fluorobenzoic Acid	F
2141	Urea	N

HIGH PURITY MATERIALS

Spectrometric Single Element Solutions

Unit Size: 50 mL

These SRMs are intended as standard solutions for use in calibrating instruments used in atomic spectrometry, including atomic absorption spectrometry, inductively coupled plasma optical spectrometry, and inductively coupled plasma mass spectrometry.

SRM	Element	Nominal Acid Concentration
3101a	Aluminum	HNO ₃ 10 %
3102a	Antimony	HNO ₃ 10 % + HF 2 %
3103a	Arsenic	HNO ₃ 15 %
3104a	Barium	HNO ₃ 1 %
3105a	Beryllium	HNO ₃ 10 %
3106	Bismuth	HNO ₃ 10 %
3107	Boron	H ₂ O
3108	Cadmium	HNO ₃ 10 %
3109a	Calcium	HNO ₃ 10 %
3110	Cerium	HNO ₃ 10 %
3111a	Cesium	HNO ₃ 1 %
3112a	Chromium	HNO ₃ 10 %
3113	Cobalt	HNO ₃ 10 %
3114	Copper	HNO ₃ 10 %
3115a	Dysprosium	HNO ₃ 10 %
3116a	Erbium	HNO ₃ 10 %
3117a	Europium	HNO ₃ 16 %
3118a	Gadolinium	HNO ₃ 10 %
3119a	Gallium	HNO ₃ 10 %
3120a	Germanium	HNO ₃ 10 % + HF 2 %
3121	Gold	HNO ₃ 5 % + HF 2 %
3122	Hafnium	HNO ₃ 10% + HF 2%
3123a	Holmium	HNO ₃ 16 %
3124a	Indium	HNO ₃ 10 %
3126a	Iron	HNO ₃ 10 %
3127a	Lanthanum	HNO ₃ 10 %
3128	Lead	HNO ₃ 10 %
3129a	Lithium	HNO ₃ 1 %
3130a	Lutetium	HNO ₃ 10 %
3131a	Magnesium	HNO ₃ 10 %
3132	Manganese	HNO ₃ 10 %

(continued)



Spectrometric Single Element Solutions (continued)

SRM	Element	Nominal Acid Concentration
3133	Mercury	HNO ₃ 10 %
3134	Molybdenum	HCl 10 %
3135a	Neodymium	HNO ₃ 10 %
3136	Nickel	HNO ₃ 10 %
3137	Niobium	HNO ₃ 10 % + HF 2 %
3138	Palladium	HCl 10 %
3139a	Phosphorus	HNO ₃ 0.8 %
3140	Platinum	HCl 10 %
3141a	Potassium	HNO ₃ 1 %
3142a	Praseodymium	HNO ₃ 10 %
3143	Rhenium	HNO ₃ 10 %
3144	Rhodium	HCl 10 %
3145a	Rubidium	HNO ₃ 1 %
3147a	Samarium	HNO ₃ 10 %
3148a	Scandium	HNO ₃ 10 %
3149	Selenium	HNO ₃ 10 %
3150	Silicon	H ₂ O
3151	Silver	HNO ₃ 10 %
3152a	Sodium	HNO ₃ 1 %
3153a	Strontium	HNO ₃ 10 %
3154	Sulfur	H ₂ SO ₄ 0.1 %
3155	Tantalum	HNO ₃ 10 % + HF 2 %
3156	Tellurium	HCl 20 %
3157a	Terbium	HNO ₃ 16 %
3158	Thallium	HNO ₃ 10 %
3159	Thorium	HNO ₃ 10 %
3160a	Thulium	HNO ₃ 10 %
3161a	Tin	HNO ₃ 5 % + HF 2 %
3162a	Titanium	HNO ₃ 10 % + HF 2 %
3163	Tungsten	HNO ₃ 7 % + HF 4 %
3164	Uranium	HNO ₃ 10 %
3165	Vanadium	HNO ₃ 10 %
3166a	Ytterbium	HNO ₃ 16 %
3167a	Yttrium	HNO ₃ 10 %
3168a	Zinc	HNO ₃ 10 %
3169	Zirconium	HNO ₃ 10 % + HF 2 %

HIGH PURITY MATERIALS

Anion Chromatography Solutions

Unit Size: 50 mL

These SRMs are single component solutions prepared gravimetrically for use in anion chromatography or any other technique that requires aqueous standard solutions for calibration of control materials.

SRM	Description	Nominal Concentration (mg/kg)
3184	Bromide	1000
3182	Chloride	1000
3183	Fluoride	1000
3185	Nitrate	1000
3186	Phosphate	1000
3181	Sulfate	1000

Stable Isotopic Materials

SRM	Description	Chemical Form	Unit Size (g)
951	Boron Isotope Standard	Boric Acid	100
952	Enriched ^{10}B Isotope Standard	Boric Acid	0.25
975a	Chlorine Isotope Standard	Sodium Chloride	0.25
976	Copper Isotope Standard	Metal	disk: 0.4
977	Bromine Isotope Standard	Sodium Bromide	0.25
978a	Silver Isotope Standard	Silver Nitrate	0.25
979	Chromium Isotope Standard	Chromium Nitrate	0.25
980	Magnesium Isotope Standard	Metal	0.25
981	Lead Isotope Standard, Natural	Metal	wire: 1.0
982	Lead Isotope Standard, $^{208}\text{Pb}/^{206}\text{Pb}$ Equal Atom	Metal	wire: 1.0
983	Lead Isotope Standard, Radiogenic	Metal	wire: 1.0
984	Rubidium Isotope Standard	Rubidium Chloride	0.25
985	Potassium Isotope Standard	Potassium Chloride	1.0
986	Nickel Isotope Standard	Metal	0.5
987	Strontium Isotope Standard	Strontium Carbonate	1.0
991	Nitrate Spike Isotope Standard, ^{206}Pb	Nitric Acid	15
994	Gallium Isotope Standard	Metal	disk: 0.25
997	Thallium Isotope Standard	Metal	rod: 0.25
3230	Iodine-129, Isotopic (low levels)	Iodine	5 × 5 mL (plus blank)
3231	Iodine-129, Isotopic (high levels)	Iodine	5 × 5 mL (plus blank)



Light Stable Isotopic Materials

These RMs are distributed by NIST on behalf of the International Atomic Energy Agency (IAEA). At the request of the IAEA, quantities of these materials are limited to *one unit of each RM per laboratory every 3 years*.

Isotopic Ratio Legend:

1. D / H

2. ^{18}O / ^{16}O

3. ^{13}C / ^{12}C

4. ^6Li / ^7Li

5. ^{30}Si / ^{28}Si

6. ^{15}N / ^{14}N

7. ^{34}S / ^{32}S

RM	Description	Isotopic Ratios	Unit Size
8535	VSMOW-Water	1,2	20 mL
8536	GISP-Water	1,2	20 mL
8537	SLAP-Water	1,2	20 mL
8538	NBS30-Biotite	1,2,3	2 g
8539	NBS22-Oil	1,2,3	1 mL
8540	PEFI-Polyethylene	1,2,3	~2 mg
8541	USGS24-Graphite	1,2,3	0.8 g
8542	Sucrose ANU-Sucrose	1,2,3	1 g
8543	NBS18-Carbonatite	2,3	0.4 g
8544	NBS18-Limestone	2,3	0.4 g
8545	LSVEC-Lithium Carbonate	3,4	0.4 g
8546	NBS28-Silica Sand (Optical)	2,5	0.4 g
8547	IAEA-N1-Ammonium Sulfate	6	0.4 g
8548	IAEA-N2-Ammonium Sulfate	6	0.4 g
8549	IAEA-N3-Potassium Nitrate	6	0.4 g
8550	USGS25-Ammonium Sulfate	6	0.4 g
8551	USGS26-Ammonium Sulfate	6	0.4 g
8552	NSVEC-Gaseous Nitrogen	6	300 μmol
8553	Soufre de Lacq - Elemental Sulfur	2,7	0.5 g
8554	IAEA-S1-Silver Sulfide	2,7	0.5 g
8555	IAEA-S2-Silver Sulfide	2,7	0.5 g
8556	NBS123-Sphalerite	2,7	0.5 g
8557	NBS127-Barium Sulfate	2,7	0.5 g
8558	USGS32-Potassium Nitrate	6	0.5 g
8559	Natural Gas Isotopic	—	1 cylinder (0.1 mole)
8560	Natural Gas Isotopic	—	1 cylinder (0.1 mole)
8561	Natural Gas Isotopic	—	1 cylinder (0.1 mole)
8562	CO ₂ -Heavy, Paleomarine Origin	2,3	2 tubes: 9 mm diameter \times 300 mm
8563	CO ₂ -Light, Petrochemical Origin	2,3	2 tubes: 9 mm diameter \times 300 mm
8564	CO ₂ -Biogenic, Modern Biomass Origin	2,3	2 tubes: 9 mm diameter \times 300 mm

HIGH PURITY MATERIALS

INDUSTRIAL MATERIALS


51 Ferrous Metals

60 Nonferrous Metals

66 Ceramics and Glasses

68 Cements

69 Lubricants





FERROUS METALS

Steels

These SRMs consist of selected steel alloys that provide a wide range of analytical values for relevant elements. Please visit our website to view the relevant certificate or report of investigation for all available certified and non-certified values. These RMs are a series of skarn deposit ores developed and certified by the Hubei Geological Research Laboratory, Hubei Province, China.

Plain Carbon Steels (chip)

Unit Size: 150 g (unless otherwise noted)

SRM	Description
178	0.4C Basic Oxygen Furnace Steel
13g	0.6 % Carbon Steel
20g	AISI 1045 Steel
14g	AISI 1078 Carbon Steel
368	AISI 1211 Steel
19h	Basic Electric Steel, 0.2 % Carbon

Basic Open-Hearth Steel

15h	0.1 % Carbon
11h	0.2 % Carbon
12h	0.4 % Carbon
152a	0.5 % Carbon (Tin-Bearing)
337a	1 % Carbon (300 g)



Low Alloy Steels (disk and rod)

Nominal Sizes for Solid Steel SRMs:


600 Series: 3.2 mm diameter \times 51 mm

1100 and 1200 Series: 31 mm diameter \times 19 mm

1700 Series: 34 mm diameter \times 19 mm

A "C" preceding the SRM number indicates a chill cast sample; 31 mm diameter \times 19 mm.

SRM	Description
1270	2-1/4 Chromium - 1 Molybdenum Low Alloy Steel, A 336 (F-22)
C1285	A242, Modified
1224	AISI 1078, Carbon Steel
C1221	AISI 1211, Modified, Resulturized/Rephosphorized
1269	AISI 1526, Modified (Line Pipe Steel)
1225	AISI 4130
661	AISI 4340
1262b	AISI 94B17 (Modified)
1254	Calcium in Low Alloy Silicon Steel
663	Chromium-Vanadium Steel, Modified
1263a	Chromium-Vanadium Steel, Modified
1265a	Electrolytic Iron
664	High Carbon Steel, Modified
1264a	High Carbon Steel, Modified
1135	High Silicon Steel
1134	High Silicon Steel
1768	High Purity Iron
1226	HY 130
1286	HY 80
1228	Basic Open Hearth Steel (0.1 % Carbon)
1227	Basic Open Hearth Steel (1 % Carbon)




(continued)



Low Alloy Steels (disk and rod) (continued)

SRM	Description
-----	-------------

1761	Low Alloy Steel
1762	Low Alloy Steel
1763	Low Alloy Steel
1764	Low Alloy Steel
1765	Low Alloy Steel
1766	Low Alloy Steel
1767	Low Alloy Steel



Low Alloy Steels (chip)

Unit Size: 150 g (unless otherwise noted)

SRM	Description
-----	-------------

72g	AISI 4130
293	AISI 8620 (Cr - Ni - Mo)
139b	AISI 8640 (Cr - Ni - Mo)
291	ASTM A213 (Cr - Mo)
163	Chromium Steel (100 g)
36b	Chromium-Molybdenum Steel
155	Chromium-Tungsten Steel
129c	SAE 112 High Sulfur
2171	HSLA 100 (6Ni - Cr - Cr - Cu - Mo)
106b	Nitralloy™ G (Cr - Mo - Al)
32e	SAE 3140 (Ni - Cr)
100b	SAE 340 (Mn)
33e	SAE 4820 (Ni)
30f	SAE 6150 (Cr - V)

Silicon Steels

179	High Silicon Steel
125b	High Silicon Steel, Calcium-Bearing
131g	Low Carbon Silicon Steel

Special Low Alloy Steels (chip and pin)


Unit Size: 150 g (unless otherwise noted)

SRM	Description
2159	Low Alloy Steel (pin - 200 g)
2160	Low Alloy Steel (pin - 200 g)
2166	Low Alloy Steel
2167	Low Alloy Steel
361	AISI 4340 Steel
362	AISI 94B17, Modified
363	Chromium-Vanadium Steel, Modified
364	High Carbon Steel, Modified
2168	High Purity Iron

High Alloy Steels (chip)

Unit Size: 150 g (unless otherwise noted)

SRM	Description
345a	Cu Precipitation Hardening Steel (15Cr - 4Ni)
344	Mo Precipitation Hardening Steel (15Cr - 7Ni)
126c	High Nickel Steel (36 % Ni)
868	High Temperature Alloy (Fe-Ni-Co) (100 g)
348a	High Temperature Alloy A286 (Ni-Cr)
862	High Temperature Alloy L605 (100 g)
346a	Valve Steel





Stainless Steels (chip)

Unit Size: 150 g (unless otherwise noted)



SRM	Description
339	SAE 303Se (17Cr - 9Ni - 0.2Se)
101g	AISI 304 L (18Cr - 10Ni)
343a	AISI 431 (16Cr - 2Ni)
123c	AISI 348 (17Cr - 11Ni - 0.6Nb)
121d	AISI 321 (17Cr - 11Ni - 0.3Ti)
160b	AISI 316 (18Cr - 12Ni - 2Mo)
166c	AISI 316L Low Carbon Stainless Steel (100 g)
893	SAE 405 (Cr)
895	SAE 201 (Cr-Mn)
73c	SAE 420 (13 % Cr)

Stainless Steels (disk)

Unit Size: 32 mm diameter × 19 mm

SRM	Description
1219	AISI 431 (16Cr - 2Ni)
1172	AISI 348 (17Cr - 11Ni - 0.6Nb)
1223	Chromium Steel
1297	SAE 201
1295	SAE 405
C1296	SAE 460
C1153a	(17Cr - 9Ni)
C1152a	(18Cr - 11Ni)
1155	AISI 316 (18Cr - 12Ni - 2Mo)
C1154a	Stainless Steel, (19Cr - 13Ni)
C1151a	Stainless Steel, (23Cr - 7Ni)
1171	AISI 321 (17 Cr - 11Ni - 0.3Ti)

Specialty Steels (disk)

SRM	Description	Unit Size
1158	High Nickel Steel, 36 % Nickel	32 mm diameter × 19 mm
1772	S-7 Tool Steel	34 mm diameter × 19 mm
1157	AISI M2, Tool Steel	32 mm diameter × 19 mm
1233	Valve Steel	35 mm diameter × 19 mm



Tool Steels (chip)

Unit Size: 150 g

SRM	Description
134a	Molybdenum - Tungsten - Chromium - Vanadium Steel
2172	S-7 Tool Steel
132b	AISI M2, Tool Steel
50c	Tungsten - Chromium - Vanadium Steel






Cast Steels, White Cast Irons, and Ductile Irons (disk)

Unit Size: 32 mm diameter × 19 mm

SRM	Description
1138a	Cast Steel (No. 1)
1139a	Cast Steel (No. 2)
C1173	Cast Steel (No. 3)
C2423	Ductile Iron A
C2423a	Ductile Iron B
C2424	Ductile Iron C
C2424a	Ductile Iron D
C1291	High Alloy White Cast Iron, Ni-Hard, Type I
C1292	High Alloy White Cast Iron, Ni-Hard, Type IV
C1290	High Alloy White Cast Iron, HC-250+V
1173	Nickel-Chromium-Molybdenum-Vanadium Steel
C1137a	White Cast Iron
C1145a	White Cast Iron

Steelmaking Alloys (fine powder)

Unit Size: 150 g




SRM	Description
57a	Silicon Metal
58a	Ferrosilicon (73 % Silicon-Regular Grade)
59a	Ferrosilicon
64c	High Carbon Ferrochromium
68c	High Carbon Ferromanganese
90	Ferrophosphorus
195	Ferrosilicon (75 % Silicon High Purity Grade)
196	Low Carbon Ferrochromium
347	Magnesium Ferrosilicon
689	Silicon Ferrochromium

Cast Irons (chip)

Unit Size: 150 g

SRM	Title
4L	Cast Iron
5m	Cast Iron
6g	Cast Iron
122i	Cast Iron
7g	High Phosphorus Cast Iron
115a	Copper-Nickel-Chromium Cast Iron
341	Ductile Cast Iron
334	Gray Cast Iron (Carbon and Sulfur)
890	High-Alloy White Cast Iron, HC 250+V
891	High-Alloy White Cast Iron, Nickel-Hard, Type I
892	High-Alloy White Cast Iron, Nickel-Hard, Type IV
82b	Nickel Chromium Cast Iron
107c	Nickel-Chromium-Molybdenum Cast Iron
342a	Nodular Cast Iron
338	White Cast Iron, Carbon and Sulfur





High Temperature Alloys (chip and disk)

SRM	Description	Unit Size
866	Incoloy™ 800	100 g
867	Incoloy™ 825	100 g
1230	High Temperature Alloy A286	disk: 32 mm diameter × 19 mm
1246	Incoloy™ 800	disk: 35 mm diameter × 19 mm
1247	Incoloy™ 825	disk: 35 mm diameter × 19 mm
1250	High Temperature Alloy (Fe - Ni - Co)	disk: 32 mm diameter × 19 mm
C2400	High Alloy Steel, ACI 17/4 PH	disk: 32 mm diameter × 19 mm
C2401	High Alloy Steel ACI-CD-4M Cu	disk: 32 mm diameter × 19 mm

Gases in Metals: Iron and Steel (rod)

These SRMs are certified for oxygen content. Materials certified for nitrogen are noted.

SRM	Description	Rod Size (mm)
1089*	Gasometric Standard, set includes: SRM 1095 AISI 4340 Steel SRM 1096 AISI 94B17 Steel, Modified** SRM 1097 Cr-V Steel, Modified SRM 1098 High Carbon Steel** SRM 1099 Electrolytic Iron	6.4 × 102 6.4 × 102 6.4 × 102 6.4 × 102 6.4 × 102
1754	AISI 4320 Oxygen in Low Alloy Steel,**	9.5 × 9.5 × 102
1090	Oxygen in Ingot Iron	6.35 × 102
1094	Oxygen in Maraging Steel	0.6 × 82
1091a	AISI 431 Oxygen in Stainless Steel	7.9 × 102
1093	Oxygen in Valve Steel	0.6 × 82

* These SRMs are sold only as a set designated SRM 1089.

** In addition to being certified for oxygen, these SRMs are also certified for nitrogen.

NONFERROUS METALS

Aluminum Base Alloys (chip and disk)

SRMs 1710 through 1715 are specially prepared to include low levels of cadmium and lead encountered in the analysis of recycled aluminum.

SRM	Description	Unit Size
87a	Silicon - Aluminum Alloy	75 g
855a	Aluminum Casting Alloy 356	30 g
856a	Aluminum Casting Alloy 380, Fine Millings	30 g
858	Alloy 6011, Modified	35 g
1258	Alloy 6011, Modified	disk: 35 mm diameter × 19 mm
859	Alloy 7075	35 g
1259	Alloy 7075	disk: 35 mm diameter × 19 mm
1710	Alloy 3004	disk: 63 mm diameter × 19 mm
1711	Alloy 3004	disk: 63 mm diameter × 19 mm
1712	Alloy 3004	disk: 63 mm diameter × 19 mm
1713	Alloy 5182	disk: 63 mm diameter × 19 mm
1714	Alloy 5182	disk: 63 mm diameter × 19 mm
1715	Alloy 5182	disk: 63 mm diameter × 19 mm

Cobalt Base Alloys (chip and disk)

SRM	Description	Unit Size
862	High Temperature Alloy L605	chip: 100 g
1242	High Temperature Alloy L605	disk: 35 mm diameter × 19 mm
1775	Refractory Alloy MP-35-N	disk: 35 mm diameter × 19 mm
2175	Refractory Alloy MP-35-N	chip: 50 g



Copper "Benchmark" (chip and rod)

Unit Size: Chip: 50 g

Rod: 6.4 mm × 103 mm

SRM		Description
Chip	Rod	
395	495	Unalloyed Copper - Cu II
396	496	Unalloyed Copper - Cu III
	457	Unalloyed Copper - Cu IV (6.6 mm diameter × 103 mm)
398	498	Unalloyed Copper - Cu V
399	499	Unalloyed Copper - Cu VI
400	500	Unalloyed Copper - Cu VII
454		Unalloyed Copper - Cu XI (35 g)

Copper Base Alloys (chip and rod)

SRM	Description	Unit Size (g)
158a	Silicon, Bronze	150
Beryllium-Copper		
458	17510	50
459	17200	50
460	17300	50
Phosphor-Bronze		
871	CDA 521	100
872	CDA 544	100
Cupro-Nickel		
874	10 % CDA 706, High-Purity	100
875	10 % CDA 706, Doped	100
Nickel-Silver		
879	CDA 762	100
880	CDA 770	100
1034	Unalloyed Copper	rod: 6.35 mm diameter × 103 mm
1035	Leaded-Tin Bronze Alloy	50

(continued)

Copper Base Alloys (block and disk)


The 1100 series SRMs are wrought disks 32 mm diameter \times 19 mm. The C1100 series SRMs are chill cast blocks 32 mm square \times 19 mm. Both forms have nearly identical elemental compositions.

SRM		Description
<i>Disk</i>	<i>Block</i>	
1104		Free-Cutting Brass
1107		Naval Brass B
1108		Naval Brass C
1110		Red Brass B
1111		Red Brass C
1112	C1112	Gilding Metal A
1113	C1113	Gilding Metal B
1114	C1114	Gilding Metal C
1115	C1115	Commercial Bronze A
1116	C1116	Commercial Bronze B
1117	C1117	Commercial Bronze C
	C1122	Beryllium-Copper
1276a		CDA 715 Cupro-Nickel



Lead Base Alloys (disk and powder forms)

SRM		Description	Unit Size (g)	
<i>Powder</i>	<i>Disk</i>		<i>Powder</i>	<i>Disk</i>
1129		Solder 63Sn - 37Pb	200	
127b	1131	Solder 40Sn - 60Pb	150	32 mm diameter \times 19 mm
53e	1132	Lead Base Bearing Metal (84Pb - 10Sb - 6Sn)	150	32 mm diameter \times 19 mm



Lead Base Materials (disk)

Unit Size: 50 mm diameter × 16 mm

SRM/RM	Description
--------	-------------

C2415	Battery Lead
C2416	Bullet Lead
C2417	Lead Base Alloy
C2418	High Purity Lead
8107	Gunpowder

Nickel Oxides (powder)

Unit Size: 25 g

SRM	Description
-----	-------------

671	Nickel Oxide No. 1
672	Nickel Oxide No. 2
673	Nickel Oxide No. 3



Nickel Base Alloys (chip and disk)

SRM	Description	Unit Size
349a	Waspaloy™	150 g
864	Inconel™ 600	100 g
865	Inconel™ 625	100 g
882	Nickel-Copper Alloy (65Ni - 31Cu - 3Al)	100 g
1159	Electronic and Magnetic Alloy Ni-Fe	disk: 31 mm diameter × 19 mm
1160	Electronic and Magnetic Alloy Ni-Mo	disk: 31 mm diameter × 19 mm
1243	Waspaloy™	disk: 34 mm diameter × 19 mm
1244	Inconel™ 600	disk: 35 mm diameter × 19 mm
C1248	Nickel-Copper Alloy (66Ni - 30Cu)	disk: 32 mm diameter × 19 mm
1249	Inconel™ 718	disk: 41 mm diameter × 19 mm
C2402	Hastelloy™ C	disk: 32 mm diameter × 19 mm

Trace Elements in Nickel Base Superalloys (chip)

Unit Size: 35 g

SRM	Description	Elemental Composition
897	"Tracealloy" A	Pb, Sc, Te, Ti
898	"Tracealloy" B	
899	"Tracealloy" C	

Tin Base Alloys (chip)

SRM	Description	Unit Size
54d	Tin Base Bearing Metal	170 g
1727	Anode Tin	30 × 30 × 30 mm

Titanium Base Alloys (chip and disk)

SRM	Description	Unit Size (g)
641	8 Mn (A)	disk: 32 mm diameter × 19 mm
642	8 Mn (B)	disk: 32 mm diameter × 19 mm
643	8 Mn (C)	disk: 32 mm diameter × 19 mm
647	6Al - 2Mo - 2Sn - 4Zr	50
648	5Al - 2Sn - 2Cr - 4Mo	50
649	15V - 3Al - 2Cr - 3Sn	50
650	Unalloyed Titanium A	30
651	Unalloyed Titanium B	30
654b	6Al - 4V	disk: 31 mm diameter × 19 mm
1128	15V - 3Al - 3Cr - 3Sn	disk: 35 mm diameter × 19 mm
2426	Galvalume	40
2431	6Al - 2Sn - 4Zr - 6Mo	50
2432	10V - 2Fe - 3Al	50
2433	8Al - 1Mo - 1V	50

**Hydrogen in Titanium (platelet)**

SRM	Description	Unit Size
352c	Hydrogen in Unalloyed Titanium	20 g
2452	Hydrogen in Titanium Alloys	1 × 10 g
2453	Hydrogen in Titanium Alloys	1 × 5 g
2454	Hydrogen in Titanium Alloys	1 × 10 g

Zirconium Base Alloys (chip)

SRM	Description	Unit Size
360b	Zircaloy-4	100 g

Zinc Base Alloys (chip and disk)

SRM	Description	Unit Size
94c	Die Casting Alloy	chip: 150 g
625	ASTM AG 40A Die Casting Alloy	disk: 44 mm diameter × 19 mm
626	ASTM AG 40A Die Casting Alloy	disk: 44 mm diameter × 19 mm
627	ASTM AG 40A Die Casting Alloy	disk: 44 mm diameter × 19 mm
628	ASTM AC 41A Die Casting Alloy	disk: 44 mm diameter × 19 mm
629	ASTM AC 41A Die Casting Alloy	disk: 44 mm diameter × 19 mm
630	ASTM AC 41A Die Casting Alloy	disk: 44 mm diameter × 19 mm
631	Zinc spelter, Modified	disk: 45 mm diameter × 19 mm
1736	Zinc-Aluminum (.31 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1737	Zinc-Aluminum (.63 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1738	Zinc-Aluminum (.10 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1739	Zinc-Aluminum (.21 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1740	Zinc-Aluminum (.42 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1741	Zinc-Aluminum (.52 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
1742	Zinc-Aluminum (.79 % Al) Die Casting Alloy	disk: 50.8 mm diameter × 12.7 mm
2139	Zinc-Aluminum (.80 % Al) Die Casting Alloy	chip: 100 g

CERAMICS AND GLASSES

Carbides (powder)


SRM	Description	Unit Size (g)
112b	Silicon Carbide	80
276b	Tungsten Carbide	75



Cemented Tungsten Carbides (powder)

Unit Size: 100 g

SRM	Description
887	Cemented Carbide (83W - 10Co)
888	Cemented Carbide (64W - 25Co - 5Ta)
889	Cemented Carbide (75W - 9Co - 5Ta - 4Ti)



Glasses (powder and solid)

SRM	Description	Unit Size (g)
81a	Glass Sand	75
89	Lead-Barium	45
92	Low-Boron Soda-Lime Powder	45
93a	High-Boron Boro-silicate	wafer: 32 mm diameter × 6 mm
165a	Glass Sand (low Iron)	75
620	Soda-Lime, Flat	3 platelets: 35 mm × 35 mm × 3 mm
621	Soda-Lime, Container	3 disks: 38 mm diameter × 5 mm
1411	Soft Borosilicate	10 platelets: 32 mm × 32 mm × 3 mm
1412	Multicomponent	8 platelets: 32 mm × 32 mm × 3 mm
1413	Glass Sand (high alumina)	75
1830	Soda-Lime, Float	3 platelets: 32 mm × 32 mm × 6 mm
1831	Soda-Lime, Sheet	3 platelets: 37 mm × 37 mm × 3 mm
1834	Fused Ore Glass	disk: 30 mm diameter × 3 mm

Trace Elements (powder and wafer)

These SRMs are for calibrating instruments and evaluating analytical techniques used to determine trace elements in inorganic matrices. SRMs 610 through 617 come in units of 6 wafers with wafer thicknesses of 3 mm for even numbered SRMs and 1 mm for odd numbered SRMs.

Also certified for isotopic ratio: $^{87}\text{Sr}/^{86}\text{Sr} = 1.20039$

SRM	Description	Certified Elements
607	Trace Elements in Potassium Feldspar (5 g)	
Trace Elements in Glass		
610/611		33 elements
612/613		33 elements
614/615		33 elements
616/617		33 elements

CEMENTS

Portland Cements (powder)

SRM	Unit Size
<i>Calcium Aluminate Cement</i>	
1882a	4 × 5 g
1883a	4 × 5 g
<i>Portland Cement</i>	
1880a	4 × 5 g
1881a	4 × 5 g
1884a	4 × 5 g
1885a	4 × 5 g
1886a	4 × 5 g
1887a	4 × 5 g
1888a	4 × 5 g
1889a	4 × 5 g

Portland Cement Clinkers (solid)

SRM	Unit Size
<i>Portland Cement Clinkers (5 phases certified)</i>	
2686	3 × 10 g
2687	3 × 10 g
2688	3 × 10 g





LUBRICANTS

Metallo-Organic Compounds

Unit Size: 5 g

These SRMs are for preparing solutions in oils of known and reproducible concentrations of metals.



SRM	Description	Elemental Composition
1075a	Aluminum 2-Ethylhexanoate	8.07 Al
1051b	Barium Cyclohexanebutyrate	28.7 Ba
1080a	Bis (1-phenyl-1,3-butanediono)copper (II)	16.37 Cu
1052b	Bis(1-phenyl-1,3-butanediono)oxovanadium (IV)	13.01 V
1053a	Cadmium Cyclohexanebutyrate	24.8 Cd
1057b	Dibutyltin bis (2-ethylhexanoate) (tin)	22.95 Sn
1059c	Lead Cyclohexanebutyrate	37.5 Pb
1060a	Lithium Cyclohexanebutyrate	4.1 Li
1065b	Nickel Cyclohexanebutyrate	13.89 Ni
1066a	Octaphenylcyclotetrasiloxane	14.14 Si
1077a	Silver 2-Ethylhexanoate	42.60 Ag
1069b	Sodium Cyclohexanebutyrate	12.0 Na
1070a	Strontium Cyclohexanebutyrate	20.7 Sr
1071b	Triphenyl Phosphate	9.48 P
1078b	Tris (1-phenyl-1,3-butanediono)chromium (III)	9.6 Cr
1079b	Tris (1-phenyl-1,3-butanediono)iron (III)	10.45 Fe
1073b	Zinc Cyclohexanebutyrate	16.66 Zn

PHYSICAL PROPERTIES

- 71 Ion Activity
- 74 Polymeric Properties
- 76 Thermodynamic Properties
- 80 Optical Properties
- 83 Electrical Properties
- 84 Optoelectronics
- 84 Metrology
- 87 Ceramics and Glasses
- 89 X-ray Spectrometry



ION ACTIVITY

pH Calibration

SRM	Description	pH(S) Values (at 25 °C)	Unit Size (g)
2193	Calcium Carbonate (In Prep)	—	—
723d	Tris(Hydroxymethyl) aminomethane (In Prep)	—	—
185h	Potassium Hydrogen Phthalate	4.006	60
188	Potassium Hydrogen Tartrate	3.557	60
189b	Potassium Tetroxalate	1.719	65
187d	Sodium Tetraborate Decahydrate (Borax)	9.182	30
Admixtures			
<i>Unit Size: 30 g (unless otherwise noted)</i>			
186Ig	Potassium Dihydrogen Phosphate	6.860*	
186IIg	Disodium Hydrogen Phosphate	7.414**	
191c	Sodium Bicarbonate (25 g)	10.015*	
192c	Sodium Carbonate		

*This pH results only when the two SRMs listed are used as an admixture in solution.

** Physiological buffer preparation.

Biological Buffer Systems

Unit Size: 60 g

SRM	Description	pH(S) Values (at 37 °C)	
		0.05 molal	0.08 molal
2181	HEPES Free Acid	7.364*	7.373*
2182	NaHEPESate		
2183	MOPSO Free Acid	6.699*	6.694*
2184	NaMOPSOate		



**This pH results only when the two SRMs listed are used as an admixture in solution.*

pD Calibration

SRM	Description	pD(S) Values (at 25°C)	Unit Size (g)
2185	Potassium Hydrogen Phthalate	4.518	60
2186I	Potassium Dihydrogen Phosphate	7.428*	30
2186II	Disodium Hydrogen Phosphate		30
2191a	Sodium Bicarbonate	10.732*	30
2192a	Sodium Carbonate		30

**This pD results only when the two SRMs listed are used as an admixture in solution.*

Ion-Selective Electrode Calibration

SRM	Description	Certified Property	Unit Size (g)
2201	Sodium Chloride	pNa, pCl	125
2202	Potassium Chloride	pK, pCl	160
2203	Potassium Fluoride	pF	125

Electrolytic Conductivity

SRM	Description	Nominal Conductivity ($\mu\text{S}/\text{cm}$)
3190	HCl in Deionized Water (In Prep)	—
<i>KCl in Deionized Water</i>		
3191		100
3192		500
3193		1000
3194		10 000
3195		100 000
<i>KCl in n-Propanol/Deionized Water</i>		
3198		5
3199		15
<i>NaCl in deionized Water</i>		
3196		—

Positive Electrophoretic Mobility


SRM	Description	Certified Property	Unit Size
1980	Goethite ($\alpha\text{-FeOOH}$)	$+\mu\text{E}$, $2.53 \mu\text{m} \cdot \text{cm}/\text{V} \cdot \text{s}$	40 mL

POLYMERIC PROPERTIES

Molar Mass/Molecular Weight (M_w)

SRM	M_w (g/mol)	Unit Size (g)
Poly(ethylene oxide)		
1924	$M_w \approx 120\,900$ ($M_w/M_n \approx 1.04$)	0.2
1923	$M_w \approx 26\,900$ ($M_w/M_n \approx 1.06$)	0.2
Poly(methylmethacrylate)		
1489*	$M_n \approx 115\,000$ ($M_w/M_n \leq 1.1$)	1.1
1488*	$M_n \approx 29\,300$ ($M_w/M_n \leq 1.1$)	2
1487*	$M_w \approx 6300$	2
Polyethylene/Polystyrene		
2887*	$M_w \approx 196\,400$	0.3
2885*	$M_w \approx 6280$	0.3
2886*	$M_w \approx 87\,000$	0.3
2888	$M_w \approx 7190$	0.3
Polyethylene, linear		
1475a*	$M_w \approx 52\,000$ ($M_w/M_n \approx 2.90$) (see also melt flow)	50
1484a*	$M_w \approx 119\,600$ ($M_w/M_n \approx 1.19$)	0.3
1482a*	$M_w \approx 13\,600$ ($M_w/M_n \approx 1.19$)	0.4
1483*	$M_w \approx 32\,100$ ($M_w/M_n \approx 1.11$)	1
Polystyrene, linear, broad molecular weight distribution		
706a	$M_w \approx 285\,000$	18
Polystyrene, linear, narrow molecular weight distribution		
1478*	$M_w \approx 37\,400$ ($M_w/M_n \approx 1.04$)	2
705a*	$M_w \approx 179\,300$ ($M_w/M_n \approx 1.07$)	5
1479	$M_w \approx 1\,050\,000$	2
Polyurethane		
1480	$M_w \approx 47\,300$	1

* Also certified for viscosity



Melt Flow Rate

SRM	Description	Melt Flow Rate (g/10 min)	Unit Size (g)
1473b	Polyethylene Resin, Low Density	1.13	50
1475a	Polyethylene, Linear	2.02	50
1474	Polyethylene Resin	5.03	60
1497	Polyethylene Gas Pipe Resin, Pigmented	0.186	9080
1496	Polyethylene Gas Pipe Resin, Unpigmented	0.26	908



Viscosity

SRM	Description	Unit Size (mL)
2490	Non-Newtonian Polymer Solution for Rheology (Polyisobutylene Dissolved in 2,6,10,14-Tetramethylpentadecane)	100
2491	Non-Newtonian Polymer Melt for Rheology	100

Biomaterials

RM	Description	Unit Size
8456	Ultra High Molecular Weight Polythylene <i>Properties:</i> - Young's Modulus - Yield Strength - Ultimate Strength - Elongation	bar: 7.62 cm diameter × 152.4 cm (3 in diameter × 60 in)
8457	Ultra High Molecular Weight Polythylene <i>Properties:</i> - Young's Modulus - Yield Strength - Ultimate Strength - Elongation	10 (0.5 cm) cubes

THERMODYNAMIC PROPERTIES

Calorimetry - Combustion

SRM	Description	Heat of Combustion (MJ/kg)*	Unit Size (g)
39j	Benzoic Acid	26.434	30
2692b	Coal, Bituminous: % S = 1.170	(32.81)**	50
2685b	Coal, Bituminous: % S = 4.730	(26.94)**	50
2682b	Coal, Sub-Bituminous: % S = 0.4917	(25.66)**	50
2151	Nicotinic Acid	22.184	25
2684b	Coal, Bituminous, Sulfur and Mercury: % S = 3.08; Hg = 97.4 µg/kg	28.56**	50
1657	Synthetic Refuse-Derived Fuel	13.87**	100
2683b	Sulfur and Mercury in Coal: % S = 1.955, Hg = 90.0 µg/kg	30.62	50
1656	Thianthrene	33.480	30
2152	Urea	10.536	25

* The calorific values (MJ/kg) may decrease upon the aging or normal oxidation of the coals. NIST will continue to monitor these calorific values and report any substantive change to the purchaser.

** Gross calorific value or HHV (Higher Heating Value).

Calorimetry - Solution

SRM	Description	Heat of Solution	Unit Size
1655	Potassium Chloride (Water Solution Calorimetry)	Absorbed (235.86 J/g)	30 g

Enthalpy and Heat Capacity

SRM	Description	Unit Size	Temperature Range (K)
RM 5	Copper	1.9 cm diameter 12 cm	25 to 300
781D2	Molybdenum	0.64 cm diameter 10 cm	273.15 to 2800
705a	Polystyrene (Molecular Weight: 170 900 g/mol)	5 g	10 to 350
720	Synthetic Sapphire	15 g	10 to 2250

Differential Scanning Calorimetry

SRM	Description	Melting Temperature (K)	Enthalpy of Fusion (J/g)	Unit Size
2222	Biphenyl (99.984 %)	342.41	120.41	1 g
2232	Indium (99.9999 %)	156.5985 °C	28.51	1 g
2225	Mercury	234.30	11.469	2.5 g
2220	Tin (99.9995 %)	505.10	60.2	(2.5 × 2.5 × 0.0127) cm
1514	Thermal Analysis Purity Set	4 levels of p-ABA (0.0 mol % to 5.0 mol %)		4 × 0.5 g

Differential Thermal Analysis

RM	Description	Temperature Range (°C)	Unit Size
GM 754	ICTA Polystyrene DTA	97.8 to 107.5	10 g
8759	ICTA Set DTA	295 to 675	5 × 10 g
8760	ICTA Set DTA	570 to 940	5 × 10 g




Defining Fixed Points, International Temperature Scale of 1990, ITS-90

SRM	Description	Temperature (°C)	Unit Size (g)
Pure Metals			
743	Mercury (Triple Point)	-38.8344	ampoule: 680
1745	Indium (Freezing Point)	156.5985	ingot: 20 × 10 g
741a	Tin (Freezing Point)	231.928	shot: 200
740a	Zinc (Freezing Point)	419.527	shot: 200
1744	Aluminum (Freezing Point)	660.323	ingot: 200
1746	Silver (Freezing Point)	961.780	shot: 300
Devices (semi-open cell)			
1747	Tin (Freezing Point), 99.9999+ %	231.928	1071
1748	Zinc (Freezing Point), 99.9999+ %	419.527	1031

Reference Points

SRM	Description	Temperature (°C)	Unit Size (g)
742	Alumina, 99.9+ % (Melting Point)	2052	powder: 10
45d	Copper (Freezing Point)	1084.6	bar: 450
49e	Lead (Freezing Point)	327.453	bar: 600



Freezing Point, Melting Point, and Triple Point Cells (sealed cell)

SRM	Description	Temperature (°C)	Unit Size (g)
1751	Gallium Melting Point	—	200
1968	Gallium (Melting Point), 99.9999+ %	29.7646	25
1972	1,3-Dioxolan-2-one (Ethylene Carbonate) (Triple Point), 99.999+ %	36.3143	60
1969	Rubidium (Triple Point), 99.9+ %	39.30	154
1973	n-Docosane (Triple Point), 99.999+ %	43.879	60
1970	Succinonitrile (Triple Point), 99.999+ %	58.0642	60
1971	Indium (Freezing Point), 99.9999+ %	156.598	100

Thermal Expansion of Metal and Glass

SRM	Description	Temperature Range (K)	Unit Size (cm)
731L1	Borosilicate Glass	80 to 680	0.64 × 5.1
731L2	Borosilicate Glass	80 to 680	0.64 × 10.2
731L3	Borosilicate Glass	80 to 680	0.64 × 15.2
736L1	Copper	20 to 800	0.64 × 5.1
738	AISI 446 Stainless Steel	293 to 780	0.64 × 5.1

Thermal Resistance of Glass, Silica, and Polystyrene

SRM	Description	Temperature Range (K)	Thermal Resistance (m ² · K · W ⁻¹)	Unit Size (cm)
1453	Expanded Polystyrene Board	285 to 310	0.381 to 0.420	66 × 93 × 1.34
1450c	Fibrous Glass Board	280 to 340	0.661 to 0.818	61 × 61 × 2.54
1449	Fumed Silica Board	297	1.195 to 1.253	60 × 60 × 2.54
1459	Fumed Silica Board	297	1.195 to 1.253	30 × 30 × 2.54

Vapor Pressure of Metals

SRM	Description	Pressure Range (Pa) (K, ITS-90)	Temperature Range	Unit Size
745	Gold	10^{-3} to 10^2	1300 to 2100	wire: 0.14 cm diameter \times 15.2 cm
746	Cadmium	10^{-6} to 10^1	350 to 594	rod: 0.64 cm diameter \times 6.4 cm

Thermal Conductivity of Graphite and Iron

RM	Conductivity Range ($\text{W}\cdot\text{m}^{-1}\cdot\text{K}^{-1}$)	Unit Size
<i>Electrolytic Iron (2 K to 1000 K)</i>		
8420	12.32 to 32.98	0.64 cm diameter \times 5.0 cm
8421	12.32 to 32.98	3.17 cm diameter \times 5.0 cm
<i>Graphite (5 K to 2500 K)</i>		
8424	0.0354 to 32.96	0.64 cm diameter \times 5.0 cm
8426	0.0354 to 32.96	2.54 cm diameter \times 5.0 cm

Laboratory Thermometer (mercury in glass)

Unit Size: 1 each

SRM	Description	Calibrated Points ($^{\circ}\text{C}$)
934	Clinical Laboratory Thermometer	-0.20 to +0.20

Thermocouple Material, Platinum

Unit Size: 1 each

SRM	Description	Temperature Range
1749	Gold vs. Platinum Thermocouple Thermometer	0 $^{\circ}\text{C}$ to 1000 $^{\circ}\text{C}$
1967	Platinum Wire, High Purity (99.999+ %)	-197 $^{\circ}\text{C}$ to 1768 $^{\circ}\text{C}$
1750	Standard Platinum Resistance Thermometer	14 K to 430 K

OPTICAL PROPERTIES

Molecular Transmittance and Absorbance



SRM	Description	Wavelength Range	Unit Size
Crystalline and Solution Forms			
935a	Crystalline Potassium Dichromate, UV Absorbance	235 nm to 350 nm	15 g
1935	Potassium Dichromate Solution, UV Absorbance	235 nm to 350 nm	10 ampoules: 5 samples, plus 5 blanks
2032	Potassium Iodide, Stray Light	240 nm to 275 nm	25 g
931f	Liquid Filters, Absorbance	302 nm to 678 nm	12 ampoules: 3 × 3 levels, plus 3 blanks
Glass Filters, Transmittance			
930e	10 %, 20 %, 30 % Transmittance	440 nm to 635 nm	3 filters, plus 1 blank
1930	1 %, 3 %, 50 % Transmittance	440 nm to 635 nm	3 filters, plus 1 blank
2030a	30 % Transmittance	465.0 nm	1 filter, plus 1 blank
2031b	Metal-on-Quartz Filters 10 %, 30 %, 90 % Transmittance	250 nm to 635 nm	3 filters, plus 1 blank
2046	Optical Density = 1	1064 nm	51 mm × 51 mm × 1.0 mm
2047	Optical Density = 2	1064 nm	51 mm × 51 mm × 2.2 mm
2048	Optical Density = 3	1064 nm	51 mm × 51 mm × 3.2 mm
2049	Optical Density = 4	1064 nm	51 mm × 51 mm × 4.2 mm
2050	Optical Density = 5	1064 nm	51 mm × 51 mm × 5.4 mm
2051	Optical Density = 6	1064 nm	51 mm × 51 mm × 6.4 mm
2053	20 nm Ni-Cr Film on Silica	2 μm to 25 μm	25 mm diameter × 250 μm
2054	90 nm Ni-Cr Film on Silica	2 μm to 25 μm	25 mm diameter × 250 μm
2055	77 nm Cu-Ni Film on Silica	2 μm to 25 μm	25 mm diameter × 250 μm
2056	97 nm Cu-Ni Film on Silica	2 μm to 20 μm	25 mm diameter × 250 μm
2930	Ultimate Range Visible Absorbance Filters	—	3 filters & 1 blank



Transmittance Wavelength Standards

SRM	Description	Wavelength Range	Unit Size
2034	Holmium Oxide Solution	240 nm to 650 nm	1 sealed cuvette
2035	Near-IR Transmission	971 nm to 1949 nm	25 mm diameter × 1.5 mm
2036	Near-IR Wavelength/Wavenumber Reflection Standard	975 nm to 1946 nm	—
2037	Red Diesel Dye	—	100 mg
2065	Transmission Wavelength/Vacuum Wavenumber	ultraviolet–visible–near-infrared	25 mm diameter × 1.5 mm
1921a	Infrared Transmission	3.2 μ m to 18.5 μ m	1 polystyrene film

Fluorescence

SRM/RM	Description	Wavelength Range	Unit Size
936a	Quinine Sulfate Dihydrate	375 nm to to 675 nm	1 g
1932	Fluorescein	488 nm to 191 nm	3 × 2 mL
8640	Fluorescein Labeled Microbead Suspension	—	—
2242	Relative Intensity Correction, Raman Spectroscopy	—	1 artifact
2241	Relative Intensity Correction, Raman Spectroscopy	785 nm	1 glass slide (10.7 × 30.4 × 2.0mm)



Specular Spectral Reflectance

SRM	Description	Wavelength Range	Unit Size
2003	First Surface, Aluminum on Glass	250 nm to 2500 nm	5.1 cm diameter × 0.65 cm
2026	Second Surface, Aluminum on Fused Quartz	250 nm to 2500 nm	5.1 cm diameter × 0.6 cm
2017	Multi-Angle White Reflectance Standard	360 nm to 780 nm	5.7 cm diameter × 1.3 cm
2040	PTFE Diffuser for Spectral Reflectance Factor	380 nm to 780 nm	5 × 26 g

Near Infrared Reflectance Wavelength Standard

SRM	Description	Wavelength Range	Unit Size
1920a	Rare Earth Oxide Mixture	740 nm to 2000 nm	5.1 cm diameter × 1.2 cm

Optical Rotation

SRM	Description	Wavelength Range	Unit Size
917b	D-Glucose (Dextrose)	546 nm to 589 nm	50 g
17e	Sucrose	546 nm to 633 nm	60 g

Liquid Refractive Index

SRM	Description	Wavelength Range	Unit Size
1922	Mineral Oil	468 nm to 589 nm	30 mL

X-ray and Photographic Imaging

SRM	Description	Unit Size
1010a	Microcopy Resolution Test Chart	5 charts
1008	Photographic Step Tablet	25.4 cm × 3.5 cm
1001	X-ray Film Step Tablet	25.4 cm × 3.5 cm



ELECTRICAL PROPERTIES

Electrical Resistivity and Conductivity of Electrolytic Iron and Graphite

Unit Size: rod: 0.64 cm diameter \times 5.0 cm

RM	Resistivity Range ($\mu\Omega \cdot m$)	Unit Size
Electrolytic Iron (2 K to 1000 K)		
8420	0.004 to 0.909	0.64 cm diameter \times 5.0 cm
8421	0.004 to 0.909	3.17 cm diameter \times 5.0 cm
Graphite (5 K to 2500 K)		
8424	28.78 to 12.59	0.64 cm diameter \times 5.0 cm
8426	28.78 to 12.59	2.54 cm diameter \times 5.0 cm



Electrical Resistivity and Conductivity of Silicon

SRM	Resistivity ($\Omega \cdot cm$)	Type
Spreading Resistance		
Unit Size: 16 chips: (5 \times 10 \times 0.625) mm		
2526	0.001 to 200	(111) p-type
2527	0.001 to 200	(111) n-type
Silicon Resistivity		
Unit Size: (100 diameter \times 0.625) mm		
2544	10	float zone n-type silicon wafer with (111) crystallographic orientation
2547	200	

OPTOELECTRONICS

SRM	Description	Unit Size
Wavelength Calibration Standards		
2514	Wavelength Calibration Reference for 1560 nm to 1595 nm - Carbon Monoxide ($^{12}\text{C}^{16}\text{O}$)	Gas Absorption Cell
2515	Wavelength Calibration Reference for 1595 nm to 1630 nm - Carbon Monoxide ($^{13}\text{C}^{16}\text{O}$)	Gas Absorption Cell
2517a	High Resolution Wavelength Calibration Reference for 1510 nm to 1540 nm - Acetylene ($^{12}\text{C}_2\text{H}_2$)	Gas Absorption Cell
2519	Wavelength Reference Absorption Cell for 1530 nm to 1560 nm Hydrogen Cyanide ($\text{H}^{13}\text{C}^{14}\text{N}$)	Gas Absorption Cell
Polarization Mode Dispersion Standards		
2518	Polarization Mode Dispersion Standard	1 each
2538	Deterministic Polarization Mode Dispersion Standard	1 each
Fiber and Fiber-Connector Geometry Standards		
2513	Mode Field Diameter Standard for Single-Mode Fiber	1 each
2520	Optical Fiber Diameter Standard	1 each
2522	Pin Gauge Standard for Optical Fiber Ferrules	1 wire-sizing bore
2523	Optical Fiber Ferrule Geometry Standard	1 ceramic connector ferrule
2553	Optical Fiber Coating Diameter ($n = 1.504$)	1 each: 250 μm diameter
2554	Optical Fiber Coating Diameter ($n = 1.515$)	1 each: 250 μm diameter

METROLOGY

Optical Microscope Linewidth Measurement



SRM	Linewidth (μm)	Pitch (μm)	Unit Size (cm)
Linewidth Measurement Standards			
475	0.9 to 10.8	2 to 36	$6.35 \times 6.35 \times 0.15$
476	0.9 to 10.8	2 to 36	$6.35 \times 6.35 \times 0.15$
2800*			$25 \times 75 \times 2.3$


* SRM 2800 is used in calibrating magnification and consists of a pattern of parallel lines whose nominal distances from the centerline range from $\pm 1 \mu\text{m}$ to $\pm 5 \text{ mm}$. Certified values are given for the center-to-center distance of each line from the centerline; the linewidths are not certified.

Scanning Electron Microscope (SEM)

SRM/RM	Description	Spacings	Unit Size (mm)
2069b	SEM Performance Standard	2 mm to 4 mm	12 mm diameter with 3 mm peg
8091	SEM Sharpness Standard		semiconductor chip: 2 mm × 2 mm
2800	Microscope Magnification Standard	1 μ m to 5 mm	25 × 75 × 2.3

Depth Profiling

SRM	Description	Value	Unit Size (cm)
2133	Phosphorus Implant in Silicon Depth Profile Standard	^{31}P : 0.04927 $\mu\text{g}/\text{cm}^2$ (9.58×10^{14} atoms/ cm^2)	crystal 1 × 1
2134	Arsenic Implant in Silicon Profile Standard	^{75}As - 7×10^{14} atoms/ cm^2	crystal: 1 × 1
2135c	Nickel-Chromium Thin-Film Depth Profile Standard	Cr: 41.3 $\mu\text{g}/\text{cm}^2$ Ni: 49.4 $\mu\text{g}/\text{cm}^2$	1 × 2.54 × 0.04
2137	Boron Implant in Silicon Depth Profile Standard	^{10}B - 1.018 v 1015 atoms/ cm^2	1 × 1



SILICON CRYSTAL

Solder Thickness for X-ray Fluorescence

Unit Size: plate: 15 mm × 15 mm


SRM	Description	Composition	Coating Mass/Area	Coating Thickness	
				(μm)	(μm)
2321	Tin-Lead Alloy	60 % Sn, 40 % Pb	6.8 mg/ cm^2	295	7.5

Coating Thickness

Unit Size: 45 mm × 45 mm

These SRMs are suitable for calibrating instruments based on magnetic induction and magnetic pull-off techniques used in the measurement of organic and non-magnetic inorganic coatings over steel.

SRM	Nominal Coating Thickness	
	(μm)	(mils)
Chromium over Copper on Steel		
1358a	80, 255, 1000	3.1, 9.8, 39
1359b	48, 140, 505, 800	2.0, 5.5, 20, 32
1361b	6, 12, 25, 48	0.2, 0.5, 1.0, 2.0
1362b	40, 80, 140, 205	1.6, 3.1, 5.5, 7.9
1363b	255, 385, 505, 635	9.8, 16, 20, 26
1364b	800, 1000, 1525, 1935	32, 39, 59, 79



Ellipsometry

Unit Size: 76 mm substrate diameter

Each unit is certified for the ellipsometric parameters delta (Δ) and psi (ψ) at the vacuum wavelength $\lambda = 633.0$ nm, and for the derived values of the thicknesses and indexes of refraction of the silicon dioxide and silicon layers.

SRM	Thickness (nm)
Thin Film Thickness Standards	
2531	50
2532	100
2533	200
2534	25
2535	14

Oxygen Concentration in Silicon

SRM	Description	Unit Size (mm)	Concentration (mg/kg)
2551	Oxygen in Silicon	4 wafers: 25 × 25 × 2	Low: 10 Medium: 13 High: 15 FZ: (<0.1)

Superconducting Critical Current (wire form)

Unit Size: wire: 8.7 cm diameter × 2.2 m

SRM	Description	Magnetic Field Range (T)	Critical Current Range (A)
1457	Niobium-Titanium Wire	2.000 to 8.000	293.30 to 69.72

CERAMICS AND GLASSES

Chemical Resistance [Durability] of Glass



SRM	Description	mL of N/50 H ₂ SO ₄	Unit Size (kg)
623	Borosilicate	0.34	2.2
622	Soda-Lime Silica	7.67	2.2

Electrical Properties of Glass

Unit Size: 5 cm × 5 cm × 2.5 cm

SRM 624 is suitable for use with ASTM C 657. SRM 774 is suitable for use with ASTM D 150.

SRM	Description	Unit Size (cm)	Value
624	Lead Silica for dc Volume Resistivity	5 × 5 × 2.5	$\log_{10} \rho \approx 9.9 \text{ } \Omega \cdot \text{cm at } 300 \text{ } ^\circ\text{C}$
774	Lead Silica for Dielectric Constant and ac Loss Characteristics	5 × 5 × 2.5	$K \approx 7.47 \text{ at } 100 \text{ Hz}$

Viscosity of Glass

SRM	Description	Unit Size (mm)
717a	Borosilicate Glass	block: 40 × 40 × 150
710a	Soda-Lime-Silica Glass	block: 100 × 100 × 40

Viscosity Fixpoints of Glass

These SRMs are for the calibration of equipment for the determination of the softening, annealing, and strain points of glass.

SRM	Description	Unit Size
714	Alkaline Earth Alumina Silicate	225 g
717a	Borosilicate	40 mm × 40 mm × 150 mm
713	Dense Barium Crown 620/603 Glass	225 g
709	Extra Dense Lead Silica	4 cm × 4 cm × 5 cm
716	Neutral Glass	250 g
710a	Soda-Lime-Silica	100 mm × 100 mm × 40 mm



Relative Stress Optical Coefficient

SRM	Description	Relative Stress Optical Coefficient (C) at $\lambda = 546.1$ nm (Value $\times 10^{-12}$ m ² /N)	Unit Size
709	Extra Dense Lead Silica	C = -1.359	bar: 4 cm × 4 cm × 5 cm

Density

SRM	Description	Density (kg/m ³)	Unit Size
1827b	Lead Silica Glass	3593.800 at 20 °C	slab: 25 cm × 25 cm × 12 cm
211d	Toluene	871.476 at 15 °C	4 × 5 mL
2214	Isooctane	695.969 at 15 °C	4 × 5 mL

Glass Liquidus Temperature

SRM	Description	Unit Size	Method	Temperature (°C)
773	Soda-Lime-Silica	2.5 cm × 2.5 cm × 0.6 cm	A (boat)	988
			B (perforated plate)	991
1416	Aluminosilicate	22 lengths of 12.7 cm tube (250 g)		1147

X-RAY SPECTROMETRY

X-ray Diffraction

SRM	Description	XRD Application	Unit Size (g)
676	Alumina (Corundum Structure)	Quantitative Analysis	20
1976	Alumina Plate, Sintered	Instrument Response	45 mm × 45 mm × 1.6 mm
2910	Calcium Hydroxyapatite	Quantitative Analysis	5
660a	Lanthanum Hexaboride Powder	Line Position, Line Shape	6
675	Mica	Low 2θ (Large d-Spacing)	7.5
1879a	Respirable Cristobalite	Quantitative Analysis	5
1878a	Respirable Quartz	Quantitative Analysis	5
656	Silicon Nitride	Quantitative Analysis	2 × 10 g
640c	Silicon Powder 2~d-Spacing	Line Position, Line Shape	7.5
674b	X-ray Powder Diffraction Intensity Set (α-Al ₂ O ₃ , CeO ₂ , Cr ₂ O ₃ , TiO ₂ , ZnO) (In Prep)	Quantitative Analysis	—

X-ray Stage Calibration

SRM	Description	Unit Size (mm)
1842	Calibration Board (X and Y dimensions)	Board: 300 × 300 × 3
1843	Calibration Board (Z dimension)	Triangular Block: 37 × 20 × 12

RADIOACTIVITY

- 91 Radioactive Solutions
- 93 Radioactive Point Sources
- 93 Radiopharmaceuticals
- 94 Beryllium Isotopic Ratio Standard
- 94 Carbon-14 Dating
- 95 Natural Matrix Materials
- 95 Neutron Density Monitor Wire
- 95 Fission Track Glass





Radioactive Solutions

Each SRM is contained in a 5 mL flame-sealed glass ampoule and, except as noted, consists of the radionuclide dissolved in an aqueous solution (usually acidic, unless noted).

SRM	Radionuclide	Massic Activity (Bq/g)	Calibration Date
4322B*	Americium-241	40	09/1991
4332D*	Americium-243	40	05/1995
4251C*	Barium-133	500 000	09/1993
4222C	Carbon-14 (as hexadecane)	50 000	09/1990
4943	Chlorine-36 (3 mL)	10 000	12/1984
4915E*	Cobalt-60	75 000	01/1995
4329*	Curium-243	70	06/1984
4320A*	Curium-244	35	02/1996
4370C*	Europium-152	90 000	02/1987
4361C	Hydrogen-3 (as water) (500 mL)	2	09/1998
4926E	Hydrogen-3 (as water) (20 mL)	5 000	09/1998
4927F	Hydrogen-3 (as water)	600 000	09/1998
4947C	Hydrogen-3 (as toluene) (4 mL)	300 000	03/1987
4949C*	Iodine-129 (dilute base solution)	3 000	03/1993
4341*	Neptunium-237	100	03/1992
4226C	Nickel-63	50 000	08/1995
4323B*	Plutonium-238	40	11/1999
4330B*	Plutonium-239	40	11/1999
4338A*	Plutonium-240	40	05/1996
4334G*	Plutonium-242	25	06/1994
4326	Polonium-209	90	03/1994

* These SRMs require a license certification.

(continued on next page)

Radioactive Solutions (continued)

SRM	Radionuclide	Massic Activity (Bq/g)	Calibration Date
4969	Radium-226	3	09/1998
4965	Radium-226	30	09/1991
4966	Radium-226	270	09/1991
4967	Radium-226	2 700	09/1991
4919H*	Strontium-90	4 000	07/1995
4234A*	Strontium-90	2 500 000	03/1995
4288A	Technetium-99	30 000	09/1996
4324B	Uranium-232	40	—
4321C	Uranium-238, Natural Uranium	250	08/1997



* These SRMs require a license certification.



Radioactive Point Sources

SRM	Radionuclide	Principal Photon Energies (keV)	Activity (Bq)	Calibration Date
4241C	Barium-133	81 to 384	60 000 to 170 000	01/1999
4218F*	Europium-152	122 to 1400	50 000 to 150 000	01/1999
4201B*	Niobium-94	702, 871	4000	04/1970

* These SRMs require a license certification.

Radiopharmaceuticals

A typical schedule of SRMs for a year:

SRM	Radionuclide	High Level†	Low Level	Month
4401*	Iodine-131	750 MBq (20 mCi)	25 MBq (700 µCi)	January
4412*	Molybdenum-99	1.5 GBq (40 mCi)	75 MBq (2 mCi)	February
4415*	Xenon-133	7.5 GBq (200 mCi)	750 MBq (20 mCi)	March
4416*	Gallium-67	375 MBq (10 mCi)	20 MBq (500 µCi)	April
	OPEN			May
4404*	Thallium-201	375 MBq (10 mCi)	35 MBq (900 µCi)	June
4425*	Samarium-153	375 MBq (10 mCi)	20 MBq (500 µCi)	July
4417*	Indium-111	375 MBq (10 mCi)	20 MBq (500 µCi)	August
4410*	Technetium-99m	7.5 GBq (200 mCi)**	**	September
4407*	Iodine-125	750 MBq (20 mCi)	6 MBq (150 µCi)	October
	OPEN			November
4427*	Yttrium-90	200 MBq (5 mCi)	20 MBq (500 µCi)	December

* These SRMs require a license certification.

† High level radiopharmaceutical SRMs are distributed through a program with the exception of technetium-99m. This program includes:

- Distribution of radioactivity Standard Reference Materials (usually a high level [10^7 Bq to 10^{10} Bq] and low level [10^6 Bq to 10^8 Bq] pair of sources), for quantification, to sponsoring participating companies ten times per year.
- Two "open" months per year when the sponsors have the opportunity to send in sources of their choice for calibration.

For more information, please contact Daniel Golas, the program's Project Manager, at NIST (301) 975-5540.

** Technetium-99m is also available to purchasers of low level radiopharmaceutical SRMs.

Beryllium Isotopic Ratio Standard

SRM	Radionuclide	Isotopic Ratio	Total Nuclide Concentration (mg/g)	Calibration Date	Volume (mL)
4325	Beryllium-10/ Beryllium-9 (in 1N HCl)	3×10^{-11}	5	08/1986	50

Carbon-14 Dating Standard

This SRM is an international standard for contemporary carbon-14 against which world-wide measurements can be compared.

SRM	Description	Unit Size
4990C	Oxalic Acid	8 × 28 g





Natural Matrix Materials

These SRMs are intended for use in testing measurements of low level, biological, and environmental radioactivity contained in matrices similar to the sample, for evaluating analytical methods, or as a generally available calibrated “real” sample matrix in interlaboratory comparisons.

SRM	Description	Unit Size (g)
4350B	Columbia River Sediment	85
4351	Human Lung	45
4352	Human Liver	45
4354	Freshwater Lake Sediment	25
4355	Peruvian Soil	75
4356	Ashed Bone	15
4357	Ocean Sediment	85



Neutron Density Monitor Wire

SRM	Description	Cobalt Composition (weight %)	Unit Size
953	Cobalt in Aluminum Wire	0.116	0.5 mm diameter × 1 m

Fission Track Glass

Each unit consists of four unirradiated glass wafers and two irradiated wafers.

SRM	Uranium Composition (µg/g)	Uranium-235 (Atom %)	Reactor Position	Neutron Fluence (× 10 ¹⁴ n/cm ²)	
				Copper Foil	Gold Foil
963a	0.823	0.2792	RT-4	39.5	43.0
			RT-3	41.2	45.8

INDUSTRIAL HYGIENE

- 97 Materials on Filter Media
- 97 Trace Constituent Elements
in Blank Filters
- 97 Respirable Silica
- 98 Lead in Paint, Dust,
and Soil
- 99 Asbestos





Materials on Filter Media

These SRMs consist of potentially hazardous materials deposited on filters to be used to determine the levels of these materials in industrial atmospheres.

SRM/RM	Description	Set Size	Elemental Composition	Diameter (mm)	Pore Size (µm)
2679a	Quartz on Filter Media	2 × 3 levels, plus 2 blanks	Quartz, Clay	47	0.45
2783	Air Particulate on Filter	2 filters, plus 2 blanks	18 certified values 9 reference values	47	0.4
RM 8785	Particulate Matter on Filters	3 filters	1 reference value 2 information values	37	—

Trace Constituent Elements in Blank Filters

SRMs 2678 and 2681 are for use in evaluating the performance of air sampling filter methods with either certified values (in µg) or limits of detection (X_D) for each of 30 constituent elements, as well as six leachable anions and cations.

SRM	Description	Diameter (mm)	Pore Size (µm)	Filter Weight (g)
2678	Cellulose Acetate Membrane	47	0.45	0.09
2681	Ashless Blank Filter	42.5	—	0.14

Respirable Silica

These SRMs are intended for use in determining, by X-ray diffraction, the levels of respirable silica in an industrial atmosphere according to the National Institute for Occupational Safety and Health (NIOSH) Analytical Method 7500 or equivalent methods.

SRM	Description	Mass Fraction/Mass Loading	Unit Size
1878a	Respirable Alpha Quartz	100.00% ± 0.21%	5 g
1879a	Respirable Cristobalite	95.6% ± 0.4%	5 g
2950	Respirable Alpha Quartz on Filter Media	(10, 20, 50, 100, 250, 500) µg/filter	set SRMs 2952-57
2951	Respirable Alpha Quartz on Filter Media	5 µg/filter	5 filters (5 blanks)
2952	Respirable Alpha Quartz on Filter Media	10 µg/filter	5 filters (5 blanks)
2953	Respirable Alpha Quartz on Filter Media	20 µg/filter	5 filters (5 blanks)
2954	Respirable Alpha Quartz on Filter Media	50 µg/filter	5 filters (5 blanks)
2955	Respirable Alpha Quartz on Filter Media	100 µg/filter	5 filters (5 blanks)
2956	Respirable Alpha Quartz on Filter Media	250 µg/filter	5 filters (5 blanks)
2957	Respirable Alpha Quartz on Filter Media	500 µg/filter	5 filters (5 blanks)
2958	Respirable Alpha Quartz on Filter Media	1000 µg/filter	5 filters (5 blanks)
2960	Respirable Alpha Cristobalite on Filter Media	(5, 10, 20, 50, 100, 250) µg/filter	set SRMs 2961-66
2961	Respirable Alpha Cristobalite on Filter Media	5 µg/filter	5 filters (5 blanks)

(continued)

Respirable Silica (continued)

SRM	Description	Mass Loading	Unit Size
2962	Respirable Alpha Cristobalite on Filter Media	10 µg/filter	5 filters (5 blanks)
2963	Respirable Alpha Cristobalite on Filter Media	20 µg/filter	5 filters (5 blanks)
2964	Respirable Alpha Cristobalite on Filter Media	50 µg/filter	5 filters (5 blanks)
2965	Respirable Alpha Cristobalite on Filter Media	100 µg/filter	5 filters (5 blanks)
2966	Respirable Alpha Cristobalite on Filter Media	250 µg/filter	5 filters (5 blanks)
2967	Respirable Alpha Cristobalite on Filter Media	500 µg/filter	5 filters (5 blanks)

Lead in Paint, Dust, and Soil

These SRMs and RM have been developed in conjunction with the U.S. EPA to monitor paint, dust, and soil sources of lead.

SRM	Lead Concentration	Unit Size
Paint Film		
2570	<0.001 mg/cm ²	1 blank film
2571	3.58 mg/cm ²	1 film, plus 1 blank
2572	1.527 mg/cm ²	1 film, plus 1 blank
2573	1.040 mg/cm ²	1 film, plus 1 blank
2574	0.714 mg/cm ²	1 film, plus 1 blank
2575	0.307 mg/cm ²	1 film, plus 1 blank
2579a (Set of 6: SRMs 2570 to 2575)	0.307 to 3.58 mg/cm ²	5 films, plus 1 blank
2576 (High Level)	5.59 mg/cm ²	1 film, plus 1 blank
Powdered Paint		
2580	4.34 %	30 g
2581	0.449 %	35 g
2582	209.8 mg/kg	20 g
2589	9.99 %	35 g
Indoor Dust, Trace Elements in (As, Cd, Cr, Hg, Pb)		
2583	85.9 mg/kg	8 g
2584	9761 mg/kg	8 g
Soil, Trace Elements in		
2586	432 mg/kg	50 g
2587	3242 mg/kg	50 g
Paint on Fiberboard		
RM 8680	1 to 2 mg/cm ²	1 sheet: (10.2 × 15.2 × 1.3) cm



Asbestos

SRM	Description	Asbestos Type	Unit Size
1866b	Common Commercial Asbestos	chrysotile grunerite (Amosite) riebeckite (Crocidolite)	3 × 4 g
1868	Quantitative Asbestos in Building Material	chrysotile grunerite	set (2) 5 – 10 g each
1876b	Chrysotile Asbestos for TEM	chrysotile	10 sections: 3 mm × 3 mm
RM 8411	Mixed Asbestos Research Filter	chrysotile asbestos grunerite (Amosite)	1 cm ²



ASBESTOS TESTING

SUBJECT INDEX

A

ABSORBANCE

- 80 See MOLECULAR...
SPECTROMETRY

ACETANILIDE

- 45 use in MICROCHEMISTRY

ACIDIMETRIC VALUE (STOICHIOMETRY)

- 44 of Benzoic Acid
- 44 of Boric Acid
- 71 of Potassium Hydrogen
Phthalate

ADHESION (TAPE ADHESION TESTING)

- 6 Linerboard for

AGRICULTURAL MATERIALS

- 11 Apple Leaves
- 9,11 Corn Kernel (Zea Mays)
- 9,11 Corn Stalk (Zea Mays)
- 11 Fluoride in Vegetation
- 11 Peach Leaves
- 10 Peanut Butter
- 11 Pine Needles
- 9,11 Spinach Leaves
- 11 Tomato Leaves
- 10 Slurried Spinach

AIR PARTICULATE

- 97 See MATERIALS
ON FILTER MEDIA

AIR POLLUTION

- 25 See PRIMARY GAS
MIXTURES

ALCOHOL

- 17 Ethanol Solutions

ALCOHOLS (FOSSIL FUELS)

- 30 Alcohol in Gasoline
- 30 Ethanol
- 30 Methanol
- 30 Methanol and t-Butanol

ALLOYS (FERROUS)

- 51 See FERROUS METALS

ALLOYS (NONFERROUS)

- 60 See NONFERROUS METALS

ALUMINA

- 33 as Bauxite (ORES)
- 35 as Burnt REFRACTORIES
- 78 REFERENCE POINT
- 2 for SURFACE AREA OF POW-
DERS
- 89 X-RAY SPECTROMETRY

ALUMINUM

- 78 Freezing Point of (DEFINING
FIXED POINT, ITS-90)
- 36 as a METALLO-ORGANIC
COMPOUND
- 46 SPECTROMETRY Solution
- 81 Specular Reflectance
(Mirrors)

ALUMINUM BASE ALLOYS

- 60 See NONFERROUS
METALS

AMERICIUM (RADIOACTIVITY)

- 91 Americium-241
- 91 Americium-243
- 95 Columbia River Sediment
- 95 Human Liver
- 95 Human Lung
- 95 Peruvian Soil

AMMONIUM DIHYDROGEN PHOSPHATE

- 11 See FERTILIZERS

ANALYZED GASES

- 25 See PRIMARY GAS
MIXTURES

ANGIOTENSIN I

- 13 See HEALTH & CLINICAL

ANISIC ACID

- 45 use in MICROCHEMISTRY

ANION CHROMATOGRAPHY

- 48 Bromide Solution
- 48 Chloride Solution
- 48 Fluoride Solution
- 48 Nitrate Solution
- 48 Phosphate Solution
- 48 Sulfate Solution

ANTICONVULSANT DRUG LEVEL ASSAY

- 13 See HEALTH & CLINICAL

ANTIEPILEPSY DRUG LEVEL ASSAY

- 13 See HEALTH & CLINICAL

ANTIMONY

- 46 SPECTROMETRY Solution

ARGILLACEOUS LIMESTONE

- 34 See ROCKS AND MINERALS

ARSENIC

- 85 Implant in Silicon (DEPTH
PROFILING)
- 46 SPECTROMETRY Solution

ARSENIC TRIOXIDE (STOICHIOMETRY)

- 44 Reductometric value

ASBESTOS

- 99 Common Commercial
- 99 Mixture on Filter

ASHED BONE (RADIOACTIVITY)

- 95 NATURAL MATRIX
MATERIALS

ATOMIC ABSORPTION SPECTROMETRY

- 46 See SPECTROMETRIC
SINGLE ELEMENTS

AUTO CATALYSTS

- 25 Recycled Monolith
- 25 Recycled Pellet

B

BALL BAR (PERFORMANCE MATERIALS)

- 6 Coordinate Measuring
Machine Probe

BARIUM

- 91 as Barium-133
(RADIOACTIVITY)
- 91 as Cesium-137 Burn-up
Standard
- 36 as a METALLO-ORGANIC
COMPOUND
- 46 SPECTROMETRY Solution

BASALT ROCK

- 34 See ROCKS AND MINERALS

BASIMETRIC VALUE (STOICHIOMETRY)

- 44 of Tris(hydroxymethyl)-aminomethane

BAUXITE (ORES)

- 33 from Arkansas
- 33 from the Dominican Republic
- 33 from Jamaica
- 33 from Surinam

BEARING METAL (PB-SB-SN)

- 62 See LEAD BASE ALLOYS

BENZOIC ACID

- 44 Acidimetric Value (STOICHIOMETRY)
- 76 Calorimetric Value (COMBUSTION CALORIMETRY)

BERYLLIUM

- 61 in COPPER BASE ALLOYS
- 46 SPECTROMETRY Solution

BET

- abbr. for Brunauer, Emmett, and Teller (method)

BET SURFACE AREA

- 2 See SURFACE AREA OF POWDERS

BILIRUBIN

- 13 See HEALTH & CLINICAL

BIOLOGICAL

- 9 See FOOD & AGRICULTURE
- 13 See HEALTH & CLINICAL

BIOLOGICAL BUFFER SYSTEMS (ION ACTIVITY)

- 13 HEPES Free Acid
- 13 MOPSO Free Acid
- 13 NaHEPESate
- 13 NaMOPSOate

BIPHENYL

- 77 for DIFFERENTIAL SCANNING CALORIMETRY

BISMUTH

- 46 SPECTROMETRY Solution

BLEACHED KRAFT PULPS

- 7 Northern Softwood
- 7 Eucalyptus Hardwood

BONE ASH

- 13 See HEALTH & CLINICAL
- 95 See NATURAL MATRIX MATERIALS

BONE MEAL

- 13 See HEALTH & CLINICAL

BORATE ORE

- 33 See ORES

BORON

- 85 Implant in Silicon (DEPTH PROFILING)
- 46 SPECTROMETRY Solution

BORIC ACID

- 44 Acidimetric/Assay Values of (STOICHIOMETRY)
- 48 Enriched in Boron-10 (STABLE ISOTOPIC MATERIALS)

BOTANICAL

- 9 See FOOD & AGRICULTURAL

BOVINE

- 9 Liver (FOOD & AGRICULTURE)
- 9 Muscle Powder
- 14 Serum Albumin (HEALTH & CLINICAL)

BRASS

- 60 See NONFERROUS METALS

BROMIDE

- 48 ANION CHROMATOGRAPHY Solution
- 48 Sodium Bromide (STABLE ISOTOPICS)

BROMINE

- 9 in FOOD & AGRICULTURE

BRONZE

- 61 See COPPER BASE ALLOYS

BUFFERS

- 71 See ION ACTIVITY

BURNT REFRACTORIES (ALUMINUM OXIDE)

- 35 See REFRACTORIES

C

CADMIUM

- 36 Cadmium
- Cyclohexanecarboxylate
- 46 SPECTROMETRY Solution
- 79 VAPOR PRESSURE OF METALS

CALCIUM

- 13 Calcium Carbonate (HEALTH & CLINICAL)
- 15 Calcium Hydroxyapatite (BIO-MATERIALS)
- 46 SPECTROMETRY Solution

CALORIMETRY (THERMODYNAMIC PROPERTIES)

- 76 COMBUSTION CALORIMETRY
- 77 DIFFERENTIAL SCANNING CALORIMETRY
- 77 DIFFERENTIAL THERMAL ANALYSIS
- 76 ENTHALPY AND HEAT CAPACITY
- 76 SOLUTION CALORIMETRY

CARBIDES (CERAMICS AND GLASSES)

- 66 Silicon CARBIDE
- 66 Tungsten CARBIDE
- 66 See CEMENTED TUNGSTON CARBIDES

CARBON

- 25 Carbon Modified Silica (INORGANICS)
- 94 Carbon-14 Dating
- 51 in PLAIN CARBON STEELS
- 51 (FERROUS METALS)

CARBON DIOXIDE (PRIMARY GAS MIXTURES)

- 26 Carbon Dioxide in Nitrogen

CARBON MONOXIDE (PRIMARY GAS MIXTURES)

- 26 Carbon Monoxide in Air
- 26 Carbon Monoxide in Nitrogen

B-CAROTENE (FAT SOLUBLE VITAMINS)

- 14 in Human Serum (HEALTH & CLINICAL)

CAST IRON

- 51 See FERROUS METALS

CAST STEEL

- 57 See FERROUS METALS

CATALYST MATERIALS

- 38 High Sulfur Gas Oil Feed (CATALYST CHARACTERIZATION MATERIAL)
- 25 Used Auto Catalysts (INORGANICS)

CEMENTS

- 2 CEMENT TURBIDIMETRY AND FINENESS(SIZING)
- 68 PORTLAND CEMENT CLINKERS
- 68 PORTLAND CEMENTS

CERAMIC MATERIALS (CERAMICS AND GLASSES)

- 66 CARBIDES
- 66 CEMENTED TUNGSTON CARBIDES
- 66 GLASSES
- 35 See REFRACTORIES
- 34 See ROCKS AND MINERALS
- 81 See SPECULAR SPECTRAL REFLECTANCE

CERIUM

- 46 SPECTROMETRY Solution

CESIUM (RADIOACTIVITY)

- 91 as Cesium-137 Burn-up Standard
- 46 SPECTROMETRY Solution

CHARPY

- 5 V-NOTCH TEST BLOCKS

CHEMICAL

- 43 See HIGH PURITY MATERIALS

CHLORIDE

- 48 ANION ION CHROMATOGRAPHY Solution

CHLORINE

- 91 as Chlorine-36 (RADIOACTIVITY)
- 37 in LUBRICATING BASE OILS
- 48 STABLE ISOTOPIC MATERIAL

CHLORO COMPOUNDS (ORGANIC CONSTITUENTS)

- 20 in Biphenyls
- 21 in Cod Liver Oil
- 20 in Halocarbons
- 45 m-Chlorobenzoic Acid (MICROCHEMISTRY)
- 20 in Pesticides
- 20 in Phenols
- 20 in Pollutants

CHOLESTEROL (HEALTH & CLINICAL)

- 10 in Coconut Oil
- 14 in freeze-dried Human Serum
- 14 in frozen Human Serum
- 10 in Whole Egg Powder

CHROMIUM

- 48 as Chromium Nitrate (STABLE ISOTOPIC MATERIALS)
- 34 in CLAYS
- 85 Cr/CrO Thin Film Depth Profile
- 37 Tris (1-phenyl-1,3-butanediol) chromium (III)
- 46 SPECTROMETRY Solution
- 53 in Steels (FERROUS METALS)

CHRYSTOLE

- 99 in ASBESTOS (INDUSTRIAL HYGIENE)

CLAYS

- 34 Brick
- 34 Flint
- 34 Plastic

CLINICAL LABORATORY MATERIALS

- 15 Amino Acids in HCl
- 13 Angiotensin I (Human)
- 14 Anticonvulsant Drug Level Assay
- 14 Antiepilepsy Drug Level Assay
- 13 Bilirubin
- 15 Bone Ash
- 15 Bone Meal
- 14 Bovine Serum Albumin
- 14 Bovine Serum (Inorganic)
- 13 Calcium Carbonate
- 13 Cholesterol
- 14 Cholesterol in Freeze-dried Human Serum
- 13 Cortisol (Hydrocortisone)
- 13 Creatinine
- 14 Electrolytes in Frozen Human Serum

- 13 d-Glucose (Dextrose)

- 14 Glucose in Frozen Human Serum

- 13 Iron Metal

- 14 Human Serum (SERUM MATERIALS)

- 13 Lead Nitrate

- 14 Lead in Blood

- 14 Lipids in Frozen Human Serum

- 13 Lithium Carbonate

- 13 Magnesium Gluconate Dihydrate

- 13 d-Mannitol

- 13 Potassium Chloride

- 13 Sodium Chloride

- 13 Sodium Pyruvate

- 13 Tripalmitin

- 13 Urea

- 13 Uric Acid

- 14 Vitamins (Fat-Soluble) and Cholesterol in Human Serum

- 13 VMA (4-hydroxy-3-methoxymandelic acid)

- 15 Cardiac Troponin

COAL

- 76 for COMBUSTION CALORIMETRY

- 32 Sulfur in (SULFUR IN FOSSIL FUELS)

- 32 TRACE ELEMENTS in

COAL FLY ASH

- 29 TRACE ELEMENTS in

COATING THICKNESS

- 86 Nonmagnetic COPPER AND CHROMIUM ON STEEL

- 85 Tin-Lead Alloy (SOLDER THICKNESS)

COBALT

- 91 as Cobalt-60

(RADIOACTIVITY)

- 46 SPECTROMETRY Solution

COBALT BASE ALLOYS

- 60 NONFERROUS METALS

COCAINE METABOLITE

- 17 See FREEZE-DRIED URINE

COCONUT OIL

- 10 Cholesterol in (FOOD & AGRICULTURE)

COD LIVER OIL

- 21 Organics in (ORGANIC CONSTITUENTS)

COLUMBIA RIVER SEDIMENT

- 95 See NATURAL MATRIX MATERIALS

CONDUCTIVITY, ELECTRICAL PROPERTIES)

- 83 of Electrolytic Iron

CONDUCTIVITY, ELECTROLYTIC (ION ACTIVITY)

- 73 Hydrochloric Acid in Water
- 72 Potassium Chloride in Water
- 72 Sodium Chloride in Water

CONDUCTIVITY, THERMAL (THERMODY- NAMIC PROPERTIES)

- 79 of Electrolytic Iron
- 79 of Graphite

COORDINATE MEASUR- ING MACHINE PROBE PERFORMANCE 6

COPPER

- 36 Bis(1-phenyl-1,3-butane-
dione)copper (II) (METALLO-
ORGANIC COMPOUNDS)
- 60 Brass (COPPER BASE
ALLOYS)
- 61 Bronze (COPPER BASE
ALLOYS)
- 62 Cupro-Nickel (COPPER BASE
ALLOYS)
- 76 ENTHALPY AND HEAT
CAPACITY of
- 51 in FERROUS METALS
- 78 Freezing Point of
(SECONDARY REFERENCE
POINTS)
- 43 High-Purity METALS
(MICROANALYSIS)
- 61 Nickel Silver (COPPER BASE
ALLOYS)
- 60 in NONFERROUS METALS
- 33 in ORES
- 46 SPECTROMETRY Solution
- 48 STABLE ISOTOPES of
- 61 as Unalloyed Copper
(COPPER γ -BENCHMARK)

COPPER BASE ALLOYS

- 60 See NONFERROUS METALS

CORN

- 10 Bran (FOOD &
AGRICULTURE)
- 9 Kernel (FOOD &
AGRICULTURAL)

- 9 Stalk (FOOD &
AGRICULTURAL)

- 10 Starch (See Nutrition
Composition)

CORROSION

- 3 Tool Steel (ABRASIVE WEAR)

CORTISOL (HYDROCORTISONE)

- 13 See HEALTH & CLINICAL

COTININE

- 17 in FREEZE-DRIED URINE

CREATININE

- 13 See HEALTH & CLINICAL

CRUDE OIL

- 29 Vanadium in (METAL
CONSTITUENTS)

CUP FURNACE (FIRE RESEARCH)

- 4 See SMOKE TOXICITY

CURIUM (RADIOACTIVITY)

- 91 as Curium-243
- 91 as Curium-244

CYSTINE

- 45 See MICROCHEMISTRY

D

DENSITY

- 88 of Lead Silica Glass
- 95 Neutron Density Monitor Wire
(RADIATION
DOSIMETRY)
- 4 of Smoke (SMOKE
DENSITY CHAMBER)

DEPTH PROFILING

- 85 Nickel/Chromium
Thin Film
- 85 Arsenic Implant in Silicon
- 85 Boron Implant in Silicon

DEXTROSE (D-GLUCOSE)

- 13 See HEALTH & CLINICAL

DIFFERENTIAL SCANNING CALORIMETRY

- 77 Biphenyl
- 77 Indium
- 77 Mercury
- 77 Thermal Analysis
Purity Set
- 77 Tin

DIFFERENTIAL THERMAL ANALYSIS 77

DIFFRACTION (X-RAY) 89

DIOXIN (IN ISOOCTANE)

- 20⁺ See ORGANIC
CONSTITUENTS

DISODIUM HYDROGEN PHOSPHATE

- 66 for pD CALIBRATION
- 65 for pH CALIBRATION

DNA

- (abbr. for Deoxyribonucleic Acid)

DNA PROFILING

- 15 See HEALTH & CLINICAL
- 17 See FORENSICS
- 17 DNA Profiling
- 17 PCR-Based DNA Profiling
- 15 DNA Mitochondrial
Sequencing

DOLOMITIC LIMESTONE

- 35 See ROCKS AND
MINERALS

DOSIMETRY (RADIOACTIVITY)

- 95 Neutron Density Monitor Wire

DRUG LEVEL ASSAY (ANTIEPILEPSY)

- 14 See HEALTH & CLINICAL

DRUGS OF ABUSE

- 17 in FREEZE-DRIED URINE

DSC

- 77 abbr. for Differential
Scanning Calorimetry

DTA

- 77 abbr. for Differential Thermal
Analysis

DUST

- 98 Urban (TRACE ELEMENTS)
- 21 Urban (ORGANIC
CONSTITUENTS)

DYE PENETRANT TEST (CRACK) BLOCK

- 5 (NONDESTRUCTIVE
EVALUATION)

DYSPROSIUM

- 46 SPECTROMETRY Solution

E

EDDY CURRENT

- 5 ARTIFICIAL FLAW FOR NDE

ELECTRICAL PROPERTIES

- 83 See ELECTRICAL RESISTIVITY AND CONDUCTIVITY OF GRAPHITE & ELECTROLITIC
- 83 See ELECTRICAL RESISTIVITY AND CONDUCTIVITY OF SILICON
- 87 See SUPERCONDUCTING CRITICAL CURRENT
- 87 of GLASS (CERAMICS AND GLASSES)

ELECTROLYTIC CONDUCTIVITY (ION ACTIVITY)

- 73 Hydrochloric Acid Solutions for
- 72 Potassium Chloride Solutions for
- 72 Sodium Chloride Solutions for

ELECTRON MICROSCOPE

- 36 THIN FILM FOR TRANSMISSION ELECTRON MICROSCOPE

ELECTROPHORETIC MOBILITY 73

ELLIPSOMETRY

- 86 Silicon Dioxide on Silicon

ENTHALPY (THERMODYNAMIC PROPERTIES)

- 76 of Copper
- 76 of Molybdenum
- 76 of Synthetic Sapphire
- 76 of Polystyrene

ENVIRONMENTAL MATRICES

- 23 See METAL CONSTITUENTS (INORGANICS)
- 95 See NATURAL MATRIX MATERIALS (RADIOACTIVITY)
- 20 See ORGANIC CONSTITUENTS (ORGANICS)
- 29 See TRACE ELEMENTS IN COALS & COKE

ERBIUM

- 46 SPECTROMETRY Solution

ESTUARINE SEDIMENT

- 23 See (SOILS, SEDIMENTS, AND SLUDGES)

ETHANOL

- 30 Ethanol (ALCOHOLS AND ETHERS IN REFERENCE FUELS)
- 17 Ethanol-Water (ETHANOL SOLUTIONS)

ETHERS (ALCOHOLS AND ETHERS IN REFERENCE FUELS)

- 30 t-Amyl Methyl Ether
- 30 Ethyl t-Butyl Ether
- 30 Methyl t-Butyl Ether

EUCALYPTUS HARDWOOD

- 7 BLEACHED KRAFT PULPS

EUROPIUM

- 91 as Europium-152 (RADIOACTIVITY)
- 46 SPECTROMETRY solution

F

FATTY ACIDS (FOOD & AGRICULTURE)

- 10 Typical Diet

FELDSPAR (ROCKS AND MINERALS)

- 34 in Potash
- 34 in Soda

FERROUS ALLOYS

- 51 See FERROUS METALS

FERTILIZERS (FOOD & AGRICULTURE)

- 11 Ammonium Dihydrogen Phosphate
- 11 Phosphate Rock (Florida & Western)
- 11 Potassium Dihydrogen Phosphate
- 11 Potassium Nitrate

FIBROUS GLASS BOARD

- 78 See THERMAL RESISTANCE OF GLASS, SILICA, AND POLYSTYRENE

FILTER MEDIA (MATERIALS ON FILTER MEDIA)

- 97 Air Particulate on Filter
- 97 Quartz on

FILTERS, OPTICAL 80

FINENESS (SIZING)

- 2 of Portland Cement (CEMENT TURBIDIMETRY AND FINENESS)

FIRE RESEARCH

- 4 FLOORING RADIANT PANEL
- 4 SMOKE DENSITY
- 4 SMOKE TOXICITY
- 3 SURFACE FLAMMABILITY

FISSION TRACK GLASS 95

FLAMMABILITY

- 3 SURFACE FLAMMABILITY (FIRE RESEARCH)

FLOORING RADIANT PANEL 4

- 3 See FIRE RESEARCH

FLOUR

- 9 Durum Wheat
- 9 Hard Red Spring Wheat
- 9 Rice
- 9 Soft Winter Wheat
- 9 Spinach Leaves
- 9 Wheat Hardness

FLUORESCENCE

- 81 Quinine Sulfate Dihydrate
- 81 Raman Spectroscopy

FLUORIDE

- 48 ANION CHROMATOGRAPHY solution
- 15 in FREEZE-DRIED URINE
- 11 in Vegetation

FLUORO COMPOUNDS

- 45 p-Fluorobenzoic Acid (MICRO-CHEMISTRY)

FLUORSPAR (ORES)

- 33 Customs Grade
- 33 High Grade

FLY ASH COAL

- 29 Coal Fly Ash (FOSSIL FUELS)
- 29 TRACE ELEMENTS

FOODS & AGRICULTURE (NUTRITION COMPOSITION)

- 10 Baking Chocolate
- 10 Baby Food Composite
- 10 Corn Bran
- 10 Corn Starch
- 10 Durham Wheat Flour
- 10 Dietary Supplement, Ephedra
- 10 Fatty Acids & Cholesterol
- 10 Infant Formula
- 10 Meat Homogenate
- 10 Peanut Butter
- 10 Typical Diet
- 10 Whole Egg Powder
- 10 Whole Milk Powder
- 10 Wheat Gluten
- 9 Bovine Liver
- 9 Non-fat Milk Powder
- 9 Oyster Tissue
- 9 Rice Flour
- 9 Wheat Flour
- 10 Slurried Spinach

FOSSIL FUELS

- 30 Alcohols & Ethers in Reference Fuels
- 70 Coal Heat of Combustion (COMBUSTION CALORIMETRY)
- 30 Ethanol (ALCOHOLS AND ETHERS IN REFERENCE FUELS)
- 29 Isooctane
- 29 n-Heptane
- 29 METAL CONSTITUENTS in Fossil Fuels
- 32 METAL CONSTITUENTS in Residual Fuel Oil
- 32 Methanol
- 32 Sulfur in Coal (SULFUR IN FOSSIL FUELS)
- 31 Sulfur in Kerosine (SULFUR IN FOSSIL FUELS)
- 32 Sulfur in Residual Fuel Oil (SULFUR IN FOSSIL FUELS)
- 70 Synthetic Refuse Derived Oil (COMBUSTION CALORIMETRY)
- 32 TRACE ELEMENTS in Coal
- 29 TRACE ELEMENTS in Coal Fly Ash
- 29 TRACE ELEMENTS in Fuel Oil

- 29 Vanadium in Crude Oil (METAL CONSTITUENTS IN FOSSIL FUELS)

FREE CUTTING BRASS

- 60 See NONFERROUS METALS

FRESHWATER LAKE SEDIMENT (RADIOACTIVITY)

- 95 Freshwater Lake Sediment (NATURAL MATRIX MATERIALS)

FREEZING POINT (THERMODYNAMIC PROPERTIES)

- 77 of Aluminum (DEFINING FIXED POINT, ITS-90)
- 78 of Copper (SECONDARY REFERENCE POINTS)
- 77 Of Indium (DEFINING FIXED POINT, ITS-90)
- 78 of Lead (REFERENCE POINTS)
- 77 of Silver (DEFINING FIXED POINT, ITS-90)
- 77 of Tin (DEFINING FIXED POINT, ITS-90)
- 77 of Zinc (DEFINING FIXED POINT, ITS-90)

FSV

- 14 abbr. for Fat Soluble Vitamins

FUELS

- 29 See FOSSIL FUELS

FUMED SILICA BOARD 78

- 78 See THERMAL RESISTANCE OF GLASS, SILICA, AND POLYSTYRENE

G

GADOLINIUM

- 46 SPECTROMETRY Solution

GALLIUM

- 24 in Buffalo River Sediment (SOILS, SEDIMENTS, AND SLUDGES)
- 32 in Coal (TRACE ELEMENTS)
- 29 in Coal Fly Ash (TRACE ELEMENTS)
- 93 as Gallium-67 (RADIOPHARMACEUTICALS)
- 67 in Glass (TRACE ELEMENTS)
- 78 Melting Point (THERMODYNAMIC PROPERTIES)

- 48 Metal (STABLE ISOTOPIC MATERIALS)

- 46 SPECTROMETRY Solution

GAS CHROMATOGRAPHY (ORGANIC CONSTITUENTS)

- 19 GC/MS System Performance
- 19 LC Selectivity

GASES (PRIMARY GAS MIXTURES)

- 25 See PRIMARY GAS MIXTURES

GASES IN METALS

- 59 in Irons (FERROUS METALS)
- 59 in Steels (FERROUS METALS)
- 61 in Unalloyed Titanium (NON-FERROUS METALS)

GASOLINE

- 29 See FOSSIL FUELS

GEOLOGICAL

- 32 See GEOLOGICAL MATERIALS AND ORES

GERMANIUM

- 46 SPECTROMETRY Solution

GILDING METAL

- 60 See NONFERROUS METALS

GLASS BEADS

- 1 See SIZING

GLASSES

- 87 Borosilicate (VISCOSITY OF GLASS)
- 67 Chemical Composition
- 87 Chemical Resistance
- 67 Fused Ore Glass
- 89 GLASS LIQUIDUS TEMPERATURE
- 67 High-Boron Borosilicate
- 67 Lead-Barium
- 87 Lead-Silica (ELECTRICAL PROPERTIES OF GLASS)
- 67 Low-Boron Soda-Lime Powder
- 79 LABORATORY THERMOMETER (MERCURY IN GLASS)
- 67 Multi Component
- 88 RELATIVE STRESS OPTICAL COEFFICIENT of
- 35 Sand (ROCKS AND MINERALS)
- 67 Soda-Lime Container
- 67 Soda-Lime Flat
- 67 Soda-Lime Float
- 67 Soda-Lime Sheet

- 88 Soda-Lime-Silica
(VISCOSITY OF GLASS)
- 67 Soft Borosilicate
- 67 SYNTHETIC GLASS (TRACE
ELEMENTS)
- 78 THERMAL EXPANSION OF
METAL & GLASS
- 78 THERMAL RESISTANCE OF
GLASS, SILICA, AND POLY-
STYRENE
- 88 VISCOSITY FIXPOINTS of

GLASS SAND

- 34 See ROCKS AND
MINERALS

GLASS SPHERES

- 1 PARTICLE SIZE (SIZING)

D-GLUCOSE

- 13 aka. Dextrose (HEALTH &
CLINICAL)
- 44 Polarimetric Value of
(STOICHIOMETRY)

GOETHITE

- 73 Aka. A-FeOOH (ELEC-
TROPHORETIC MOBILITY)

GOLD

- 43 METALS (HIGH PURITY MET-
ALS)
- 33 Ore Refractories
- 46 SPECTROMETRY Solution
- 79 VAPOR PRESSURE OF
METALS
- 44 Royal Canadian Mint
Reference Materials
(HIGH PURITY MATERIALS)

GRAPHITE

- 79 THERMAL CONDUCTIVITY OF
GRAPHITE AND IRON

GRAVITY SEDIMENTATION

- 1 Zirconium Oxide
(PARTICLE SIZE)

H

HAFNIUM

- 46 SPECTROMETRY Solution
- 65 in Zircaloy (ZIRCONIUM BASE
ALLOYS)

HARDNESS (FOOD AND AGRICULTURE)

- 9 WHEAT HARDNESS

HARDNESS (SURFACE FINISH)

- 6 of Bright Copper (MICRO-
HARDNESS)
- 6 of Bright Nickel (MICROHARD-
NESS)
- 6 Of Ceramic
(MICROHARDNESS)
- 5 ROCKWELL HARDNESS

HASTELLOY

- 63 NICKEL BASE ALLOYS

HEALTH, NUTRITION COMPOSITION

- 10 Baby Food Composite
- 10 Cholesterol.....in Coconut
Oil
- 10 Fatty Acids.....Frozen Diet
Composite
- 10 Infant Formula (milk-based)
- 10 Typical Diet
- 10 Whole Egg Powder
- 10 Whole Milk

HEAT (THERMODYNAMIC PROPERTIES)

- 76 COMBUSTION
CALORIMETRY
- 77 DEFINING FIXED POINT,ITS-90
- 77 DEFINING FIXED POINT
CELLS,ITS-90
- 77 DIFFERENTIAL SCANNING
CALORIMETRY
- 77 DIFFERENTIAL THERMAL
ANALYSIS
- 76 ENTHALPY AND HEAT
CAPACITY
- 78 FREEZING POINT, MELTING
POINT, AND TRIPPLE POINT
CELLS
- 79 LABORATORY
THERMOMETER
- 78 REFERENCE POINTS
- 76 SOLUTION CALORIMETRY
- 79 THERMAL CONDUCTIVITY
OF GRAPHITE AND IRON
- 78 THERMAL EXPANSION
OF METAL & GLASS
- 78 THERMAL RESISTANCE
OF GLASS, SILICA, AND
POLYSTYRENE

- 79 THERMOCOUPLE
MATERIAL, PLATINUM
- 79 VAPOR PRESSURE OF
METALS

HEPES (BIOLOGICAL BUFFERS)

- 72 abbr. for N-2-Hydroxyethyl-
piperazine-N-2-ethanesulfonic
Acid
- 72 HEPES Free Acid
- 72 NaHEPESate

N-HEPTANE (FOSSIL FUELS)

- 29 REFERENCE LIQUIDS
FOR EVALUATING FUELS

HIGH ALLOY STEELS (FERROUS METALS)

- 54 Chromium Nickel (Copper
Precipitation Hardening)
- 52 Chromium Nickel
(Molybdenum Precipitation
Hardening)
- 56 High Nickel
- 54 High Temperature Alloy
(A286) Nickel-Chromium
- 54 High Temperature Alloy L605
- 58 High Temperature Alloy Iron-
Nickel-Cobalt
- 54 Valve Steel

HIGH PURITY METALS

- 43 High Purity Gold
- 43 High Purity Platinum
- 43 High Purity Zinc
- 43 Refined Copper
- 43 Selenium Intermediate Purity
- 43 Zinc Intermediate Purity
- 43 Zinc Metal

HIGH TEMPERATURE ALLOYS

- 51 See FERROUS METALS

HOLMIUM

- 81 Holmium Oxide Solution
Wavelength
- 46 SPECTROMETRY Solution

HUMAN

- 13 See HEALTH & CLINICAL
- 95 LIVER (NATURAL
MATRIX MATERIALS)
(RADIOACTIVITY)
- 95 LUNG (NATURAL
MATRIX MATERIALS)
(RADIOACTIVITY)
- 14 Serum (SERUM
MATERIALS)

HUMAN SERUM (HEALTH & CLINICAL)

- 14 Cholesterol in Human Serum
- 14 Electrolytes in (SERUM MATERIALS)
- 14 Fat Soluble Vitamins in
- 14 Glucose in Frozen (SERUM MATERIALS)
- 14 Lipids in Frozen (SERUM MATERIALS)
- 14 SERUM MATERIALS

HUMAN SERUM (ORGANICS)

- 22 Polychlorinated Biphenyls in (ORGANIC CONSTITUENTS)

HYDROGEN

- 91 as Hydrogen-3 (RADIOACTIVITY SOLUTIONS)
- 64 Unalloyed Titanium for (GASES IN METALS)

HYDROXYAPATITE

- 15 See Calcium Hydroxyapatite

4-HYDROXY-3-METHOXY- DL-MANDELIC ACID (VMA) 13

I

ICTAC

- 77 abbr. for International Confederation of Thermal Analysis and Calorimetry
- 77 X-RAY AND PHOTOGRAPHY

INCONEL

- 63 NICKEL BASE ALLOYS (NON-FERROUS METALS)

INDIUM

- 93 as Indium-111 (RADIOPHARMACEUTICALS)
- 77 DEFINED FIXED POINT, ITS-90
- 78 FREEZING POINT, MELTING POINT, AND TRIPLE POINT CELLS
- 46 SPECTROMETRY Solution

INDUSTRIAL HYGIENE

- 97 See INDUSTRIAL HYGIENE

INFRARED, NEAR

- 81 INFRARED REFLECTANCE

IODINE (RADIOACTIVITY)

- 93 as Iodine-125 (RADIOPHARMACEUTICALS)
- 48 Iodine, Isotopic

- 91 as Iodine-129 (RADIOACTIVE LUTIONS)
- 93 as Iodine-131 (RADIOPHARMACEUTICALS)

ION ACTIVITY

- 13,72 BIOLOGICAL BUFFER SYSTEMS
- 73 ELECTROLYTIC CONDUCTIVITY
- 72 ION-SELECTIVE ELECTRODE CALIBRATION
- 72 pD CALIBRATION
- 71 pH CALIBRATION

IRON

- 79 Electrolytic Iron (THERMAL CONDUCTIVITY OF GRAPHITE AND IRON)
- 51 See FERROUS METALS
- 13 Iron Metal (HEALTH & CLINICAL)
- 46 SPECTROMETRY Solution
- 37 Tris(1-phenyl-1-3 butaine-diono)-iron(III) (METALLO-ORGANIC COMPOUNDS)

ISOTOPE(S)

- 49 See LIGHT STABLE ISOTOPIC MATERIALS
- 43 See HIGH PURITY MATERIALS
- 94 See RADIOACTIVITY

K

KEROSINE

- 31 Sulfur in (SULFUR IN FOSSIL FUELS)

KNOOP MICROHARD- NESS (SURFACE FINISH)

- 6 Bright Copper
- 6 Bright Nickel
- 6 Silicon Nitride

L

LANTHANUM

- 46 SPECTROMETRY Solution

LAKE SEDIMENT (RADIOACTIVITY)

- 95 Freshwater Lake Sediment (NATURAL MATRIX MATERIALS)

LEAD

- 36 Lead Cyclohexanecarboxylate (METALLO-ORGANIC COMPOUNDS)
- 14 Lead in Blood (HEALTH & CLINICAL)
- 13 Lead Nitrate (HEALTH & CLINICAL)
- 48 Metal Equal Atom (STABLE ISOTOPIC MATERIALS)
- 48 Metal, Natural (STABLE ISOTOPIC MATERIALS)
- 48 Metal, Radiogenic (STABLE ISOTOPIC MATERIALS)
- 98 In Paint Film
- 98 In Powdered Paint
- 98 In Indoor Dust, Trace Elements
- 98 In Paint on Fiberboard
- 98 In Soil, Trace Elements
- 60 See NONFERROUS METALS
- 98 Powdered Lead Base Paint (LEAD IN PAINT, DUST AND SOIL)
- 29 in Reference Fuel (METAL CONSTITUENTS IN FOSSIL FUELS)
- 46 SPECTROMETRY Solution

LEAD BASE ALLOYS/MATERIALS

- 60 See NONFERROUS METALS

LEAVES (FOOD & AGRICULTURE)

- 11 Apple
- 11 Peach
- 11 Pine Needles
- 11 Spinach
- 11 Tomato

LIMESTONE (ROCKS AND MINERALS)

- 35 Argillaceous
- 35 Dolomitic

LINERBOARD

- 6 for TAPE ADHESION TESTING

LINEWIDTH (METROLOGY)

- 84 OPTICAL MICROSCOPE LINEWIDTH MEASUREMENT

LIPIDS

- 14 in Human Serum (SERUM MATERIALS)

LIQUID CHROMATOGRAPHY

- 19 GS/MS AND LC SYSTEM PERFORMANCE

LIQUIDUS TEMPERATURE

- 89 Soda-Lime Silica
- 89 Aluminosilicate

LITHIUM

- 49 Carbonate (LIGHT STABLE ISOTOPIC MATERIALS)
- 13 Carbonate (HEALTH & CLINICAL)
- 36 Lithium Cyclohexanebutyrate (METALLO-ORGANIC COMPOUNDS)
- 33 Ore, Lepidolite
- 33 Ore, Petalite (ORES)
- 33 Ore, Spodumene (ORES)
- 46 SPECTROMETRY Solution

LIVER

- 9 Bovine (FOODS AND BEVERAGES)
- 95 Human (NATURAL MATRIX MATERIALS) (RADIOACTIVITY)

LUBRICATING BASE OIL

- 37 Total Chlorine
- 37 Total Nitrogen
- 37 Total Sulfur
- 38 WEAR-METALS IN OIL

LUNG (RADIOACTIVITY)

- 95 Human (NATURAL MATRIX MATERIALS)

LUTETIUM

- 46 SPECTROMETRY Solution

M

MAGNETIC MOMENT

- 7 Nickel Disk
- 7 Nickel Sphere
- 7 Yttrium Garnet Sphere

MAGNESIUM

- 13 Magnesium Gluconate Dihydrate (HEALTH & CLINICAL)
- 48 Magnesium Metal (STABLE ISOTOPIC MATERIALS)
- 46 SPECTROMETRY Solution

MAGNIFICATION

- 85 SCANNING ELECTRON MICROSCOPE (SEM)

MANGANESE

- 85 SEM Performance Standard
- 85 SEM Sharpness Standard
- 46 SPECTROMETRY Solution

D-MANNITOL (HEALTH & CLINICAL) 13

MARIJUANA METABOLITE

- 17 THC-9-COOH (DRUGS OF ABUSE IN URINE)

MARINE MATERIALS

- 24 Buffalo River Sediment (METAL CONSTITUENTS IN NATURAL MATRICES)
- 24 Estuarine Sediment (METAL CONSTITUENTS IN NATURAL MATRICES)
- 24 Marine Sediment
- 35 Limestone Argillaceous
- 34 (ROCKS AND MINERALS)
- 35 Limestone Dolomitic (ROCKS AND MINERALS)
- 21 Organics in Marine Sediment (ORGANIC CONSTITUENTS)
- 21 Organics in Mussel Tissue (ORGANIC CONSTITUENTS)
- 21 Organics in Whale Blubber (ORGANIC CONSTITUENTS)
- 9 Oyster Tissue (FOOD & AGRICULTURE)
- 22 Polychlorinated Biphenyls (Congeners) in River Sediment A (ORGANIC CONSTITUENTS)
- 24 Sediment for Solid Sampling

MASS SPECTROMETRY

- 19 GC/MS AND LC SYSTEM PERFORMANCE (ORGANICS)
- 19 GC/MS SYSTEM
- 19 Lc Chiral Selectivity
- 19 LC Performance
- 19 LC Selectivity
- 45 See LIGHT STABLE ISOTOPIC MATERIALS
- 91 See RADIOACTIVITY
- 48 See STABLE ISOTOPIC MATERIALS

MATERIALS ON FILTER MEDIA

- 97 Quartz on Filter Media
- 97 Air Particulate on Filter
- 97 Cellulose Acetate Membrane
- 97 Ashless Blank Filter
- 97 Respirable Alpha Quartz
- 97 Respirable Cristobalite

MELTING POINT AND TRIPLE POINT (THERMODYNAMIC PROPERTIES) 78

MERCURY

- 15 Mercury (TOXIC SUBSTANCES IN URINE)
- 77 Mercury (Triple Point) (DEFINING FIXED POINT ITS-90)
- 46 SPECTROMETRY Solution
- 29 TRACE ELEMENTS (FOSSIL FUELS)
- 29 Trace Mercury in Coal (TRACE ELEMENTS)
- 24 in Water (METAL CONSTITUENTS IN NATURAL MATRICES)

METALLO-ORGANICS (ENGINE WEAR MATERIALS) 36

METALS ON FILTER MEDIA

- 97 See MATERIALS ON FILTER MEDIA

METHANE (PRIMARY GAS MIXTURES)

- 27 Methane in Air

METROLOGY 84

MICROANALYSIS 35

MICROCHEMISTRY (HIGH PURITY MATERIALS)

- 45 Acetanilide
- 45 Anisic Acid
- 45 m-Chlorobenzoic Acid
- 45 Cystine
- 45 p-Fluorobenzoic Acid
- 45 Nicotinic Acid
- 45 Urea

MICROCOPY

- 82 Microcopy Resolution Test Chart(X-RAY AND PHOTOGRAPHY)

MICROHARDNESS (SURFACE FINISH)

- 6 of Bright Copper
- 6 of Bright Nickel
- 6 of Ceramic

MICROSCOPY (METROLOGY)

- 85 DEPTH PROFILING
- 86 ELLIPSOMETRY
- 84 OPTICAL MICROSCOPE LINEWIDTH MEASUREMENT
- 85 SCANNING ELECTRON MICROSCOPE (SEM)

MICROSPHERE (SIZING)

- 1 Glass Spheres (PARTICLE SIZE)
- 1 Polystyrene Spheres (PARTICLE SIZE)

MILK (FOOD AND AGRICULTURE)

- 10 Infant Formula
- 9 Non-fat Milk Powder

MINERALS

- 34 See ROCKS AND MINERALS

MIXTURES AND POLLUTANTS (PRIMARY GAS MIXTURES)

- 25 Ambient Non-Methane Organics in Nitrogen
- 26 Carbon Dioxide in Nitrogen
- 25 Carbon Monoxide in Air
- 26 Carbon Monoxide in Nitrogen
- 27 Hydrogen Sulfide in Nitrogen
- 27 Methane in Air
- 27 Nitric Oxide in Nitrogen
- 27 Oxides of Nitrogen in Air
- 28 Oxygen in Nitrogen
- 28 Propane in Air
- 28 Sulfur Dioxide in Nitrogen

MOLECULAR WEIGHT AND MELT FLOW (POLYMERIC PROPERTIES)

- 75 Polyethylene Gas Pipe Resin
- 74 Polyethylene Linear
- 74 Poly(ethylene oxide)
- 75 Polyethylene Resin
- 74 Poly(methylmethacrylate)
- 74 Polystyrene

MOLYBDENUM

- 76 ENTHALPY AND HEAT CAPACITY
- 93 as Molybdenum-99-Technetium-99m (RADIO-PHARMACEUTICALS) 87
- 46 SPECTROMETRY Solution

N

NAVAL BRASS

- 60 See NONFERROUS METALS

NDE

- 5 abbr. for Nondestructive Evaluation

NEODYMIUM

- 46 SPECTROMETRY Solution

NEUTRON MONITOR (RADIOACTIVITY)

- 95 Neutron Density Monitor Wire (RADIATION DOSIMETRY)

NICKEL

- 91 as Nickel-63 (RADIOACTIVE SOLUTION)
- 36 Nickel
- 36 Cyclohexanecarboxylate(METAL-LO-ORGANIC COMPOUNDS)
- 48 Nickel (STABLE ISOTOPIC MATERIALS)
- 85 Nickel-Chromium Thin Film(DEPTH PROFILING)
- 63 NICKEL BASE ALLOYS (NON-FERROUS METALS)
- 63 NICKEL OXIDES (NONFERROUS METALS)
- 7 Nickel Disk (MAGNETIC MOMENT)
- 7 Nickel Sphere (MAGNETIC MOMENT)
- 46 SPECTROMETRY Solution

NICOTINIC ACID

- 45 MICROCHEMISTRY (HIGH PURITY MATERIALS)

NIOBIUM

- 93 as Niobium-94 (GAMMA RAY POINT SOURCES)
- 46 SPECTROMETRY Solution

NITRATE

- 48 ANION CHROMATOGRAPHY Solution

NITRIC OXIDE (PRIMARY GAS MIXTURES)

- 27 Nitric Oxide in Nitrogen

NITRIDE

- 2 Silicon Nitride (SURFACE AREA OF POWDERS)
- 6 (MICROHARDNESS)

NITROGEN (PRIMARY GAS MIXTURES)

- 37 Total Nitrogen (LUBRICATING BASE OILS)

NONDESTRUCTIVE EVALUATION

- 5 ARTIFICIAL FLAW FOR EDDY CURRENT NDE

NONFERROUS ALLOYS

- 60 See NONFERROUS METALS

NORTHERN SOFTWOOD

- 7 BLEACHED KRAFT PULPS

NUCLEAR MATERIALS (RADIOACTIVITY)

- 94 Carbon-14 DATING
- 95 FISSION TRACK GLASS
- 95 NATURAL MATRIX MATERIALS
- 91 RADIOACTIVE SOLUTIONS
- 93 RADIOPHARMACEUTICALS
- 94 RADON EMANATION

NUTRITION

- 9 See FOOD & AGRICULTURE

O

OBSIDIAN ROCK

- 34 ROCKS AND MINERALS

OCEAN MATERIALS (RADIOACTIVITY) (NATURAL MATRIX MATERIALS)

- 95 Ocean Sediment

OIL

- 37 Chlorine in (LUBRICATING BASE OILS)
- 31 Fuel Oil (FOSSIL FUELS)
- 38 High Sulfur Gas Oil Feed (CATALYST CHARACTERIZATION MATERIALS)
- 32 Moisture in Oils (FOSSIL FUELS)
- 37 Nitrogen (LUBRICATING BASE OILS)
- 21 Organics in Cod Liver Oil (ORGANIC CONSTITUENTS)
- 21 Petroleum Crude Oil (ORGANIC CONSTITUENTS)
- 22 Polychlorinated Biphenyls in (ORGANIC CONSTITUENTS)
- 21 Shale Oil (ORGANIC CONSTITUENTS)
- 37 Sulfur in (LUBRICATING BASE OILS)
- 32 Sulfur in Residual Fuel Oil (SULFUR IN FOSSIL FUELS)
- 29 Vanadium in Crude Oil (METAL CONSTITUENTS IN FOSSIL FUELS)
- 38 WEAR-METALS IN OIL (ENGINE WEAR MATERIALS)

OPTOELECTRONICS (METROLOGY)

- 84 Optical Fiber Coating
- 84 Optical Fiber Diameter
- 84 Optical Fiber Ferrule Geometry
- 84 Pin Gauge for Optical Fiber Ferrules
- 84 Polarization Mode Dispersion
- 84 Wavelength Reference Absorption Cell

ORES (GEOLOGICAL MATERIALS AND ORES)

- 33 Alumina (Reduction Grade)
- 33 Bauxite, Arkansas
- 33 Bauxite, Dominican
- 33 Bauxite, Jamaican
- 33 Bauxite, Surinam
- 33 Borate Ore
- 33 Chinese Ores
- 33 Copper Ore Mill Heads
- 33 Copper Ore Mill Tails
- 33 Fluorspar, Customs Grade
- 33 Fluorspar, High Grade
- 33 Gold Ore, Refractory
- 33 Iron Ore, Canada
- 33 Iron Ore, Labrador
- 33 Iron Ore, Nimba
- 33 Iron Oxide Reduced
- 33 Lithium Ore (Petalite)
- 33 Lithium Ore (Spodumene)
- 33 Lithium Ore (Lepidolite)
- 33 Manganese Ore
- 33 Phosphate Rock Florida
- 33 Phosphate Rock Western
- 33 Pyrite Ore (ORE BIOLEACHING SUBSTRATE)
- 33 Rutile Ore
- 33 Scheelite Ore
- 33 Tungsten Concentrate
- 33 Zinc

ORGANICS

- 22 EPA: ORGANIC COMPOUNDS RELATED TO (WATER ANALYSIS)
- 20 ORGANIC CONSTITUENTS
- 19 GC/MS AND LC SYSTEM PERFORMANCE

OXALIC ACID (RADIOACTIVITY)

- 94 Carbon-14 Dating

OXYGEN (PRIMARY GAS MIXTURES)

- 28 Oxygen in Nitrogen

OXYGENATES

- 30 ALCOHOLS...IN REFERENCE FUELS

OYSTER TISSUE

- 9 FOOD & AGRICULTURE

P

PAINT

- 98 LEAD IN PAINT, DUST AND SOIL

PALLADIUM

- 46 SPECTROMETRY Solution

PARTICAL COUNT MATERIALS

- 2 For suspensions

PARTICLE SIZE (SIZING)

- 1 Glass Spheres
- 1 Polystyrene Spheres
- 1 Silicon Nitride
- 1 Zirconium Oxide

PARTICULATES

- 21 Diesel Particulate Matter (ORGANIC CONSTITUENTS)
- 97 MATERIALS ON FILTER MEDIA
- 22 Urban Dust/Organics (ORGANIC CONSTITUENTS)
- 23 Urban Particulate Matter (INORGANICS)
- 66 pH CALIBRATION (ION ACTIVITY)
- 72 Disodium Hydrogen Phosphate
- 72 Potassium Dihydrogen Phosphate
- 72 Potassium Hydrogen Phthalate
- 72 Sodium Bicarbonate
- 72 Sodium Carbonate

PERUVIAN SOIL (RADIOACTIVITY) 89

PESTICIDES (ORGANIC CONSTITUENTS)

- 20 Chlorinated Pesticides in Hexane
- 20 Chlorinated Pesticides in Isooctane

PH CALIBRATION (ION ACTIVITY)

- 71 Calcium Carbonate
- 71 Disodium Hydrogen Phosphate
- 71 Potassium Dihydrogen Phosphate
- 71 Potassium Hydrogen Phthalate
- 71 Potassium Hydrogen Tartrate
- 71 Potassium Tetroxalate
- 71 Sodium Bicarbonate
- 71 Sodium Carbonate
- 71 Sodium Tetraborate Decahydrate
- 72 See BIOLOGICAL BUFFER SYSTEMS

PHOSPHATE

- 72 See pH CALIBRATION
- 71 See pH CALIBRATION
- 48 ANION CHROMATOGRAPHY Solution
- 33 Phosphate Rock (ORES)
- 37 Triphenyl Phosphate (METAL-LO-ORGANIC COMPOUNDS)

PHOSPHORUS

- 46 SPECTROMETRY Solution

PHOTOGRAPHY

- 82 See X-RAY AND PHOTOGRAPHY

PINE NEEDLES

- 11 See FOOD & AGRICULTURAL

PLASTIC

- 74 See POLYMERIC PROPERTIES

PLATINUM (HIGH PURITY METALS)

- 43 High Purity Platinum
- 46 SPECTROMETRY Solution

PLUTONIUM (RADIOACTIVITY)

- 95 Ashed Bone
- 95 Columbia River Sediment
- 95 Human Liver
- 95 Human Lung
- 95 Ocean Sediment
- 95 Peruvian Soil
- 91 Plutonium-238
- 91 Plutonium-239
- 91 Plutonium-240
- 91 Plutonium-241
- 91 Plutonium-242

POLLUTANTS

- 29 METAL CONSTITUENTS IN FOSSIL FUELS
- 25 PRIMARY GAS MIXTURES
- 20 ORGANIC CONSTITUENTS (ORGANICS)

POLONIUM (RADIOACTIVITY)

- 91 Polonium-209 (RADIOACTIVE SOLUTIONS)

POLYCHLORINATED BIPHENYLS PCBS

- 20 Chlorinated Biphenyls
- 20 Chlorinated Biphenyl Congeners in Isooctane
- 22 Polychlorinated Biphenyl Congeners in Isooctane
- 22 Polychlorinated Biphenyls in River Sediment

POLYETHYLENE (MOLECULAR WEIGHT AND MELT FLOW)

- 75 Polyethylene Gas Pipe Resin
- 74 Polyethylene Linear
- 74 Poly(ethylene Oxide)
- 75 Polyethylene Resin

POLYMER

- 74 See POLYMERIC PROPERTIES

POLY(METHYL-METHACRYLATE) (POLYMERIC PROPERTIES)

- 75 MOLECULAR WEIGHT AND MELT FLOW

POLYSTYRENE

- 76 ENTHALPY AND HEAT CAPACITY
- 76 (THERMODYNAMIC PROPERTIES)
- 74 MOLECULAR WEIGHT AND MELT FLOW
- 74 (POLYMERIC PROPERTIES)

POTASSIUM

- 46 SPECTROMETRY Solution

POTASSIUM CHLORIDE

- 13 See PURE CRYSTALLINE STANDARDS
- 73 ELECTROLYTIC CONDUCTIVITY
- 72 ION-SELECTIVE ELECTRODE CALIBRATION
- 48 STABLE ISOTOPIC MATERIALS
- 76 SOLUTION CALORIMETRY
- 44 STOICHIOMETRY

POTASSIUM DICHROMATE

- 80 MOLECULAR ABSORPTION
- 44 STOICHIOMETRY

POTASSIUM DIHYDROGEN PHOSPHATE

- 11 FERTILIZERS
- 72 pD CALIBRATION
- 71 pH CALIBRATION

POTASSIUM FLUORIDE

- 72 ION-SELECTIVE ELECTRODE CALIBRATION

POTASSIUM HYDROGEN PHTHALATE

- 72 pD CALIBRATION
- 71 pH CALIBRATION
- 44 STOICHIOMETRY

POTASSIUM HYDROGEN TARTRATE

- 71 pH CALIBRATION

POTASSIUM IODIDE

- 80 MOLECULAR ABSORPTION

POTASSIUM NITRATE

- 11 FERTILIZERS
- 49 LIGHT STABLE ISOTOPIC MATERIALS

POTASSIUM TETROXALATE

- 71 pH CALIBRATION

POWDERED LEAD BASE PAINT

- 98 LEAD IN PAINT, DUST, AND SOIL

PRASEODYMIUM

- 46 SPECTROMETRY Solution

PRIMARY CHEMICALS

- 44 STOICHIOMETRY

PRIORITY

POLLUTANT PAH

- 21 ORGANIC CONTAMINANTS

PYRITE ORE

- 34 ORE BIOLEACHING SUBSTRATE

Q

QUARTZ

- 97 MATERIALS ON FILTER MEDIA

R

RADIOACTIVITY

- 95 FISSION TRACK GLASS
- 95 NATURAL MATRIX MATERIALS
- 91 RADIOACTIVE SOLUTIONS
- 93 RADIOPHARMACEUTICALS
- 94 RADON EMANATION
- 94 Carbon-14 DATING

RADIUM (RADIOACTIVITY)

- 92 Radium-226 (RADIOACTIVE SOLUTIONS)
- 94 Radium-226 (RADON EMANATION)

REFERENCE FUELS

- 29 See FOSSIL FUELS

REFLECTANCE (OPTICAL PROPERTIES)

- 81 DIFFUSE SPECTRAL REFLECTANCE
- 81 INFRARED REFLECTANCE
- 81 SPECULAR SPECTRAL REFLECTANCE

REFRACTORIES (GEOLOGICAL MATERIALS AND ORES)

- 35 Burnt Refractory

REFORMULATED GASOLINES

- 29 See FOSSIL FUELS

RESIDUAL RESISTIVITY RATIO (ELECTRICAL PROPERTIES) 83

RESISTANCE (THERMODYNAMIC PROPERTIES)

- 78 THERMAL RESISTANCE OF
GLASS, SILICA, AND POLY-
STYRENE

RESISTIVITY (ELECTRICAL PROPERTIES)

- 83 ELECTRICAL RESISTIVITY
AND CONDUCTIVITY OF
METALS
- 83 ELECTRICAL RESISTIVITY
AND CONDUCTIVITY OF
SILICON

RHENIUM

- 46 SPECTROMETRY Solution

RHODIUM

- 46 SPECTROMETRY Solution

RICE FLOUR (FOOD & AGRICULTURE) 9

RIVER SEDIMENT (INORGANICS)

- 24 SOILS, SEDIMENTS, AND
SLUDGES
- 24 Buffalo River Sediment
- 24 Estuarine Sediment

RIVER SEDIMENT (ORGANICS)

- 22 Polychlorinated Biphenyls in
River Sediment A

RIVER SEDIMENT (RADIOACTIVITY)

- 95 Columbia River Sediment
(NATURAL MATRIX
MATERIALS)

ROCKS

- 34 Basalt Rock (ROCKS AND
MINERALS)
- 35 Obsidian Rock (ROCKS AND
MINERALS)
- 11 Phosphate Rock (Florida)
(FERTILIZERS)
- 11 Phosphate Rock (Western)
(FERTILIZERS)

ROYAL CANADIAN MINT REFERENCE MATERIALS 44

RUBIDIUM

- 78 Rubidium (FREEZING POINT,
MELTING POINT AND TRIPLE
POINT CELLS)
- 48 Rubidium Chloride (STABLE
ISOTOPIC MATERIALS)
- 46 SPECTROMETRY Solution

S

SAMARIUM

- 46 SPECTROMETRY Solution

SAND (GLASS)

- 34 See ROCKS AND MINERALS

SCANDIUM

- 46 SPECTROMETRY Solution

SCANNING ELECTRON MICROSCOPE (METROLOGY)

- 85 SEM Performance Standard
- 85 SEM Sharpness Standard

SHEELITE ORE

- 33 ORES

SEDIMENT

- 23 METAL CONSTITUENTS IN
NATURAL MATRICES
- 95 NATURAL MATRIX MATERIALS
(RADIOACTIVITY)

SELENIUM

- 47 Selenium Intermediate Purity
(HIGH PURITY METALS)
- 46 SPECTROMETRY Solution

SERUM MATERIALS

- 14 Bovine Serum Albumin
- 14 Electrolytes in Frozen Human
Serum
- 14 Glucose in Frozen Human
Serum
- 14 Human Serum
- 14 Lipids in Frozen Human
Serum

SHELLFISH

- 21 Mussel Tissue (ORGANIC
CONSTITUENTS)
- 10 Oyster Tissue (FOOD & AGRI-
CULTURE)

SILICA

- 25 Carbon Modified Silica
(INORGANICS)
- 78 Fumed Silica Board (THER-
MAL RESISTANCE OF GLASS,
SILICA, AND POLYSTYRENE)
- 78 THERMAL EXPANSION OF
GLASS AND SILICA)
- 89 Lead Silica Glass (DENSITY
AND REFRACTIVE INDEX)
- 97 Respirable Alpha Quartz
(RESPIRABLE SILICA)
- 97 Respirable Cristobalite
(RESPIRABLE SILICA)
- 35 Silica Brick (REFRACTORIES)

SILICON

- 83 ELECTRICAL RESISTIVITY
AND CONDUCTIVITY OF
SILICON
- 36 Octaphenylcyclotetrasiloxane
(METALLO-ORGANIC
COMPOUNDS)
- 57 Silicon Metal (STEELMAKING
ALLOYS)
- 89 Silicon Powder (X-RAY DIF-
FRACTION)
- 46 SPECTROMETRY Solution
- 57 See STEELMAKING ALLOYS

SILICON DIOXIDE

- 86 Thin Film Thickness
(ELLIPSOMETRY)

SILICON NITRIDE (SIZING) (SURFACE FINISH)

- 1 PARTICLE SIZE
- 2 SURFACE AREA OF POWDERS
- 6 MICROHARDNESS

SILVER

- 35 Alloy (METALS)
(MICROANALYSIS)
- 44 Royal Canadian Mint
Reference Materials
- 36 Silver 2-ethylhexanoate (MET-
ALLO-ORGANIC MATERIALS)
- 48 Silver Nitrate (STABLE
ISOTOPIC MATERIALS)
- 46 SPECTROMETRY Solution

SINUSOIDAL ROUGHNESS

- 3 SURFACE ROUGHNESS
(SURFACE FINISH)

SIZING

- 2 CEMENT TURBIDIMETRY AND
- 2 FINENESS
- 1 PARTICLE SIZE
- 2 SURFACE AREA OF POWDERS

SLUDGE

- 24 Domestic Sludge (METAL CONSTITUENTS IN NATURAL MATRICES)
- 24 Industrial Sludge (METAL CONSTITUENTS IN NATURAL MATRICES)
- 23 SOILS, SEDIMENTS, AND SLUDGES (GEOLOGICAL MATERIALS AND ORES)

SMOKE

(FIRE RESEARCH)

- 4 SMOKE DENSITY CHAMBER
- 4 SMOKE TOXICITY

SODA LIME GLASS (CERAMICS AND GLASSES)

- 67 Soda-Lime, Container (GLASSES)
- 67 Soda-Lime, Flat (GLASSES)
- 67 Soda-Lime, Float (GLASSES)
- 67 Soda-Lime, Sheet (GLASSES)

SODIUM

- 72,71 Disodium Hydrogen Phosphate
- 72,71 Sodium Bicarbonate (ION ACTIVITY)
- 44 Sodium Carbonate (STOICHIOMETRY)
- 71 Sodium Carbonate (ION ACTIVITY)
- 13 Sodium Chloride (HEALTH & CLINICAL)
- 37 Sodium Cyclohexanecarboxylate (METALLO-ORGANIC MATERIALS)
- 44 Sodium Oxalate (STOICHIOMETRY)
- 13 Sodium Pyruvate (HEALTH & CLINICAL)
- 71 Sodium Tetraborate Decahydrate
- 46 SPECTROMETRY Solution

SOILS

- 23 METAL CONSTITUENTS IN NATURAL MATRICES
- 23 SOILS, SEDIMENTS, AND SLUDGES

SOLDER (METROLOGY)

- 85 Tin-Lead Alloy (SOLDER THICKNESS)

SPECTRAL REFLECTANCE (OPTICAL PROPERTIES)

- 81 SPECULAR SPECTRAL REFLECTANCE

SPHERES (SIZING)

- 1 PARTICLE SIZE

SPECTROMETRY

- 46 SINGLE ELEMENT solutions
- 80 See MOLECULAR ABSORPTION

SPECULAR SPECTRAL REFLECTANCE

- 81 First Surface, Aluminum on Glass

STAINLESS STEEL

- 51 See FERROUS METALS

STEEL COATINGS

- 86 CHROMIUM OVER COPPER ON STEEL

STEELS (FERROUS METALS) 51

STRONTIUM

- 92 Strontium-90 (RADIOACTIVITY)
- 44 Strontium Carbonate (STOICHIOMETRY)
- 48 Strontium Carbonate (STABLE ISOTOPIC MATERIALS)
- 37 Strontium Cyclohexanecarboxylate (METALLO-ORGANIC MATERIALS)
- 46 SPECTROMETRY Solution
- 78 SUCCINONITRILE (THERMODYNAMIC PROPERTIES)

SUCROSE

- 82 OPTICAL ROTATION
- 44 STOICHIOMETRY

SULFATE

- 48 ANION CHROMATOGRAPHY Solution

SULFIDE (PRIMARY GAS MIXTURES)

- 27 Hydrogen Sulfide in Nitrogen

SULFUR

- 38 CATALYST CHARACTERIZATION MATERIALS
- 46 SPECTROMETRY Solution
- 31 SULFUR IN FOSSIL FUELS
- 38 WEAR-METALS IN OIL

SULFUR DIOXIDE (PRIMARY GAS MIXTURES)

- 28 Sulfur Dioxide in Nitrogen
- 3 SURFACE FINISH
- 3 ABRASIVE WEAR
- 6 MICROHARDNESS
- 3 SURFACE ROUGHNESS

SURFACE FINISH

- 3 Abrasive Wear
- 3 Surface Roughness

SURFACE FLAMMABILITY (FIRE RESEARCH)

- 3 Hardboard Sheet

T

TANTALUM

- 46 SPECTROMETRY Solution

TAPE ADHESION TESTING

- 6 Linerboard for Tape Adhesion Testing

TECHNETIUM

- 92 Technetium-99 (RADIOACTIVE SOLUTIONS)
- 93 Technetium-99m (RADIO-PHARMACEUTICALS)

TELLURIUM

- 46 SPECTROMETRY Solution

TERBIUM

- 46 SPECTROMETRY Solution
- 17 TETRAHYDROCANNABINOL (Marijuana Metabolite)
- 17 DRUGS OF ABUSE IN URINE, SINGLE ANALYTE
- 17 DRUGS OF ABUSE IN URINE, MULTIANALYTE

THALLIUM

- 46 SPECTROMETRY Solution
- 93 Thallium-201 (RADIOPHARMACEUTICALS)

THERMAL ANALYSIS (THERMODYNAMIC PROPERTIES)

- 76 COMBUSTION CALORIMETRY
- 77 DIFFERENTIAL SCANNING CALORIMETRY
- 77 DIFFERENTIAL THERMAL ANALYSIS
- 76 ENTHALPY AND HEAT CAPACITY
- 76 SOLUTION CALORIMETRY

THERMAL CONDUCTIVITY OF GRAPHITE AND METALS

- 76 (THERMODYNAMIC PROPERTIES)
- 79 Electrolytic Iron
- 79 Graphite

THERMAL EXPANSION OF METAL GLASS AND SILICA

- 78 Borosilicate Glass
- 78 Copper
- 78 Stainless Steel (AISI 446)

THERMAL RESISTANCE OF GLASS, SILICA, AND POLYSTYRENE

- 78 Expanded Polystyrene Board
- 78 Fibrous Glass Board
- 78 Fumed Silica Board

THERMOMETER (THERMODYNAMIC PROPERTIES)

- 79 Laboratory Thermometer

THERMOMETRIC FIXED POINTS (THERMODYNAMIC PROPERTIES)

- 77 DEFINING FIXED POINT, ITS-90
- 77 DEFINING FIXED POINT CELLS, ITS-90
- 78 FREEZING POINT, MELTING POINT AND TRIPLE POINT

THIANTHRENE

- 76 COMBUSTION CALORIMETRY

THICKNESS (METROLOGY)

- 86 CHROMIUM OVER COPPER ON STEEL
- 86 ELLIPSOMETRY
- 85 SOLDER THICKNESS

THORIUM

- 46 SPECTROMETRY Solution

THORIUM (RADIOACTIVITY)

- 92 RADIOACTIVE SOLUTIONS

THULIUM

- 46 SPECTROMETRY Solution

TIN

- 77 DEFINING FIXED POINT, ITS-90
- 77 DEFINING FIXED POINT CELLS, ITS-90
- 36 Dibutyltin bis (2-ethylhexanoate) (METALLO-ORGANIC COMPOUNDS)

- 77 DIFFERENTIAL SCANNING CALORIMETRY
- 46 SPECTROMETRY Solution

TIN BASE ALLOYS

- 60 See NONFERROUS METALS

TITANIUM

- 64 GASES IN METALS (NONFERROUS METALS)
- 46 SPECTROMETRY Solution
- 64 TITANIUM BASE ALLOYS (NONFERROUS METALS)

TITANIUM DIOXIDE

- 35 REFRACTORIES

TOXIC METALS

- 17 TOXIC SUBSTANCES IN URINE

TRACE ELEMENTS

- 66 See CERAMICS AND GLASSES
- 29 See FOSSIL FUELS
- 64 See TRACE ELEMENTS IN NICKEL BASE SUPERALLOYS

TRANSMISSION ELECTRON MICROSCOPE

- 36 See THIN FILM FOR TRANSMISSION
- 99 See ASBESTOS

TRANSMITTANCE

- 80 See MOLECULAR ABSORPTION

TRIPLE POINT

- 78 (THERMODYNAMIC PROPERTIES)
- 29 REFERENCE LIQUIDS FOR RATING FUELS

TRIPALMITIN

- 13 HEALTH & CLINICAL
- 44 TRIS(HYDROXYMETHYL)-AMINOMETHANE
- 44 STOICHIOMETRY

TUNGSTEN

- 46 SPECTROMETRY Solution
- 66 Tungsten Carbide (CARBIDES)
- 6 Tungsten Carbide (MICRO-HARDNESS)
- 33 Tungsten Concentrate (ORES)

TURBIDIMETRY (SIZING)

- 2 Portland Cement (CEMENT TURBIDIMETRY AND FINENESS)

U

UNIVERSITY OF PITTSBURGH I (FIRE RESEARCH)

- 4 See SMOKE TOXICITY

URANIUM

- 46 SPECTROMETRY Solution

URANIUM (RADIOACTIVITY)

- 95 Fission Track Glass
- 91 RADIOACTIVE SOLUTIONS
- 95 NATURAL MATRIX MATERIALS

UREA

- 13 HEALTH & CLINICAL
- 76 COMBUSTION CALORIMETRY (THERMODYNAMIC PROPERTIES)
- 14 in Human Serum (SERUM MATERIALS)
- 45 MICROCHEMISTRY

URIC ACID

- 13 HEALTH & CLINICAL

URINE FREEZE-DRIED (FORENSICS)

- 17 Cocaine and Metabolites in
- 17 Cotinine in
- 15 Fluoride in
- 15 Mercury in
- 17 Morphine and Codeine in
- 17 Morphine and Glucuronide in
- 17 Multi-drugs of Abuse in
- 17 THC (Marijuana Metabolite) in
- 17 Toxic Metals in

USA/CANADA COLLABORATIVE MATERIALS

- 9 Bovine Muscle
- 9 Corn Kernel
- 9 Corn Stalk
- 9 Corn Starch
- 10 Durum Wheat Flour
- 9 Hard Red Spring Wheat Flour
- 9 Soft Winter Wheat Flour
- 10 Wheat Gluten
- 10 Whole Egg
- 10 Whole Milk

V

VANADIUM

- 36 Bis(1-phenyl-13-butane-diono)oxovanadium (IV)
- 46 SPECTROMETRY Solution
- 29 Vanadium in Crude Oil (METAL CONSTITUENTS IN FOSSIL FUELS)

VAPOR PRESSURE OF METALS (THERMODYNAMIC PROPERTIES)

- 79 Cadmium
- 79 Gold

VICKERS (MICROHARDNESS) (SURFACE FINISH)

- 6 Bright Copper
- 6 Bright Nickel
- 6 Tungsten Carbide

VISCOSITY OF GLASS (CERAMICS AND GLASSES)

- 88 VISCOSITY FIXPOINTS
- 88 VISCOSITY OF GLASS
- 75 VISCOSITY OF POLYMERS

VITAMINS

- 10 Baby Food Composite (NUTRITION COMPOSITION)
- 10 Cholesterol and FSV in Coconut Oil
- 14 Fat Soluble Vitamins in Human Serum (HEALTH & CLINICAL)
- 10 Infant Formula

VMA

- 13 aka. 4-hydroxy-3-methoxymandelic acid

W

WASPALLOY

- 63 NICKEL BASE ALLOYS (NONFERROUS METALS)
- 22 WATER ANALYSIS (INORGANICS)
- 24 Mercury in Water (METAL CONSTITUENTS IN NATURAL MATRICES)
- 24 Natural Water (METAL CONSTITUENTS IN NATURAL MATRICES)
- 24 Trace Elements in Water (METAL CONSTITUENTS IN NATURAL MATRICES)

WAVELENGTH (OPTICAL PROPERTIES)

- 81 Holmium Oxide Solution

WEAR (SURFACE FINISH)

- 3 D-2 Tool Steel (ABRASIVE WEAR)

WEAR-METALS (ENGINE WEAR MATERIALS)

- 38 WEAR METALS IN OIL

WHALE BLUBBER (ORGANICS) 21

WHEAT FLOUR (FOOD AND AGRICULTURE)

- 9 USA/CANADA COLLABORATIVE MATERIALS

X

XENON (RADIOACTIVITY)

- 93 as Xenon-133 (RADIOPHARMACEUTICALS)

X-RAY

- 89 X-RAY DIFFRACTION
- 89 X-RAY STAGE CALIBRATION

X-RAY FILM

- 82 X-Ray Film Step Tablet (X-RAY AND PHOTOGRAPHY)

Y

YTTERBIUM

- 46 SPECTROMETRY Solution

YTTRIUM

- 46 SPECTROMETRY Solution

Z

ZINC

- 77 DEFINING FIXED POINT, ITS-90
- 77 DEFINING FIXED POINT CELLS, ITS-90
- 77 DIFFERENTIAL SCANNING CALORIMETRY
- 43 METALS (HIGH PURITY METALS)
- 97 Metals on Filter Media (MATERIALS ON FILTER MEDIA)
- 46 SPECTROMETRY Solution
- 37 Zinc Cyclohexanecarboxylate (METALLO-ORGANIC COMPOUNDS)
- 33 Zinc Concentrate (ORES)

ZIRCONIUM

- 46 SPECTROMETRY Solution
- 65 Zircaloy-4 (ZIRCONIUM BASE ALLOYS)

NUMERIC INDEX

SRM	Descriptor	Page	SRM	Descriptor	Page
1c	Limestone, Argillaceous	35	83d	Arsenic Trioxide (Reductometric)	44
4l	Cast Iron	58	84k	Potassium Hydrogen Phthalate	44
5m	Cast Iron	58	87a	Aluminum-Silicon Alloy	60
6g	Cast Iron	58	88b	Dolomitic Limestone	35
7g	Cast Iron (High Phosphorus)	58	89	Glass, Lead Barium	67
11h	Carbon Steel, 0.2 C	51	90	Ferrophosphorus	57
12h	Carbon Steel, 0.4 C	51	92	Low-Boron, Soda-Lime Powder	67
13g	Carbon Steel, 0.6 C	51	93a	High-Boron Borosilicate	67
14g	Carbon Steel, 0.8 C	51	94c	Zinc-Base Die Casting Alloy	65
15h	Carbon Steel, 0.1 C	51	97b	Flint Clay	34
17e	Sucrose (Polarimetric)	44, 82	98b	Plastic Clay	34
19h	Carbon Steel, 0.2 C	51	99a	Feldspar, Soda	34
20g	Carbon Steel	51	100b	LA Steel, Manganese (SAE T340)	53
25d	Manganese Ore	33	101g	Stainless Steel (AISI 304L)	55
30f	LA Steel, Cr-V (SAE 6150)	53	106b	LA Steel, Cr-Mo-Al (Nitalloy rG)	53
32e	LA Steel, Ni-Cr (SAE 3140)	53	107c	Cast Iron (Ni-Cr-Mo)	58
33e	LA Steel, Ni-Mo (SAE 4820)	53	112b	Silicon Carbide	66
36b	LA Steel, Cr-Mo	53	113b	Zinc Concentrate	33
39j	Benzoic Acid (Calorimetric Standard)	76	114p	Portland Cement Fineness Standard	2
45d	Cu Freezing Point	78	115a	Cast Iron (Cu-Ni-Cr)	58
49e	Lead Freezing Point	78	120c	Phosphate Rock (Florida)	11, 33
50c	Tungsten-Chromium-Vanadium Steel	56	121d	Stainless Steel Cr-Ni-Ti (AISI 321)	55
53e	Bearing Metal (84Pb-10Sb-6Sn)	62	122i	Cast Iron	58
54d	Bearing Metal (Tin Base)	64	123c	Stainless Steel Cr-Ni-Nb (AISI 348)	55
57a	Silicon Metal	57	125b	High Silicon Steel – Calcium Bearing	53
58a	Ferrosilicon (73 % Si)	57	126c	High Alloy Steel, High Nickel	54
59a	Ferrosilicon	57	127b	Solder, 40Sn-60Pb	62
64c	Ferrochromium, High Carbon	57	129c	LA Steel, High Sulfur (SAE 112)	53
68c	Ferromanganese, High Carbon	57	131g	LA Steel, High Silicon	53
69b	Bauxite (Arkansas)	33	132b	Tool Steel (AISI M2)	56
70a	Feldspar, Potash	34	134a	Tool Steel, Mo-W-Cr-V	56
72g	LA Steel (AISI 4130)	53	136e	Potassium Dichromate (oxidimetric standard)	44
73c	Stainless Steel, Cr (SAE 420)	55	139b	LA Steel, Cr-Ni-Mo (AISI 8640)	53
76a	Burnt Refractory (Al2O3-40 %)	35	141d	Acetanilide	45
77a	Burnt Refractory, (Al2O3-60 %)	35	142	Anisic Acid	45
78a	Burnt Refractory, (Al2O3-70 %)	35	143d	Cystine	45
79a	Fluorspar, Customs Grade	33	148	Nicotinic Acid	45
81a	Glass Sand	35, 67	152a	Carbon Steel, 0.5 C	51
82b	Cast Iron (Ni-Cr)	58			

SRM	Descriptor	Page	SRM	Descriptor	Page
154c	Titanium Dioxide	35	345a	HA Steel, (Cu Precipitation Hardening)	54
155	LA Steel, Cr-W	53	346a	Valve Steel	54
158a	Bronze, Silicon	61	347	Magnesium Ferrosilicon	57
160b	Stainless Steel Cr-Ni-Mo (AISI 316)	55	348a	Hi Temp. Alloy, (A286) Ni-Cr	54
163	LA Steel, 1.0 C	53	349a	Waspalloy	63
165a	Glass Sand (Low Iron)	35, 67	350a	Benzoic Acid	44
166c	Stainless Steel, Carbon Only	55	351	Sodium Carbonate	44
178	Carbon Steel, 0.4 C	51	352c	Unalloyed Titanium, Hydrogen	65
179	LA Steel, High Silicon	53	360b	Zircaloy 4, Zr-Base Alloy	65
180	Fluorspar, High Grade	33	361	LA Steel (AISI 4340)	54
181	Lithium Ore (Spodumene)	33	362	LA Steel (AISI 94B17) (mod.)	54
182	Lithium Ore (Petalite)	33	363	LA Steel, Cr-V (mod.)	54
183	Lithium Ore (Lepidolite)	33	364	LA Steel, High C (mod.)	54
185h	Potassium Hydrogen Phthalate, pH	71	368	Carbon Steel (AISI 1211)	51
186lg	Potassium Dihydrogen Phosphate	71	395	Unalloyed Copper II (chips)	61
186llg	Disodium Hydrogen Phosphate	71	396	Unalloyed Copper III (chips)	61
187d	Sodium Tetraborate (Borax), pH	71	398	Unalloyed Copper V (chips)	61
188	Potassium Hydrogen Tartrate, pH	71	399	Unalloyed Copper VI (chips)	61
189b	Potassium Tetroxalate, pH	71	400	Unalloyed Copper VII (chips)	61
191c	Sodium Bicarbonate, pH	71	454	Unalloyed Copper XI (chips)	61
192c	Sodium Carbonate, pH	71	457	Unalloyed Copper	61
193	Potassium Nitrate	11	458	Beryllium-Copper (17510)	61
194	Ammonium Dihydrogen Phosphate	11	459	Beryllium-Copper (17200)	61
195	Ferrosilicon (75 % Si-HP Grade)	57	460	Beryllium-Copper (17300)	61
196	Ferrochromium, Low Carbon	57	475	Optical Linewidth	84
198	Silica Brick	35	476	Optical Linewidth	84
199	Silica Brick	35	480	Tungsten-Molybdenum EPMA	35
200a	Potassium Dihydrogen Phosphate	11	481	Gold-Silver EPMA	35
211d	Toluene Liquid Density	89	482	Gold-Copper EPMA	35
276b	Tungsten Carbide	66	494	Unalloyed Copper I (solid)	57
277	Tungsten Concentrate	33	495	Unalloyed Copper II (solid)	61
278	Obsidian Rock	35	496	Unalloyed Copper III (solid)	61
291	LA Steel, Cr-Mo (ASTM A 213)	53	498	Unalloyed Copper V (solid)	61
293	LA Steel, Cr-Ni-Mo (AISI 8620)	53	499	Unalloyed Copper VI (solid)	61
330	Copper Ore Mill Heads	33	500	Unalloyed Copper VII (solid)	61
331	Copper Ore Mill Tails	33	600	Bauxite, Australian	33
334	Gray Cast Iron (Carbon & Sulfur)	58	607	Potassium Feldspar	67
337a	Basic Open Hearth Steel, 1 % Carbon	51	610	Trace Elements in Glass	67
338	White Cast Iron (Carbon & Sulfur)	58	611	Trace Elements in Glass	67
339	Stainless Steel, Cr-Ni-Se (SAE 30)	55	612	Trace Elements in Glass	67
341	Ductile Cast Iron	58	613	Trace Elements in Glass	67
342a	Nodular Cast Iron	58	614	Trace Elements in Glass	67
343a	Stainless Steel (AISI 431)	55	615	Trace Elements in Glass	67
344	HA Steel, (Mo Precipitation Hardening)	54	616	Trace Elements in Glass	67

NUMERIC INDEX

SRM	Descriptor	Page	SRM	Descriptor	Page
617	Trace Elements in Glass	67	688	Basalt Rock	34
620	Soda Lime, Flat	67	689	Ferrochromium Silicon	57
621	Soda-Lime Container	67	690	Iron Ore (Canada)	33
622	Soda-Lime Silica (Durability)	87	691	Iron Oxide, Reduced	33
623	Borosilicate (Durability)	87	692	Iron Ore, Labrador	33
624	Lead-Silica Glass for dc Resistivity	87	693	Iron Ore, Nimba	33
625	Zinc-Base A	65	694	Phosphate Rock, Western	33, 11
626	Zinc-Base B	65	696	Bauxite, Surinam	33
627	Zinc-Base C	65	697	Bauxite, Dominican	33
628	Zinc-Base D	65	698	Bauxite, Jamaican	33
629	Zinc-Base E-ASTM AC 41A	65	699	Alumina (Reduction Grade)	33
630	Zinc-Base F	65	705a	Polystyrene 179k Mol/Wt	74, 76
631	Zinc Spelter (mod)	65	706a	Polystyrene 258k mol/wt	74
640c	Silicon Line Position (XRD)	89	709	Extra Dense Lead	88
641	Titanium Alloy, 8 Mn (A)	64	710a	Soda-Lime Silica Glass	88
642	Titanium Alloy, 8 Mn (B)	64	713	Barium Glass Anneal Pt	88
643	Titanium Alloy, 8 Mn (C)	64	714	Alumina Glass Anneal Pt	88
647	Titanium Alloy, Al-Mo-Sn-Zr	64	716	Neutral Glass Anneal Pt	88
648	Titanium Alloy, Al-Sn-Zr-Cr-Mo	64	717a	Hi Boron Glass Viscosity	88
649	Titanium Alloy V-Al-Cr-Sn	64	720	Sapphire Heat Capacity	76
650	Unalloyed Titanium A	64	723d	Tris (hydroxymethy) amionmethane	44, 71
651	Unalloyed Titanium B	64	726	Selenium, Inter-Purity	43
654b	Titanium Alloy, Al-V	64	728	Zinc, Intermediate Purity	43
656	Silicon Nitride Quantitative Analysis	89	731L1	Borosilicate Glass - Thermal Expansion	78
659	Silicon Nitride, Particle Size	1	731L2	Borosilicate Glass - Thermal Expansion	78
660a	Line Profile LaB6	89	731L3	Borosilicate Glass - Thermal Expansion	78
661	LA Steel (AISI 4340)	52	736L1	Copper Thermal Expansion	78
663	LA Steel, Cr-V (mod.)	52	738	Stainless Steel - Thermal Expansion	78
664	LA Steel, High Carbon, (mod.)	52	740a	Zinc (Freezing Point)	77
670	Rutile Ore	33	741a	Tin (Freezing Point)	77
671	Nickel Oxide 1	63	742	Alumina (Reference Point)	78
672	Nickel Oxide 2	63	743	Mercury (Triple Point)	77
673	Nickel Oxide 3	63	745	Gold-Vapor Pressure	79
674b	X-Ray Powder Diffraction Intensity, set	89	746	Cadmium-Vapor Pressure	79
675	Line Position, Mica (XRD)	89	762	Magnetic Moment Standard Nickel Disk	7
676	Quantitative Analysis, Alumina (XRD)	89	772a	Nickel Sphere for Magnetic Moment	7
679	Brick Clay	34	773	Soda-Lime Silica (Glass Liquidus)	89
680L1a	High Purity Platinum	43	774	Lead-Silica (Dielectric Constant)	87
680L2a	High Purity Platinum	43	781D2	Molybdenum (Heat Capacity)	76
682	High Purity Zinc	43	855a	Aluminum Casting Alloy 356	60
683	Zinc, Metal	43	856a	Aluminum Casting Alloy 380	60
685R	High Purity Gold	43	858	Aluminum Alloy 6011	60
685W	High Purity Gold	43	859	Aluminum Alloy 7075	60

SRM	Descriptor	Page	SRM	Descriptor	Page
861	Nickel-based Superalloy	50	921	Cortisol (Hydrocortisone)	13
862	High Temperature Alloy L-605	54, 60	924a	Lithium Carbonate (Clinical)	13
864	Inconel 600	63	925	VMA (Clinical)	13
865	Inconel 625	63	927c	Bovine Serum Albumin (7 % solution)	14
866	Incoloy, 800	59	928	Lead Nitrate (Clinical)	13
867	Incoloy, 825	59	929	Magnesium Glutamate Dihydrate	13
868	High Temp Alloy Fe-Ni-Co	54	930e	Glass Filters Transmittence	80
869a	LC Column Selectivity	19	931f	Liquid Absorbance Filters UV-VIS	80
870	LC Column Performance	19	934	Clinical Thermometer	79
871	Bronze, Phosphor (CDA521)	61	935a	Potassium Dichromate, UV Absorbance	80
872	Bronze, Phosphor (CDA 544)	61	936a	Quinine Sulfate	81
874	Cupro-Nickel, 10 % (CDA 706) "H-P"	61	937	Iron Metal Clinical	13
875	Cupro-Nickel, 10 % (CDA 706)	61	951	Boric Acid, Assay and Isotopic	44, 48
877	LC Chiral Selectivity	19	952	Boric Acid 95 % enr 10B	44, 48
879	Nickel Silver (CDA 762)	61	953	Cobalt in Aluminum Wire	95
880	Nickel Silver (CDA 770)	61	955b	Lead in Blood	15
882	Alloy Ni-Cu-Al	63	956a	Electrolytes in Frozen Human Serum	14
885	Refined Copper	43	963a	Fission Track Glass U-1 mg/g	95
886	Gold, Ore Refractory	33	965a	Glucose in Human Serum	14
887	Cemented Carbide (W-83,Co-10)	66	966	Toxic Metals in Bovine Blood	14
888	Cemented Carbide (W-64,Co-25,Ta-5)	66	968c	Fat-Sol Vit,Caroten,Cholest in Hum Serum	14
889	Cemented Carbide(W-75,Co-9,Ta-5,Ti-4)	66	970	Ascorbic Acid in Frozen Human Serum	14
890	Cast Iron HC250+V	58	975a	Chlorine (Isotopic)	48
891	Cast Iron, Ni-Hard Type 1	58	976	Copper (Isotopic)	48
892	Cast Iron, Ni-Hard, Type IV	58	977	Bromine (Isotopic)	48
893	Stainless Steel (SAE 405)	55	978a	Silver (Isotopic)	48
895	Stainless Steel (SAE 201)	55	979	Chromium (Isotopic)	48
897	Tracealloy A	64	980	Magnesium (Isotopic)	48
898	Tracealloy B	64	981	Natural Lead (Isotopic)	48
899	Tracealloy C	64	982	Equal Atom Lead (Isotopic)	48
900	Antiepilepsy Drug (4) Level	14	983	Radiogenic Lead (Isotopic)	48
909b	Human Serum	14	984	Rubidium Assay (Isotopic)	48
910	Sodium Pyruvate	13	985	Potassium (Isotopic)	48
911b	Cholesterol	13	986	Nickel (Isotopic)	48
912a	Urea	13	987	Strontium Assay and Isotopic	44, 48
913a	Uric Acid	13	991	Lead-206 Spike Assay and Isotopic	48
914a	Creatinine	13	994	Gallium (Isotopic)	48
915a	Calcium Carbonate (Clinical)	13	997	Thallium (Isotopic)	48
916a	Bilirubin	13	998	Angiotensin I (Human)	13
917b	D-Glucose (Dextrose-Clinical)	13, 44, 82	999a	Potassium Chloride(Assay)	44
918a	Potassium Chloride (Clinical)	13	1001	X-ray Film Step Tablet	82
919a	Sodium Chloride (Clinical)	13	1002d	Hard Board (Surface Flammability)	3
920	D-Mannitol	13	1003c	Glass Spheres (Particle Size)	1

SRM	Descriptor	Page	SRM	Descriptor	Page
1004b	Glass Beads - Particle Size Distribution	1	1111	Red Brass C	62
1006d	Smoke Density, Cellulose	4	1112	Gilding Metal A (disk)	62
1007b	Plastic (Smoke Density)	4	C1112	Gilding Metal A (block)	62
1008	Photographic Step Tablet	82	1113	Gilding Metal B (disk)	62
1010a	Microcopy Test Chart	82	C1113	Gilding Metal B (block)	62
1012	Flooring Radiant Panel	4	1114	Gilding Metal C (disk)	62
1017b	Glass (Particle Size)	1	C1114	Gilding Metal C (block)	62
1018b	Glass (Particle Size)	1	1115	Commercial Bronze A (disk)	62
1019b	Glass (Particle Size)	1	C1115	Commercial Bronze A (block)	62
1021	Glass Beads, Soda Lime	1	1116	Commercial Bronze B (disk)	62
1034	Unalloyed Copper	61	C1116	Commercial Bronze B (block)	62
1035	Leaded-Tin Bronze Alloy	61	1117	Commercial Bronze C (disk)	62
1048	Smoke Toxicity (Cup Furnace)	4	C1117	Commercial Bronze C (block)	62
1049	Smoke Toxicity (Univ of Pittsburgh)	4	C1122	Beryllium-Copper (block)	62
1051b	Barium (Metallo-Organic)	69	1128	Ti Alloy (15V-3AL-3CR-3SN)	64
1052b	Vanadium (Metallo-Organic)	69	1129	Solder (63Sn-37Pb)	62
1053a	Cadmium (Metallo-Organic)	69	1134	LA Steel, High Silicon	52
1057b	Tin (Metallo-Organic)	69	1135	LA Steel, High Silicon	52
1059c	Lead (Metallo-Organic)	69	C1137a	White Cast Iron	57
1060a	Lithium (Metallo-Organic)	69	1138a	Cast Steel (No 1)	57
1065b	Nickel (Metallo-Organic)	69	1139a	Cast Steel (No 2)	57
1066a	Silicon (Metallo-Organic)	69	C1145a	White Cast Iron	57
1069b	Sodium (Metallo-Organic)	37, 69	C1151a	Stainless Steel 23Cr-7Ni	55
1070a	Strontium (Metallo-Organic)	37, 69	C1152a	Stainless Steel 18Cr-11Ni	55
1071b	Phosphorus (Metallo-Organic)	37, 69	C1153a	Stainless Steel 17Cr-9Ni	55
1073b	Zinc (Metallo-Organic)	37, 69	C1154a	Stainless Steel 19Cr-13Ni	55
1075a	Aluminum (Metallo-Organic)	36, 69	1155	Stainless Steel Cr18-Ni12-Mo2 (AISI 316)	55
1077a	Silver (Metallo-Organic)	69	1157	Specialty Steel, Tool (AISI M2)	56
1078b	Chromium (Metallo-Organic)	37, 69	1158	Specialty Steel, High Nickel (36 % Ni)	56
1079b	Iron (Metallo-Organic)	37, 69	1159	Elec/Mag Ni-Fe	63
1080a	Copper (Metallo-Organic)	69	1160	Elec/Mag Ni-Mo-Fe	63
1083	Wear Metals (Base Oil)	38	1171	Stainless Steel Cr17-Ni11-Ti0.3 AISI 321	55
1084a	Wear Metals in Oil, 100 mg/kg	38	1172	Stainless Steel, Cr17-Ni11-Nb.6 AISI 348	55
1085b	Wear Metals in Oil, 300 mg/kg	38	1173	Ni-Cr-Mo-V Steel	57
1089	Steels, Set (consists of SRMs 1095-1099)	59	C1173	Cast Steel 3	57
1090	Ingot Iron, Oxygen	59	1216	Carbon Modified Silica	25
1091a	Stainless Steel (AISI 431)	59	1219	Stainless Steel Cr-Ni (AISI 431)	55
1093	Valve Steel, Oxygen	59	C1221	Carbon Steel	52
1094	Maraging Steel	59	1223	Chromium Steel	55
1104	Fire Cutting Brass	62	1224	LA Steel, Carbon (AISI 1078)	52
1107	Naval Brass B	62	1225	LA Steel AISI 4130	52
1108	Naval Brass C	62			
1110	Red Brass B	62			

SRM	Descriptor	Page	SRM	Descriptor	Page
1226	LA Steel	52	1416	Aluminosilicate Glass for Liquidus Temp	89
1227	LA Steel, Basic Open Hearth, 1 %C	52	1449	Fumed Silica Board	78
1228	LA Steel 0.1 % C	52	1450c	Fibrous Glass Board	78
1230	High Temp Alloy A286	59	1453	Thermal Resis Expanded Polystyrene Board	78
1233	Specialty Steel, Valve Steel	56	1457	Superconducting Nb-Ti Wire	87
1242	High Temp Alloy L-605	60	1459	Fumed Silica Board	78
1243	Waspalloy	63	1473b	Low Density Polyethylene Resin	75
1244	Inconel 600	63	1474	Polyethylene Resin	75
1246	Incoloy 800	59	1475a	Polyethylene, Linear	74, 75
1247	Incoloy 825	59	1478	Polystyrene Narrow Mol Wt	74
C1248	Nickel-Copper Alloy	63	1479	Polystyrene, Narrow Mol Wt	74
1249	Inconel 718	63	1480	Polyurethane	74
1250	High Temp Alloy Fe-Ni-Co	59	1482a	Polyethylene, 14K Molecular Weight	74
1254	LA Steel (Ca only)	52	1483	Polyethylene, Linear	74
1258	Aluminum Alloy 6011	60	1484a	Polyethylene, Linear	74
1259	Aluminum Alloy 7075	60	1486	Bone Meal	15
1262b	LA Steel (AISI 94B17)	52	1487	Poly (methyl methacrylate)	74
1263a	Cr Steel Cr-V (mod)	52	1488	Poly (methyl methacrylate)	74
1264a	LA Steel, High Carbon (mod)	52	1489	Poly (methyl methacrylate)	74
1265a	Electrolytic Iron	52	1491	Arom Hydro/Hexane Toluene	20
1269	Line Pipe (AISI 1521 mod)	52	1492	Chlor Pesticides/Hexane	20
1270	LA Steel, Cr-Mo (A336) (F-22)	52	1493	PCB Congeners	20
1271	LA Steel (HSLA-100)	48	1496	Polyethylene Gas Pipe Resin	75
1276a	Cupro-Nickel (CDA 715)	62	1497	Polyethylene Gas Pipe Resin	75
C1285	LA Steel (A242) (mod)	52	1507b	THC-COOH in Freeze-Dried Urine	17
1286	Low Alloy Steel (HY 80)	52	1508	Benzoylcgonine(Cocaine Meta) FR-DR URINE	17
C1290	High Alloy (HC-250 + V)	57	1511	Multi Drugs of Abuse in Frez-Dried Urine	17
C1291	High Alloy (Ni-Hard, Type I)	57	1514	Thermal Analysis Purity Set (DSC)	77
C1292	High Alloy (Ni-Hard, Type IV)	57	1515	Apple Leaves	11
1295	Stainless Steel (SAE 405)	55	1543	GC/MS System Performance	19
C1296	Stainless Steel	55	1544	Fatty Acids & Chol in Froz Diet Composit	10
1297	Stainless Steel (SAE 201)	55	1546	Meat Homogenate	10
1358a	Cu & Cr Coating on Steel	86	1547	Peach Leaves	11
1359b	Cu & Cr Coating on Steel	86	1548a	Typical Diet	10, 9
1361b	Cu & Cr Coating on Steel	86	1549	Non-Fat Milk Powder	9
1362b	Cu & Cr Coating on Steel	86	1563	Cholesterol & Fat Soluble Vitamins in Coconut Oil	10
1363b	Cu & Cr Coating on Steel	86	1566b	Oyster Tissue	10, 9
1364b	Cu & Cr Coating on Steel	86	1567a	Wheat Flour	9
1400	Bone Ash	15			
1411	Soft Borosilicate Glass	67			
1412	Multicomponent Glass	67			
1413	Glass Sand (High Alumina)	35, 67			

SRM	Descriptor	Page	SRM	Descriptor	Page
1568a	Rice Flour	9	1656	Thianthrene Combustion Calorimeter	76
1570a	Trace Elements in Spinach Leaves	11, 9, 10	1657	Synthetic Refuse Derived Fuel	76
1573a	Tomato Leaves	11	1658a	CH ₄ /Air, 1umol/mol	27
1575a	Trace Elements in Pine Needles	11	1659b	CH ₄ /Air, 10 umol/mol	27
1577b	Bovine Liver	9	1660a	CH ₄ /C ₃ H ₈ /Air 1 umol/mol	27, 28
1580	Shale Oil	21	1661a	SO ₂ /N ₂ 500 umol/mol	28
1582	Petroleum Crude Oil	21	1662a	SO ₂ /N ₂ 1000 umol/mol	28
1584	Phenols in Methanol	20	1663a	SO ₂ /N ₂ 1500 umol/mol	28
1586	Isotope Label Pollutants	20	1664a	SO ₂ /N ₂ 2500 umol/mol	28
1587	Nitro PAH in Methanol	20	1665b	C ₃ H ₈ /Air 3 umol/mol	28
1588a	Organics in Cod Liver Oil	21	1666b	Propane in Air 10 umol/mol	28
1589a	PCBs,Pesti,Dioxins/ Furans in Human Serum	21, 10, 14	1667b	Propane in Air 50 umol/mol	28
1595	Tripalmitin	13	1668b	Propane in Air 100 umol/mol	28
1596	Dinitropyrene Imrs,1Nitropyrene Meth-Chl	20	1669b	Propane in Air 500 umol/mol	28
1597	Complex PAH Mix	21	1671a	CO ₂ /Air, 340 umol/mol	25
1598	Inorganic Constituents in Bovine Serum	14	1672a	CO ₂ /Air, 350 umol/mol	25
1599	2 Anticonvulsant Drugs	14	1674b	CO ₂ /N ₂ mol 7%	26
1614	Dioxin in Isooctane	20	1675b	CO ₂ /N ₂ mol 14%	26
1616a	Sulfur in Kerosene	31	1676	CO ₂ /Air, 365 umol/mol	25
1617a	Sulfur in Kerosene	31	1677c	CO/N ₂ 10 ppm	26
1619b	Sulfur in Residual Fuel Oil 0.7 %	32	1678c	CO/N ₂ 50 umol/mol	26
1620c	Sulfur in Residual Fuel Oil 4 %	32	1679c	CO/N ₂ 100 umol/mol	26
1621e	Sulfur in Residual Fuel Oil 1 %	32	1680b	CO/N ₂ 500 umol/mol	26
1622e	Sulfur in Residual Fuel Oil 2 %	32	1681b	CO/N ₂ 1000 umol/mol	26
1623c	Sulfur in Residual Fuel Oil 0.3 %	32	1683b	NO/N ₂ 50 umol/mol	27
1632c	Trace Elements in Coal	29, 32	1684b	NO/N ₂ 100 umol/mol	27
1633b	Trace Elements in Coal Fly Ash	29	1685b	NO/N ₂ 250 umol/mol	27
1634c	Trace Elements in Fuel Oil	29	1686b	NO/N ₂ 500 umol/mol	27
1635	Trace Elements in Coal (Subbitumimous)	29, 32	1687b	NO/N ₂ 1000 umol/mol	27
1639	Halocarbons (in Methanol)	20	1690	Polystyrene (Particle Size)	1
1640	Natural Water	24	1691	Polystyrene (Particle Size)	1
1641d	Mercury in Water	24	1692	Polystyrene (Particle Size)	1
1643e	Trace Elements in Water	24	1693a	SO ₂ /N ₂ 50 umol/mol	28
1646a	Estuarine Sediment	24	1694a	SO ₂ /N ₂ 100 umol/mol	28
1647d	Priority Pollutant PAHs	20	1696a	SO ₂ /N ₂ , 3500 umol/mol	28
1648	Urban Particulate Matter	22, 23	1710	Aluminum Alloy 3004	60
1649a	Urban Dust/Organics	22	1711	Aluminum Alloy 3004	60
1650b	Diesel Particulate Matter	21	1712	Aluminum Alloy 3004	60
1655	KCl Solution Calorimetry	76	1713	Aluminum Alloy 5182	60
			1714	Aluminum Alloy 5182	60
			1715	Aluminum Alloy 5182	60
			1727	Anode Tin	64

SRM	Descriptor	Page	SRM	Descriptor	Page
1736	Zinc-Aluminum Alloy	65	1836	Nitrogen in Lub Base Oil	37
1737	Zinc-Aluminum Alloy	65	1837	Methanol and Butanol (in Gasoline)	30
1738	Zinc-Aluminum Alloy	65	1838	Ethanol (in Gasoline)	30
1739	Zinc-Aluminum Alloy	65	1839	Methanol (in Gasoline)	30
1740	Zinc-Aluminum Alloy	65	1842	X-Ray Stage Calibration Board (X,Y Dim)	89
1741	Zinc-Aluminum Alloy	65	1843	X-Ray Stage Calibration Board (Z Dim)	89
1742	Zinc-Aluminum Alloy	65	1845	Whole Egg Powder	10
1744	Aluminum (Freezing Point)	77	1846	Infant Formula (milk-based)	10
1745	Indium (Freezing Point)	77	1848	Lubricating Oil Additive Pkg	38
1746	Silver (Freezing Point)	77	1857	Tool Steel for Abrasive Wear	3
1747	Tin Freezing Point Cell	77	1866b	Common Commercial Asbestos	41, 99
1748	Zinc Freezing Point Cell	77	1868	Asbestos in Building Materials	41, 99
1749	Gold vs. Platinum Thermocouple Thermometer	79	1872	Synthetic Glass	35
1750	Standard Platinum Resistance Thermometer	79	1873	Synthetic Glass	35
1751	Gallium Melting Point	78	1876b	Chrysotile Asbestos	41, 99
1754	Steel (AISI 4320)	59	1878a	Respirable Alpha Quartz	39, 89, 97
1755	Low Alloy Steel	55	1879a	Respirable Cristobalite	39, 89, 97
1761	Low Alloy Steel	53	1880a	Portland Cement (Formerly Black)	68
1762	Low Alloy Steel	53	1881a	Portland Cement	68
1763	Low Alloy Steel	53	1882a	Calcium Aluminate Cement	68
1764	Low Alloy Steel	53	1883a	Calcium Aluminate Cement	68
1765	Low Alloy Steel	53	1884a	Portland Cement	68
1766	Low Alloy Steel	53	1885a	Portland Cement	68
1767	Low Alloy Steel	53	1886a	Portland Cement	68
1768	High-Purity Iron	52	1887a	Portland Cement	68
1772	Tool Steel (S-7)	56	1888a	Portland Cement	68
1775	MP 35N Refractory Alloy	60	1889a	Portland Cement	68
1800	Organic Compounds/N ₂	25	1893	Microhardness Cu-Knoop	6
1804c	Organic Compounds/N ₂	25	1894a	Microhardness Ni-Vickers	6
1810a	Linerboard	6	1895	Microhardness Ni-Knoop	6
1815a	n-Heptane (Fuel Rating)	29	1896a	Microhardness Ni-Vickers	6
1816a	Isooctane (Fuel Rating)	29	1897	Specific Surface Area	2
1818a	Chlorine in Lub Base Oil	37	1899	Specific Surface Area for BET	2
1819a	Sulfur in Lub Base Oil	37	1900	Specific Surface Area for BET	2
1827b	Lead Silica Glass Density	89	1905	Microhardness, Ni-Knoop	6
1828b	Ethanol-Water Solution	17	1906	Microhardness, Ni-Knoop	6
1829	Alcohols in Reference Fuel	30	1907	Microhardness, Ni-Knoop	6
1830	Soda Lime Float (Glass)	67	1908	Microhardness, Ni-Vickers	6
1831	Soda Lime Sheet (Glass)	67	1909	Microhardness, Ni-Vickers	6
1834	Fused Ore (Glass)	67	1917	Mercury Porosimeter Instrusion	2
1835	Borate Ore	33	1918	Mercury Porosimeter Instrusion	2

SRM	Descriptor	Page	SRM	Descriptor	Page
1920a	Near IR Reflectance	81	2031b	Metal-on-Quartz Filters	80
1921a	IR Transmiss Wavelength Polystyrene film	81	2032	Potassium Iodide, Stray Light	80
1922	Liquid Refractive Index - Mineral Oil	82	2034	Holmium Oxide Wavelength	81
1923	Poly(ethylene oxide)	74	2035	Near Infrared Transmission Wavelength	81
1924	Poly(ethylene oxide)	74	2036	Near-IR Wavelength/Wavenumber Reflection	81
1930	Glass Filters, Transmittance	80	2037	Red Diesel Dye	81
1932	Fluorescein	81	2040	Bidirectional White Diffuser	81
1935	Potassium Dichromate Soln/ UV Absorbance	80	2046	Transmission Filter	80
1939a	PCBs in River Sediment A	22	2047	Transmission Filter	80
1941b	Organics in Marine Sediment	21	2048	Transmission Filter	80
1944	New York/New Jersey Waterway Sediment	21, 24	2049	Transmission Filter	80
1945	Organics in Whale Blubber	21	2050	Transmission Filter	80
1946	Lake Superior Fish Tissue	21, 24	2051	Transmission Filter	80
1951b	Lipids in Frozen (Liquid) Human Serum	14	2053	IR Transmission Filter	80
1952a	Cholesterol in Human Serum	14	2054	IR Transmission Filter	80
1960	Polystyrene (10 um)	1	2055	IR Transmission Filter	80
1961	Polystyrene (30 um)	1	2056	IR Transmission Filter	80
1963	Polystyrene Spheres	1	2063a	Mineral Glass (Thin Film)	36
1965	Polystyrene (on Slide) (Particle Size)	1	2065	UV-Vis-NIR Transmission Wavelength	81
1967	PT Thermocouple Wire	79	2066	K-411 Glass Microspheres	35
1968	Gallium Melting Point	78	2069b	SEM Performance	85
1969	Rubidium Triple Point	78	2071b	Sinusoidal Roughness	3
1970	Succinonitrile Triple Point	78	2073a	Sinusoidal Roughness	3
1971	Indium Freezing Point	78	2074	Sinusoidal Roughness	3
1972	1, 3-Dioxolan-2-one Triple Point	78	2075	Sinusoidal Roughness	3
1973	N-Docosane Triple Point	78	2084	CMM Probe Performance Standard	6
1974a	Organics-Mussel Tissue (Mytilus edulis)	10, 21	2084R	CMM Probe (10-mm sphere)	6
1975	Diesel Particulate Extract	21	2085	CMM Probe Performance Standard	6
1976	Instrument Sens.for Xray Pwder Diffraction	89	2092	Low-Energy Charpy V-Notch	5
1978	Zirconium Oxide (Particle Size)	1	2096	High-Energy Charpy V-Notch	5
1980	Geothite	73	2098	Super High-Energy Charpy V-Notch	5
1982	Zirconia Thermal Spray Powder	1	2100	Fracture Toughness of Ceramic	7
1984	Thermal Spray Pwder Particle Size Distribution	1	2133	Phosphorus Implant in Silicon Depth Profile	85
1985	Thermal Spray Pwder Particle Size Distribution	1	2134	Arsenic in Silicon	85
2003	First Surface Aluminum on Glass	81	2135c	Ni-Cr Thin Film Depth Profile	85
2017	Multi-Angle White Reflectance	81	2137	B Implant in Si Depth Profile	85
2026	First Surface, Black Glass	81	2139	Zinc-Aluminum Alloy	65
2030a	30% Transmittance	80	2141	Urea	45
			2143	p-Fluorobenzoic Acid	45
			2144	m-Chlorobenzoic Acid	45
			2151	Nicotinic Acid (Combustion Calorimetric Standard)	76

SRM	Descriptor	Page	SRM	Descriptor	Page
2152	Urea (Combustion & Calorimetric Standard)	76	2286	Ethanol (in Gasoline)	30
2159	LA Steel, Carbon & Sulfur Only	54	2287	Ethanol (in Gasoline)	30
2160	LA Steel, Carbon & Sulfur only	54	2288	t-Amyl-methyl-Ether (in Gasoline)	30
2166	LA Steel, F	54	2289	t-Amyl-methyl-Ether (in Gasoline)	30
2167	LA Steel, G	54	2290	Ethyl-t-butyl Ether (in Gasoline)	30
2168	High Purity Iron	54	2291	Ethyl-t-butyl Ether (in Gasoline)	30
2171	LA Steel, (HSLA-100)	53	2292	Methyl-t-Butyl Ether (in Gasoline)	30
2172	S-7 Tool Steel	56	2293	Methyl-t-Butyl Ether (in Gasoline)	30
2175	MP 35N Refractory Alloy	60	2294	Reformulated Fuels (Nominal 11 % MTBE)	30, 31
2181	HEPES Free Acid	13, 72	2295	Reformulated Fuel (Nominal 15 % MTBE)	30, 31
2182	NaHEPESate	13, 72	2296	Reformulated Fuel (Nominal 13 % ETBE)	30, 31
2183	MOPSO Free Acid	13, 72	2297	Reformulated Fuel (Nominal 10 % ETOH)	30, 31
2184	NaMOPSOate	13, 72	2298	Sulfur in Gasoline	31
2185	Pot. Hydrogen Phthalate	72	2299	Sulfur in Gasoline	31
2186I	Potassium Dihydrogen Phosphate	72	2321	Sn-Pb Alloy Coating	85
2186II	Disodium Hydrogen Phosphate	72	2379	Cocaine in Human Hair Segments I	17
2191a	Sodium Bicarbonate	72	2380	Codeine in Human Hair Segments II	17
2192a	Sodium Carbonate	72	2381	Morphine and Codeine in Urine	17
2193	Calcium Carbonate	71	2382	Morphine Glucoronide in Urine	17
2201	Sodium Chloride (Ion-Selective)	72	2383	Baby Food Composite	10
2202	Potassium Chloride (Ion-Selective Electr)	72	2384	Baking Chocolate	10
2203	Potassium Fluoride (Ion-Selective Electr)	72	2385	Slurried Spinach	10
2214	Isooctane Liquid Density	89	2387	Peanut Butter	10
2220	Tin (99.9995%)	77	2389	Amino Acids in 0.1 mol/L Hydrochlor Acid	15
2222	Biphenyl (Differen Scanning Calorimeter)	77	2390	DNA Profiling	17
2225	Mercury (Differen Scanning Calorimeter)	77	2391b	PCR-Based DNA Profiling	17
2232	Indium DSC Calibr Std Temp & Enth of Fus	77	2392	DNA Mitochondrial Sequencing	15
2241	Relative Intensity Correction Standard	81	C2400	HA Steel ACI (17/4 PH)	59
2242	Relative Intensity Correction Standard	81	C2401	HA Steel (ACI-C-4M-Cu)	59
2260a	Aromatic Hydrocarbon in Toluene	20	C2402	Hastelloy 7C	63
2261	Chlorinated Pesticides in Hexane	20	C2415	Battery Lead	63
2262	Chlorinated Biphenyls in Isooctan	20	C2416	Bullet Lead	63
2269	Perdeuterated PAH I	20	C2417	Lead-Base Alloy	63
2270	Perdeuterated PAH II	20	C2418	High-Purity Lead	63
2273	DDT and Metabolites	20	C2423	Ductile Iron A	57
2274	PCB Congeners II	20	C2423a	Ductile Iron B	57
2275	Chlorinated Pesticide II	20	C2424	Ductile Iron C	57
2276	Coplanar PCBs	20	C2424a	Ductile Iron D	57
2285	Arson Text Mixture	30	2426	Galvalume	64
			2430	Scheelite Ore	33

NUMERIC INDEX

SRM	Descriptor	Page	SRM	Descriptor	Page
2431	Titanium Base Alloy	64	2573	Lead Paint Film (Red) Nominal 1.0 mg/cm ²	40, 98
2432	Titanium Base Alloy	64	2574	Lead Paint Film (Gold) Nominal .7 mg/cm ²	40, 98
2433	Titanium Alloy	64	2575	Lead Paint Film (Green) Nominal .3 mg/cm ²	40, 98
2452	Hydrogen in Titanium Alloys	65	2576	Lead Paint Film, High Level	40, 98
2453	Hydrogen in Titanium Alloys	65	2579a	Lead Paint Films for Portable XRF Analyz	40, 98
2454	Hydrogen in Titanium Alloys	65	2580	Powdered Paint Nominal 4 % Lead	40, 98
2490	Non-Newtonian Polymer Solution/Rheology	75	2581	Powdered Paint Nominal 0.5 % Lead	40, 98
2491	Non-Newtonian Polymer Melt for Rheology	75	2582	Powdered Paint Nominal 200 mg/kg L	40, 98
2513	Mode-Field Diameter of Single-Mode Fiber	84	2583	Trace Elements in Indoor Dust	23, 40, 98
2514	Wavelength Reference Absorption Cell-12CO	84	2584	Trace Element in Indoor Dust	23, 40, 98
2515	Wavelength Reference Absorption Cell-13CO	84	2586	Trace Elements in Soil w/lead from paint	24, 40, 98
2517a	Wavelength Reference Absorption Cell	84	2587	Trace Elements in Soil w/Lead from Paint	24, 40, 98
2518	Polarization Mode Dispersion	84	2589	Powdered Paint Nominal 10 % Lead	40, 98
2519	Wavelength Reference Absorption	84	2612a	CO/Air 10 umol/mol	26
2520	Optical Fiber Geometry Standard	84	2613a	CO/Air 20 umol/mol	26
2522	Pin Gage for Optical Fiber Ferrul	84	2614a	CO/Air 45 umol/mol	26
2523	Optical Fiber Ferrule Geometry	84	2619a	Carbon Dioxide in Nitrogen .5 % mol/mol	26
2526	111 p-Type Silicon Resistivity Specimens	83	2620a	Carbon Dioxide in Nitrogen 1.0 % mol/mol	26
2527	111 n-Type Silicon Resistivity Specimens	83	2621a	Carbon Dioxide in Nitrogen 5 % mol/mol	26
2531	Si/SiO ₂ Thickness-50 nm	86	2622a	Carbon Dioxide in Nitrogen 2.0 % mol/mol	26
2532	Si/SiO ₂ Thickness-100 nm	86	2623a	Carbon Dioxide in Nitrogen 2.5 % mol/mol	26
2533	Si/SiO ₂ Thickness-200 nm	86	2624a	Carbon Dioxide in Nitrogen 3.0 % mol/mol	26
2534	Si/SiO ₂ Thickness-25 nm	86	2625a	Carbon Dioxide in Nitrogen 3.5 % umol/mol	26
2535	Si/SiO ₂ Thickness-14 nm	86	2626a	Carbon Dioxide in Nitrogen 4.0 % umol/mol	26
2538	Deterministic Polarization Mode Dispersion	84	2629a	NO/N ₂ , 20 umol/mol	27
2544	Silicon Resistivity	83	2630	NO/N ₂ , 1500 umol/mol	27
2547	Silicon Resistivity	83	2631a	NO/N ₂ , 3,000 umol/mol	27
2551	Oxygen in Silicon	87	2635a	CO/N ₂ 25 umol/mol	26
2553	Optical Fiber Coating Standard	84	2636a	CO/N ₂ 250 umol/mol	26
2554	Optical Fiber Coating Standard	84	2637a	CO/N ₂ 2500 umol/mol	26
2556	Recycled Pellet (Autocatalyst)	25	2638a	CO/N ₂ 5000 umol/mol	26
2557	Recycled Monolith (Autocatalyst)	25	2639a	CO/N ₂ 1.0 % mol/mol	26
2570	Lead Paint Film White/Blank .001 mg/cm ²	40, 98			
2571	Lead Paint Film (Yellow) Nominal 3.5 mg/cm ²	40, 98			
2572	Lead Paint Film (Orange) Nominal 1.6 mg/cm ²	40, 98			

SRM	Descriptor	Page	SRM	Descriptor	Page
2640a	CO/N ₂ 2.0 % mol/mol	26	2721	Moisture & Sulfur in Crude Oil (Yeates Sour)	32
2641a	CO/N ₂ 4 % mol/mol	26	2722	Moisture & Sulfur in Crude Oil (Rufrio Sweet)	32
2642a	CO/N ₂ 8 % mol/mol	26	2723a	Sulfur in Diesel Fuel Oil	31
2643a	Propane in Nitrogen 100 umol/mol	28	2724b	Sulfur in Diesel Fuel Oil, 0.04 %	31
2644a	Propane in Nitrogen 250 umol/mol	28	2730	H ₂ S/N ₂ , 5 umol/mol	27
2645a	Propane in Nitrogen 500 umol/mol	28	2731	H ₂ S/N ₂ , 20 umol/mol	27
2646a	C ₃ H ₈ /N ₂ , 1000 umol/mol	28	2735	NO/N ₂ , 800 umol/mol	27
2647a	C ₃ H ₈ /N ₂ , 2500 umol/mol	28	2736a	NO/N ₂ , 2000 umol/mol	27
2648a	C ₃ H ₈ /N ₂ , 5000 umol/mol	28	2737	NO/N ₂	27
2657a	O ₂ /N ₂ 2 % mol/mol	28	2738	NO/N ₂	27
2658a	O ₂ /N ₂ 10 % mol/mol	28	2740a	CO/N ₂ , 10 % mol/mol	27
2659a	O ₂ /N ₂ , 21 % mol/mol	28	2741a	CO/N ₂ , 13 % mol/mol	27
2660a	Total Oxides of Nitr in Air 100 umol/mol	27	2745	CO ₂ /N ₂ , 16 % mol/mol	26
2670a	Toxic Elements in Urine	15	2750	CH ₄ /Air 50 umol/mol	27
2671a	Fluoride in Freeze-Dried Urine	15	2751	CH ₄ /Air 100 umol/mol	27
2672a	Mercury in Urine	15	2764	C ₃ H ₈ /Air .25 umol/mol	28
2678	Membrane Blank Filter	39, 97	2775	Foundry Coke	31
2679a	Quartz on Filter Media	38, 97	2776	Furnace Coke	31
2681	Ashless Blank Filter	39, 97	2780	Hard Rock Mine Waste	24
2682b	Sulfur & Mercury in Coal	32, 76	2781	Domestic Sludge	24
2683b	Sulfur in Coal, 2 %	32, 76	2782	Industrial Sludge	24
2684b	Sulfur & Mercury in Coal	32, 76	2783	Air Particulate on Filter Media	23, 38, 97
2685b	Sulfur & Mercury in Coal	32, 76	2798a	Microhardness Ni-Vickers	6
2686	Portland Cement Clinker	68	2800	Microscope Magnification Standard	84, 85
2687	Portland Cement Clinker	68	2806	Medium Test Dust(MTD) in Hydraulic Fluid	2
2688	Portland Cement Clinker	68	2810	Rockwell C Hardness, Low	5
2689	Coal Fly Ash	29	2811	Rockwell C Hardness, Mid	5
2690	Coal Fly Ash	29	2812	Rockwell C Hardness, High	5
2691	Coal Fly Ash	29	2830	Microhardness, Ceramic-Knoop	6
2692b	Sulfur & Mercury in Coal	32, 76	2831	Microhardness, Ceramic-Vickers	6
2693	Low Sulfur/Mercury Coal	32	2853	Magnetic Moment Standard - Yttrium Iron Garnet	7
2695	Fluoride in Vegetation	11	2885	Polyethylene (Molar Mass 6,280 g/mol)	74
2702	Marine Sediment	24	2886	Polyethylene (Molar Mass 87,000 g/mol)	74
2703	Sediment for Solid Sampling	24	2887	Polyethylene (Molar Mass 196,400 g/mol)	74
2709	San Joaquin Soil	24	2888	Polyethylene/Polystyrene	74
2710	Montana I Soil	24	2890	Water Saturated Octanol	32
2711	Montana II Soil	24	2910	Calcium Hydroxyapatite	15, 89
2713	Lead in Reference Fuel	29	2921	Cardiac Troponin	15
2714	Lead in Reference Fuel	29			
2717a	Sulfur in Residual Fuel Oil	32			
2718	Green Petroleum Coke	29, 31			
2719	Calcined Petroleum Coke	29, 31			

SRM	Descriptor	Page	SRM	Descriptor	Page
2930	Ultimate Range Visible Absorbance Filters	80	3006	Carbon Tetrachloride in Methanol	22
2950	Respirable Alpha Quartz on Filter Media	39, 97	3008	Methylene Chloride in Methanol	22
2951	Respirable Alpha Quartz on Filter Media	39, 97	3009	1,2 Dichloropropane in Methanol	22
2952	Respirable Alpha Quartz on Filter Media	39, 97	3010	Tetrachloroethylene in Methanol	22
2953	Respirable Alpha Quartz on Filter Media	39, 97	3011	1,1,1 Trichloroethane in Methanol	22
2954	Respirable Alpha Quartz on Filter Media	39, 97	3012	1,2-Dichloroethane in Methanol	22
2955	Respirable Alpha Quartz on Filter Media	39, 97	3014	1,2,3 Trichloropropane in Methanol	22
2956	Respirable Alpha Quartz on Filter Media	39, 97	3015	Isopropylbenzene in Methanol	22
2957	Respirable Alpha Quartz on Filter Media	39, 97	3016	sec-Butylbenzene in Methanol	22
2958	Respirable Alpha Quartz on Filter Media	39, 97	3063	Dioxin in Methanol	22
2960	Respirable Alpha Cristobalite on Filter Media	39, 97	3064	Endothall in Water	22
2961	Respirable Alpha Cristobalite on Filter Media	39, 97	3067	Toxaphene in Methanol	23
2962	Respirable Alpha Cristobalite on Filter Media	39, 98	3068	Chlordane in Methanol	23
2963	Respirable Alpha Cristobalite on Filter Media	39, 98	3071	Glyphosate	23
2964	Respirable Alpha Cristobalite on Filter Media	39, 98	3072	Diquat Dibromide Monohydrate in Water	23
2965	Respirable Alpha Cristobalite on Filter Media	39, 98	3075	Aroclor 1016 in Transformer Oil	23
2966	Respirable Alpha Cristobalite on Filter Media	39, 98	3076	Aroclor 1232 in Transformer Oil	23
2967	Respirable Alpha Cristobalite on Filter Media	39, 98	3077	Aroclor 1242 in Transformer Oil	23
2975	Diesel Partic.Matter (Indus.Forklift)	21	3078	Aroclor 1248 in Transformer Oil	23
2976	Mussel Tissue T.E. & Methylmercury Frz-Dr	21	3079	Aroclor 1254 in Transformer Oil	23
2977	Mussel Tissue Organic Contaminants &T.E.	21	3080	Aroclor 1260 in Transformer Oil	23
2978	Mussel Tissue Org.Contam Raritan Bay, NJ	21	3081	Aroclor 1016 in Methanol	23
3000	Benzene in Methanol	22	3082	Aroclor 1232 in Methanol	23
3001	Toluene in Methanol	22	3083	Aroclor 1242 in Methanol	23
3002	Ethylbenzene in Methanol	22	3084	Aroclor 1248 in Methanol	23
3003	o-Xylene in Methanol	22	3085	Aroclor 1254 in Methanol	23
3004	m-Xylene in Methanol	22	3086	Aroclor 1260 in Methanol	23
3005	p-Xylene in Methanol	22	3090	Aroclors in Transformer Oil (set SRMs 3075-3080)	23
			3091	Aroclors in Methanol (set SRMs 3081 - 3086)	23
			3101a	Aluminum Standard Solution	46
			3102a	Antimony Standard Solution	46
			3103a	Arsenic Standard Solution	46
			3104a	Barium Standard Solution	46
			3105a	Beryllium Standard Solution	46
			3106	Bismuth Standard Solution	46
			3107	Boron Standard Solution	46
			3108	Cadmium Standard Solution	46
			3109a	Calcium Standard Solution	46
			3110	Cerium Standard Solution	46
			3111a	Cesium Standard Solution	46
			3112a	Chromium Standard Solution	46
			3113	Cobalt Standard Solution	46

SRM	Descriptor	Page	SRM	Descriptor	Page
3114	Copper Standard Solution	46	3159	Thorium Standard Solution	47
3115a	Dysprosium Standard Solution	46	3160a	Thulium Standard Solution	47
3116a	Erbium Standard Solution	46	3161a	Tin Standard Solution	47
3117a	Europium Standard Solution	46	3162a	Titanium Standard Solution	47
3118a	Gadolinium Standard Solution	46	3163	Tungsten Standard Solution	47
3119a	Gallium Standard Solution	46	3164	Uranium Standard Solution	47
3120a	Germanium Standard Solution	46	3165	Vanadium Standard Solution	47
3121	Gold Standard Solution	46	3166a	Ytterbium Standard Solution	47
3122	Hafnium Standard Solution	46	3167a	Yttrium Standard Solution	47
3123a	Holmium Standard Solution	46	3168a	Zinc Standard Solution	47
3124a	Indium Standard Solution	46	3169	Zirconium Standard Solution	47
3126a	Iron Standard Solution	46	3181	Sulfate Anion Solution	48
3127a	Lanthanum Standard Solution	46	3182	Chloride Anion Solution	48
3128	Lead Standard Solution	46	3183	Fluoride Anion Solution	48
3129a	Lithium Standard Solution	46	3184	Bromide Anion Solution	48
3130a	Lutetium Standard Solution	46	3185	Nitrate Anion Solution	48
3131a	Magnesium Standard Solution	46	3186	Phosphate Anion Solution	48
3132	Manganese Standard Solution	46	3190	Aqueous Electrolytic Conductivity 25 uS/cm	73
3133	Mercury Standard Solution	47	3191	Aqueous Electrolytic Conductivity 100 uS/cm	73
3134	Molybdenum Standard Solution	47	3192	Aqueous Electrolytic Conductivity 500 uS/cm	73
3135a	Neodymium Standard Solution	47	3193	Aqueous Electrolytic Conductivity 1000 uS/cm	73
3136	Nickel Standard Solution	47	3194	Aqueous Electrolytic Conductivity 10,000 uS/cm	73
3137	Niobium Standard Solution	47	3195	Aqueous Electrolytic Conductivity 100,000 uS/cm	73
3138	Palladium Standard Solution	47	3196	Aqueous Electrolytic Conductivity 20,000 uS/cm	73
3139a	Phosphorus Standard Solution	47	3198	Aqueous Electrolytic Conductivity 5 uS/cm	73
3140	Platinum Standard Solution	47	3199	Aqueous Electrolytic Conductivity 15 uS/cm	73
3141a	Potassium Standard Solution	47	3230	Iodine-129, Isotopic (low levels)	48
3142a	Praseodymium Standard Solution	47	3231	Iodine-129, Isotopic (high levels)	48
3143	Rhenium Standard Solution	47	3240	Ephedra-Related Dietary Supplement	10
3144	Rhodium Standard Solution	47	3241	Ephedra-Related Dietary Supplement	10
3145a	Rubidium Standard Solution	47	3242	Ephedra-Related Dietary Supplement	10
3147a	Samarium Standard Solution	47	3243	Ephedra-Related Dietary Supplement	10
3148a	Scandium Standard Solution	47	3244	Ephedra-Related Dietary Supplement	10
3149	Selenium Standard Solution	47	4201B	Niobium-94 Point Source	93
3150	Silicon Standard Solution	47	4203D	Colbalt-60	93
3151	Silver Standard Solution	47	4218F	Europium-152 Point Source	93
3152a	Sodium Standard Solution	47	4222C	Carbon-14 (as hexadene)	91
3153a	Strontium Standard Solution	47			
3154	Sulfur Standard Solution	47			
3155	Tantalum Standard Solution	47			
3156	Tellurium Standard Solution	47			
3157a	Terbium Standard Solution	47			
3158	Thallium Standard Solution	47			

SRM	Descriptor	Page	SRM	Descriptor	Page
4226C	Nickel-63 Solution	91	4427	Yttrium-90 Solution (Lot 5)	93
4234A	Strontium/Yttrium-90	92	4915E	Cobalt-60 Solution	91
4241C	Barium-133 Point Source	93	4919H	Strontium-90 Solution	92
4251C	Barium-133 Solution	91	4926E	Hydrogen-3 Water	91
4288A	Technetium-99	92	4927F	Hydrogen-3 Water	91
4320A	Curium-244 Solution	91	4941	Neptunium-237	91
4321C	Natural Uranium Solution	92	4943	Chlorine-36 Solution	91
4322B	Americium-241 Solution	91	4947C	Hydrogen-3 Toluene	91
4323B	Plutonium-238 Solution	91	4949C	Iodine-129 Solution	91
4324B	Uranium-232	92	4965	Radium-226 Solution	92
4325	Beryllium-10/9 Solution	94	4966	Radium-226 Solution	92
4326	Polonium-209 Solution	91	4967	Radium-226 Solution	92
4328C	Thorium-299	92	4968	Radium-226/Radon-222 Capsule	94
4329	Curium-243 Solution	91	4969	Radium-226 Solution	92
4330B	Plutonium-239 Solution	91	4990C	Oxalic Acid Powder	94
4332D	Americium-243 Solution	91	8010	Sand for Sand Sieve Analysis	1
4334G	Plutonium-242 Solution	91	8040	Sodium Oxalate	44
4338A	Plutonium-240 Solution	91	8091	SEM Sharpness Standard	85
4339B	Radium-228 Solution	92	8107	Gun Powder	63
4340B	Plutonium-241 Solution	91	8407	Buffalo River Sediment	24
4341	Neptunium-237 Solution	91	8411	Mixed Asbestos Research Filter	41, 99
4342A	Thorium-230	92	8412	Corn Stalk (Zea Mays)	9, 11
4350B	River Sediment (Radioactivity)	95	8413	Corn Kernel (Zea Mays)	9, 11
4351	Human Lung Powder	95	8414	Bovine Muscle Powder (Beef)	9
4352	Human Liver Powder	95	8415	Whole Egg Powder	10
4354	Lake Sediment Powder	95	8418	Wheat Gluten	10
4355	Peruvian Soil Powder	95	8420	Iron Electrolytic	79, 83
4356	Ashed Bone (Radioactivity)	95	8421	Iron Electrolytic	79, 83
4357	Ocean Sediment Powder	95	8424	Graphite Thermal Conductivity	79, 83
4361C	Hydrogen-3 Water	91	8426	Graphite Thermal Conductivity	79, 83
4370C	Europium-152 Solution	91	8432	Corn Starch	10
4401	Iodine-131 Solution	93	8433	Corn Bran	10
4404	Thallium-201	93	8435	Whole Milk Powder	10
4407	Iodine-125 Solution	93	8436	Durum Wheat Flour	9, 10
4410	Technetium-99m	93	8437	Hard Red Spring Wheat Flour	9
4412	Molybdenum-99 Solution	93	8438	Soft Winter Wheat Flour	9
4415	Xenon-133 Solution	93	8441	Wheat Hardness	9
4416	Gallium-67 Solution	93	8443	GC/MS System Performance	19
4417	Indium-111	93	8444	Cotinine in Freeze Dried Human Urine	17
4425	Samarium-153	93	8455	Pyrite Ore	34

SRM	Descriptor	Page			
8456	Ultra-hi Molecular Wt. Polyethylene Bar	15, 75	8549	IAEA-N3-Potassium Nitrate	49
8457	Ultra-hi Molecular Wt. Polyethylene Bar	75	8550	USGS25-Ammonium Sulfate	49
8458	Artificial Flaw for Eddy Current	5	8551	USGS26-Ammonium Sulfate	49
8466	Y-HCH (Lindane)(neat)	20	8552	NSVEC-Gaseous Nitrogen	49
8467	4, 4'-DDE (neat)	20	8553	Soufre de Lacq-Elemental Sulfur	49
8469	Pesticide, 4,4'-DDT (neat)	20	8554	NZ1-Silver Sulfide	49
8480	Secondary Ferrite # Standard - Low Range	7	8555	NZ2-Silver Sulfide	49
8481	Secondary Ferrite # Standard - High Range	7	8556	NBS123-Sphalerite	49
8491	Sugar Cane Bagasse	11	8557	NBS127-Barium Sulfate	49
8492	Eastern Cottonwood	11	8558	USGS32-Potassium Nitrate	49
8493	Monterey Pine	11	8559	Natural Gas Isotopic	49
8494	Wheat Straw	11	8560	Natural Gas Isotopic	49
8495	Northern Softwood	7	8561	Natural Gas Isotopic	49
8496	Eucalyptus Hardwood	7	8562	CO2-Heavy, Paleomarine Origin	49
8505	Vanadium in Crude Oil	29	8563	CO2-Light, Paleomarine Origin	49
8506a	Transformer Oil	32	8564	CO2-Biogenic, Modern Biomass Origin	49
8507	Mineral Oil	32	8590	High Sulfur Gas Oil Feed	38
8509	Moisture in Methanol, 93 mg/kg	32	8600	Chinese Copper Ore	34
8510	Moisture in Methanol, 325 mg/kg	32	8601	Chinese Copper Ore	34
8535	Vsmow-Water	49	8602	Chinese Lead Ore	34
8536	GISP-Water	49	8603	Chinese Lead Ore	34
8537	SLAP-Water Light Stable Isotopic Std	49	8604	Chinese Zinc Ore	34
8538	NBS30-Biotite	49	8605	Chinese Molybdenum Ore	34
8539	NBS22-Oil	49	8606	Chinese Molybdenum Ore	34
8540	PEFI-Polyethylene Foil	49	8607	Chinese Tungsten Ore	34
8541	USGS24-Graphite	49	8608	Chinese Tungsten Ore	34
8542	Sucrose ANU-Sucrose	49	8631	Medium Test Dust (MTD)	2
8543	NBS18-Carbonatite	49	8632	Ultrafine Test Dust	2
8544	NBS19-Limestone	49	8640	Fluorescein Labeled Microbead	81
8545	LSVEC-Lithium Carbonate	49	8680	Paint on Fiberboard	40, 98
8546	NBS28-Silica Sand	49	8704	Buffalo River Sediment	24
8547	IAEAN1-Ammonium Sulfate	49	8759	ICTA Set DTA	77
8548	IAEAN2-Ammonium Sulfate	49	8760	ICTA Set DTA	77
			8785	Particulate Matter on Filters	24, 38, 97
			GM754	ICTA Polystyrene DTA	77
			RM5	Cu Low Temperature Heat Capacity	76

NIST develops and promotes measurements, standards, and technology to enhance productivity, facilitate trade, and improve the quality of life. As the U.S. National Metrology Institute, NIST continually strives to meet the nation's measurement needs with Standard Reference Materials, Calibration Services, and Standard Reference Data. Please visit our website at www.nist.gov for further information.



**National Institute of
Standards and Technology**
Technology Administration
U.S. Department of Commerce